

## Product datasheet for **MG222098**

### **Hoxa1 (NM\_010449) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Hoxa1 (NM\_010449) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Hoxa1  
**Synonyms:** ERA1; Hox-1.6  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG222098 representing NM\_010449  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGACAATGCAAGAATGAACTCCTTTCTGGAATACCCATCCTTGGCAGTGGCGACTCTGGACCTGCT  
CGGCGCGAGCTTACCCCTCTGACCATGGGATTACAACCTTCCAATCCTGCGCGGTCAGTGCCAACAGCTG  
CGGCGGCGACGACCGCTTCTAGTGGGCAGGGGGTGCAGATCAGCTCGCCCCACCACCACCACCACC  
CACCACCACCATCACCACCCAGACGGCTACTACCAGACTTCTGAAACCTTGGGATTTCTTATTCCC  
ACTCGAGTTGTGGTCCAAGCTATGGCGCGCAGAACTTCAGTGCGCCTTATGGCCCTATGGATTAATCA  
GGAAGCAGACGTAAGTGGTGGGTACCCCGTGCCTCCCGCTGTTACTCTGAAACCTCTCGACTCCC  
ATGGTCCAGCATCACCACCACCAGGGTTATGCTGGGGCAGTGGGCTCGCCCCAGTACATTCACC  
ACTCATATGGACAAGAGCAGCAGACTCTGGCCCTGGCCACGTATAATAACTCCTTATCCCTCTCCACGC  
CAGCCACCAAGAAGCCTGTGTTCCCTGCTTCAGAGACGCTTCTCCAGCGCAGACCTTGGACTGGATG  
AAAGTTAAAAGAAACCTCCAAAACAGGGAAAGTTGGAGAGTACGGCTACGTGGGTCAACCAACGCAG  
TGCGCACCAATTCACCACCAAGCAGCTCACAGAGCTGGAGAAGGAGTTCACCTCAACAAGTACCTTAC  
ACGAGCGCGCAGGGTGGAGATTGCCGCTCCCTACAGCTCAATGAGACCCAGGTGAAGATCTGGTTCCAG  
AATCGCCGATGAAGCAGAAGAAGCGTGAAGAAGGAGGGCTCCTGCCATCTCCCTGCCACTCCTCCTG  
GCAGCGATGAGAAAACGGAAGAATCATCTGAGAAATCTAGCCCTCGCCAGTGCCCTTCTCCGGCATC  
GTCTACCTCAGACTCTGACTACCTCCAC

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG222098 representing NM\_010449  
 Red=Cloning site Green=Tags(s)

```
MDNARMNSFLEYPILGSGDSGTCSARAYPSDHGITTFFQSCAVSANS CGGDDRFLVGRGVQISSPHHHHHH
HHHHHHHPQTATYQTSNGLGISYSHSSCGPSYGAQNF SAPYGPYGLNQEADVSGGYPPCAPAVYSNLSSTP
MVQHHHHHQGYAGGTVGSPOYIHHSYGGQEQQTLALATYNNLSPLHASHQEACRSPASETSSPAQTFDWM
KVKRNPPKTGKVG EYGYVGPNAVRTNFTTKQL TELEKEFHFNKYLTRARRVEIAASLQLNETQVKIWFQ
NRRMKQKKREKGLLLPISPATPPGSDKTEESSEKSSSPSPASPASSTSDTLTTSH
```


TRTRPLE - GFP Tag - V

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja2354\\_b02.zip](https://cdn.origene.com/chromatograms/ja2354_b02.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_010449


**ORF Size:** 1011 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_010449.4</a> , <a href="#">NP_034579.3</a>
<b>RefSeq Size:</b>	2474 bp
<b>RefSeq ORF:</b>	1011 bp
<b>Locus ID:</b>	15394
<b>Cytogenetics:</b>	6 25.4 cM

### Product images:



Circular map for MG222098