

Product datasheet for **MG222122**

Snap23 (NM_001177792) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Snap23 (NM_001177792) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Snap23
Synonyms:	23kDa; AA408749; SNAP-23; Sndt; Syndet
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG222122 representing NM_001177792 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGC**C

ATGGATAATCTGTCCCCAGAGGAAGTTCAGCTGCGGGCTCACCAGGTTACTGATGAGTCTCTGGAAAGCA
 CAAGGAGAATCCTGGGTTTAGCCATTGAGTCTCAGGATGCAGGAATCAAGACCATCACTATGCTGGATGA
 GCAAGGGGAACAATAATCGCATAGAAGAAGGCATGGACCAATAAATAAGACATGAGAGAGGCAGAG
 AAGACTTTAACAGAACTCAACAAGTGTGTGGCCTCTGCATCTGCCCTTGAATAGATTCTCAGTTGGTG
 ACTGTTTCTTTGAAACCAGGACAAAGAAGTCTGAGTCTGGAAAGAACTATAAGGCAACATGGGGAGATGG
 TGGCGACAACCTACCTAGCAATGTAGTATCTAAGCAACCGAGCCGATTACAAATGGTCAGCCTCAGCAA
 ACTACAGGAGCAGCCAGTGGTGGATACATTAAACGTATACTAATGATGCCAGAGAAGATGAGATGGAAG
 AGAACCTGACTCAAGTGGGCAGCATCCTAGGCAACCTAAAGAACATGGCTCTGGATATGGGCAATGAAAT
 TGATGCTCAAAACCAGCAATACAGAAGATCACAGAAAAGGCTGACACCAACAAGAATCGTATTGACATT
 GCCAATACAAGAGCAAAGAACTCATTGATAGT

ACGCGTACGCGGCGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


Protein Sequence: >MG222122 representing NM_001177792
Red=Cloning site Green=Tag(s)

MDNLSPEEVQLRAHQVTDESLESTRILGLAIESQDAGIKTITMLDEQGEQLNRIEEGMDQINKDMREAE
KTLTELNKKCCGLCICPCNRFVSGDCFFETRTKNFESGKNYKATWGDGDNPSNVVSKQPSRITNGQPQQ
TTGAASGGYIKRITNDAREDEMEENLTQVGSILGNLKNMALDMGNEIDAQNNQIQKITEKADTNKNRIDI
ANTRAKKLIDS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001177792

ORF Size: 663 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001177792.1, NP_001171263.1</u>
RefSeq Size:	2223 bp
RefSeq ORF:	666 bp
Locus ID:	20619
UniProt ID:	<u>O09044</u>
Cytogenetics:	2 E5
Gene Summary:	Essential component of the high affinity receptor for the general membrane fusion machinery and an important regulator of transport vesicle docking and fusion.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG222122