

Product datasheet for **MG222138**

Lpar5 (NM_001163268) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lpar5 (NM_001163268) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Lpar5
Synonyms:	Gm1072; Gpr92; GPR93; LPA5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG222138 representing NM_001163268 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTTGCCAATTCCTCAGCCAACACGACTTCTACCAACAGCTCTGTGCTCCAGTGCCCTGACTATCGAG
ATACACATCGTTTGCATATGGTGGTCTACAGCCTGGTATTGGCGACTGGTCTCCCTCTCAACGCTCTGGC
TCTCTGGGTCTTCTGCGTGTACTGCGGTACACTCAGTGGTGGCGGTACATGTGCAACCTGGCAGCC
AGCGACTTGCTCTTACCCTGTCACTTCCCTGCGCCTCTCTACTATGCACAGCACCCTGGCCTTTTC
CAGGCTTCTGTGCCAGACGTGCGGCCATCTCCAGATGAACATGTACGGCAGCTGTCTTTTCTGAT
GCTCATCAACGTGGACCGCTATGCGGCCATCGTGCACCCGCTGAGACTGCGCCACCTACGGCGGCCCGT
GTGGCACGGCGGCTCTGCCTGGGCGTGTGGGCTCTCATCTGCTGTTTCGCTGTGCCCGCCCGCGGTGC
ACAGCCCGTCCCCTGCACGTACAAGAATCACTGTGCGCCTGTGCTTCGAGAGCTTCAGCGATGAACT
GTGGAAGGGCAGGCTGCTGCCGCTCCTGCTGCTGGCCGAGATACTAGGCTTTCTGCTGCCCTGGCGGT
GTCGTCTATTGCTGCGCAGAGTCTTCTGGACTGGCGAGGCCGACTCAGAGCCAACGGCGAC
GGAAGACCGTGGCCTCCTGCTGGCCAATCTCATCTTCTGCTGTGCTTCGTCGCCCTATAACTCCAC
GCTGGCTGATATGGGTGCTACGGGCCAACTTGGTGAAGAAGCAGCATTACAGACCGCGATCAGGTGGC
GGGTGCTGATGATAATGGTGTGCTGGCGGCCAACTGCGTGGTCCACTGGTTTACTACTTCA
GTGCCGAGGGTTTCCGTAACACCCTTCGCAACCTGGGCGCCCGCTGAATACCAGGCTTTGGCTACCAA
TGGGGTGCAGGCGTGTCCACCGAACTACCCCTCAGAAAAGCACCCAAAACACTGGGCGAGGATGCCACAAGT
CAGGTTCTACTCCAGCCTGCCACTCTGGGTACACCCCGGACAACCTGCTCCAGGATTTCGGCTCTCTG

ACGGTACGGCGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG222138 representing NM_001163268
 Red=Cloning site Green=Tags(s)

MFANSSANTTSTNSSVLQCPDYRDTHRLHMVVYSLVLTATGLPLNALALWVFLRVLRVHVSVSVYMCNLAA
 SDLLFTLSLPLRLSYAQHHWPFPGFLCQTSGAIFQNMNYGSCLFLMLINVDRYAAIVHPLRLRHLRRPR
 VARRLCLGVWALILLFAVPAARVHSPSHCTYKNIIVRLCFESFSDELWKGRLPLLLLLAEILGFLPLAA
 VVYSSGRVFWTLARPDATQSQRRTKTVRLLLANLIIFLLCFVPYNSTLAVYGLLRANLVKNSIQDRDQVR
 GVLMIIVLLAGANCVLDPLVYVYFSAEGFRNTRLNLAGPLNTRPLATNGAAGVLTPELSESTQNTGQDQATS
 QVLLQPATLGTTPDNCSQDSAL


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001163268

ORF Size:	1152 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001163268.1 , NP_001156740.1
RefSeq Size:	1458 bp
RefSeq ORF:	1155 bp
Locus ID:	381810
UniProt ID:	Q149R9
Cytogenetics:	6 F2
Gene Summary:	Receptor for lysophosphatidic acid (LPA), a mediator of diverse cellular activities. [UniProtKB/Swiss-Prot Function]