

Product datasheet for **MG222198**

Tnnt3 (NM_001163666) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tnnt3 (NM_001163666) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tnnt3
Synonyms:	fTnT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG222198 representing NM_001163666 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTGACGAGGAACTGAACAAGTTGAGGAACAATACGAAGAGGAAGAGGAGGCCAAGAGGAAGAGG
AAGTCCAAGAAGATGCTGTCGCCGAGGAGGACGGGAGGAGGACGAGGAAGAGGAGAAACCAAGACCCAA
ACTTACTGCTCCTAAGATCCCGAAGGGGAGAAAGTAGACTTTGATGACATCCAGAAGAAGCGTCAAAAC
AAGGACCTCATGGAGCTCAAGCCCTCATTGACAGCCACTTTGAAGCTAGAAAGAAAGAAGAAGAGGAGC
TAATCGCACTTAAGAAAGAATTGAGAAGCGCCGTGCAGAGAGGGCTGAGCAGCAAGAATTCCGCGCTGA
GAAGGAGCGGGAACGCCAGAACAGATTGGCGGAGGAGAAGGCCAGAAGAGAGGAGGAGGATGCCAAGAGG
AGAGCTGAAGATGACATGAAGAAGAAAAAGGCTCTGTCCTCCATGGGTGCCAACTACAGCAGCTACCTGG
CCAAGGCTGACCAGAAGAGAGGCAAGAAACAGACAGCCCGGAAAATGAAAAAGAAGATTCTTGCCGAAA
GCGCAAGCCTCTGAACATTGACCATCTTAGCGATGACAAGCTGAGGGACAAGGCCAAGGAACTCTGGGAT
ACCTTGTACCAACTGGAGACTGACAAATTCGAGTTTGGGGAGAAGCTGAAACGTCAGAAATATGATATCA
CCACCCCTCAGGAGCCGCATTGACCAAGCCAGAAGCACAGCAAGAAGGCCGGTGCCACAGCCAAGGGCAA
AGTCGGTGGGCGCTGGAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG222198 representing NM_001163666
 Red=Cloning site Green=Tags(s)

MSDEETEQVEEQYEEEEEAQEEEEVQEDAVAEEREEDDEEEKPRPKLTAPKIPEGEKVDFDDIQKKRQN
 KDLMELQALIDSHFEARKKEEEELIALKERIEKRRAERAQQRIRAEKERERQNRLAEEKARREEDAKR
 RAEDDMKKKKALSSMGANYSSYLAKADQKRGKKQTAREMKKKILAERRKPLNIDHLSDDKLRDKAKELWD
 TLYQLETDKFEFGKELKRQKYDITTLRSRIDQAQKHSKKAGATAKGVGGRWK

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001163666

ORF Size: 789 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001163666.1](#), [NP_001157138.1](#)

RefSeq Size: 1081 bp


RefSeq ORF: 792 bp

Locus ID: 21957

Cytogenetics: 7 87.94 cM

Gene Summary: Troponin T is the tropomyosin-binding subunit of troponin, the thin filament regulatory complex which confers calcium-sensitivity to striated muscle actomyosin ATPase activity. [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG222198