

## Product datasheet for **MG222316**

### Gpx4 (NM\_001037741) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Gpx4 (NM\_001037741) Mouse Tagged ORF Clone  
**Symbol:** Gpx4  
**Synonyms:** GPx-4; GSHPx-4; mtPHG; mtPHGPx; PHG; PHGPx; sn; snGPx  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG222316 representing NM\_001037741  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGGGCCGCGCGGCCGCCGCAAGCGGGGACGCTGCAGACAGCGCGCGGATCCCCGAGAGGCCGGCGAC  
 GCCGTGGACCTGGACGCCAAAGTCCTAGGAAACGCCCGGGCCCTCGGCGAAGGAAAGCGCGCGCGGCCG  
 CCGCAGGAGGGCGCGCCCTCGCCGGATGGAGCCATTCTGAACCTTTCAACCCGGGGCCTCTGTGCAA  
 GAGCCTCCCACTACTGCAACAGCTCCGAGTTCCTGGGCTTGTGTGCATCCCGGATGATTGGCGCTGTG  
 CGCGCTCCATGCACGAATCTCAGCCAAGGACATCGACGGGCACATGGTCTGCCTGGATAAGTACAGGGG  
 TTTCTGTGCATCGTACCAACGTGGCCTCGCAATGAGGCAAACTGACGTAAACTACACTCAGCTAGTC  
 GATCTGCATGCCGATATGCTGAGTGTGGTTTACGAATCCTGGCCTTCCCTGCAACCAAGTTTGGGAGGC  
 AGGAGCCAGGAAGTAATCAAGAAATCAAGGAGTTTGCAGCCGGCTACAACGTCAAGTTTGCATGTACAG  
 CAAGATCTGTGTAATGGGGACGATGCCACCCACTGTGGAAATGGATGAAAGTCCAGCCCAAGGGCAGG  
 GGCATGCTGGGAAATGCCATCAAATGGAACCTTACCAAGTTTCTATTGATAAGAACGGCTCGTGGTGA  
 AGCGCTATGGTCCATGGAGGAGCCCCAGGTGATAGAGAAGGACCTGCCGTGCTATCTC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG222316 representing NM\_001037741  
Red=Cloning site Green=Tags(s)

MGRAAARKRGRCRQRGGSPRRRRRGPGRQSPRKRPGPRRRRKARARRRRRARRRMEPIPEFPNPGPLLQ  
 EPPQYCNSSSEFLGLCASRDDWRCARSMHEFSAKDIDGHMVCLDKYRGFVCIVTNVASQ\*GKTDVNYTQLV  
 DLHARYAECGLRILAFPCNQFGRQEPGSNQEIKEFAAGYNVVKFDMYSKICVNGDDAHLWKWMKVQPKGR  
 GMLGNAIKWNFTKFLIDKNGCVVKRYGPMEEPQVIEKDLPCYL


TRTRPLE - GFP Tag - V



[View online »](#)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_001037741

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#) The expression of this clone is not guaranteed due to the nature of selenoproteins.

**OTI Annotation:** This clone encodes a selenoprotein containing the rare amino acid selenocysteine (Sec). Sec is encoded by UGA codon, which normally signals translational termination. Expression of this clone is not guaranteed due to the nature of selenoproteins.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001037741.4](#)

**RefSeq Size:** 1012 bp

**RefSeq ORF:** 762 bp


**Locus ID:** 625249

UniProt ID: [O70325](#)

Cytogenetics: 10 39.72 cM

**Gene Summary:** The protein encoded by this gene belongs to the glutathione peroxidase family, members of which catalyze the reduction of hydrogen peroxide, organic hydroperoxides and lipid hydroperoxides, and thereby protect cells against oxidative damage. Several isozymes of this gene family exist in vertebrates, which vary in cellular location and substrate specificity. This isozyme has a high preference for lipid hydroperoxides and protects cells against membrane lipid peroxidation and cell death. It is also required for normal sperm development; thus, it has been identified as a 'moonlighting' protein because of its ability to serve dual functions as a peroxidase, as well as a structural protein in mature spermatozoa. Disruption of this gene in mouse spermatocytes is associated with male infertility. This isozyme is also a selenoprotein, containing the rare amino acid selenocysteine (Sec) at its active site. Sec is encoded by the UGA codon, which normally signals translation termination. The 3' UTRs of selenoprotein mRNAs contain a conserved stem-loop structure, designated the Sec insertion sequence (SECIS) element, that is necessary for the recognition of UGA as a Sec codon, rather than as a stop signal. Transcript variants resulting from alternative splicing or use of alternate promoters have been described to encode isoforms with different subcellular localization. Pseudogenes of this locus have been identified on chromosomes 10 and 17. [provided by RefSeq, Jan 2019]

### Product images:



Circular map for MG222316