

Product datasheet for **MG222459**

Arhgef37 (NM_177828) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgef37 (NM_177828) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Arhgef37
Synonyms:	4933429F08Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG222459 representing NM_177828
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCTGACTTTGGAGCAGATGAGGCCTCTTCCAAGTCAGAGAGTCCTGAGCAGGAAGATCAGGGCTCTG
AGGACAGATCGCTGCTCCATCAGAGGCTGGCCATCAGAGAGCTCATTGACACAGAAGTCTCTACTTACA
CATGCTCCGACTCTGTGCCCTGTACATCAGAAGCCACCTACAACAGCTGCCACCAGGCGATCTGGACATC
CTGTTCTCCAACATCGATGACATCATCCAAGTGAGCAGCAGGTTCCCTCGAGGGTGCCAGGAGACAGCCT
GCAGGGAAGAGGAGCAAGCGCACCTGATTGGTAACTTATTCTTGGAAATCCAAGAGGAGTTTGAACAAGT
CTATAAGTCTACTGTGCCAACTATGACCAGGCCCTGCTGCTGGTCAAGGCGTACCAGAAGGAGCCAGAA
CTGCAGAAGGAGATCCAGGGCATCATCGAAGCTGCCATGCCACAAGCTGGACCTTCGGGTCTCAGTTTCT
TCCTTGTGATCCCTCTACAGAGGATTACCAAGTACCCCTCCTGCTGCAGAAAATCTGGAGAACACACC
TGCAGATGCCAGCGCCCACCCTGCCCTCAGAGAGCCACCTCTGCCCTCCAGGATGTGAACGGCAATATC
AACGAATACAAGATGCGCAAGGAAGTGGCTTAAAGTATACCAAAGTGGAGCAGCTGAGCCTTCGGGAGC
GGCTGGCCCGCATCAACACACATACTCTTTCGAAGAAGACCACCAGACTGAGCCAGATGTGAAGCAGGA
GGCGGGGCTTGTATCCCGGACAGAAGACAAGGAATTTGATGACCTGGAGGAGAGATTCCAGTGGGTGCT
CTGTGTGTGAAGGAGCTCAGAAACAATGTGGTCGATTTTATGGTAATCTGGAGGTTTTCTCAGCTTCC
GGCCACACGAACGGAGCCTGGATATCCCTGGGGACCTGCGGAGCAGTACTGCAGCCTAGCGAGGGACCT
TCAGCTCCATGCCTTCTGCAGTTTAAAGCAGCGGTTAAGAAGCCTGGTGTGGCAACCATTGTGCAGTCTG
GCCAGAGCCCTGGCTGGTCCCTCAGAACTTGATCAAGAAGCGCCTGGACAACTCCTGGACTTCGATCGAG
TAGAGGAGAAGCTGTTAGATGTGGCAGCGTGACCTGGGAGGAGGAGGCAGCCCGCAGCAATACCAGGC
GCTCAACTCCTTACTAGTGGCCGAGCTCCACAGTTTAAACCAACTCGTCATGCAGTGGCTGGGTGAGATC
CTGTGCACATTTGTGGTCTCCAGAGAGACCTCGCCGATCAAGTGTACGAAGGGCAGAGAGCACCATGG
CTCTGCTGCCTCATCACCACGTCTCAGAGCCGGACTTCCAGAAGCTGCTGGAGCACACACTAGGCCAGAG
CAGCAGTCAGCTCCGCTCTTTCAGGAGAGCTTGGAGAAGGTGCTTCCACCCTCCACCCTCAGCCACTC
CTTCCAGGGTCTGAACACCAGATGCAGGCTCTCCTTACCAGATATGGCCCTGAGAAGCTATACCAGGTGT
CAAGCAACATCAACGGGACTGGGACTCTGGACCTGACCCTGCTGCGGGGCCAAATCGTGGCCCTTCTACA
AAACAAGGATACCAAAGGCAACAACAGTCGCTGGCTGGTAGACACAGGGGACACCGGGGCTATGTGCA
GCTGGAAAATTCAGTTATACCACCCATCATTCCCAGTGAAGGAGCTCAGAGGCCAGACAGGGACAC
ATGAAGACTCCTGGCTTCCAGCCCCAGAGCCCACCCAGCTCTCTGTTCCCACTGTCCAACCATGAGTCA
GGTGGTGGCAGTATACCTTTTGTGGCCAGAAGCACCCACGAAGTGAAGTCTACAGGCAGGCCAGCCTGTG
ACCATCTAGAGGCCAAGACAAGAAGGGGAACCCAGAATGGAGCCTGGTGGAAAGCAAACGGACAGAGGG
GCTATGTGCCTTCCAACCTTCTAGCCAGGACCCCAAGCCCAACTCCCTGGGGCTGGAATCTTCCCTCT

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG222459 representing NM_177828
 Red=Cloning site Green=Tags(s)


MADFGADEASSKSESPEDQGSSEDRSLHQLAIRELIDTEVSYLHMLRLCASDIRSHLQQLPPGDLDI
 LFSNIDDIQVSSRFLRGLQETACREEQAHLIGNLFLEFQEEFEQVYKVCANYDQALLLVKAYQKEPE
 LQKEIQGIEEAAMPQAGPSGLSFFLVIPLQRITKYPLLLQKILENTPADASAHPALQRATSALQDVNGNI
 NEYKMRKEVALKYTKVEQLSLRERLARINTHTLSKKTTRL SQMLKQEAGLVSRTEDEKFDLEERFQWVS
 LCVKELRNNVDFMGNLEVFLSFRPHERSLDIPGGPAEQYCSLARDLQLHAF LQFKQRLRSLVWQPLCSL
 ARALAGPQNL IKKRLDKLLDFDRVEEKLLDVGSVTWEEEAARHEYQALNSLLVAELPQFNQLVMQWLGQI
 LCTFVVLQRDLADQVLRRAESTMALLPHHVSEPDFQKLEHTLGQSSQLRFLQESFEKVLPPSTTQPL
 LPGSEHQMQALLTRYGPEKLYQVSSNINGTGTLDLTLRQIVALLQNKDTKGNNSRWLVD TGGRGYVP
 AGKLQLYHPIIPSEKELRGQTGTHEDSWLPAPEPTQLSVPTVPTMSQVVAVYPPVARSTHEL SLQAGQPV
 TIL EAQDKKGNPEWSLVEANGQRGYVPSNFLARTPSPTPWGWNLPS


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_177828

ORF Size: 2028 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_177828.4](#), [NP_808496.2](#)

RefSeq Size: 3234 bp

RefSeq ORF: 2031 bp

Locus ID: 328967

UniProt ID: [A1IGU4](#)

Cytogenetics: 18 E1

Gene Summary: May act as a guanine nucleotide exchange factor (GEF).[UniProtKB/Swiss-Prot Function]