

Product datasheet for **MG222529**

Klrk1 (NM_001083322) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Klrk1 (NM_001083322) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Klrk1
Synonyms:	D6H12S2489E; NKG2-D; Nkg2d
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG222529 representing NM_001083322 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCAAATGCCATAATTACGACCTCAAGCCAGCAAAGTGGGATACTTCTCAAGAACACAGAAACAAA
GATTAGCACTAACTACCAGTCAACCTGGAGAAAATGGTATCATAAGAGGAAGATACCCTATAGAAAACT
CAAAATATCTCCAATGTTTCGTTGTTTCGAGTCCTTGCTATAGCCTTGGCAATTCGATTCACCCTTAACACA
TTGATGTGGCTTGCCATTTTCAAAGAGACGTTTCAGCCAGTATTGTGCAACAAGGAAGTCCCAGTTTCCT
CAAGAGAGGGCTACTGTGGCCCATGCCCTAACAACTGGATATGTCACAGAAACAAGTGTACCAATTTTT
TAATGAAGAGAAAACCTGGAACAGAGCCAAGCTTCTGTTTGTCTCAAATTCAGCCTTCTGAAGATA
TACAGTAAAGAAGAACAGGATTTCTTAAAGCTGGTTAAGTCTATCACTGGATGGGACTGGTCCAGATCC
CAGCAAATGGCTCCTGGCAGTGGGAAGATGGCTCCTCTCTCATACAATCAGTTAACTCTGGTGGAAAT
ACCAAAAGGATCCTGTGCTGTCTATGGCTCAAGCTTAAAGGCTTACACAGAAGACTGTGCAAATCTAAAC
ACGTACATCTGCATGAAAAGGGCGGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG222529 representing NM_001083322
Red=Cloning site Green=Tags(s)

MSKCHNYDLKPAKWDTSQEQKQRLALTTSPGENGIIRGRYPKIEKLKISPMFVVRVLAIALAIRFTLNT
 LMWLAIFKETFPVLCNKEVPVSSREGYCGPCPNWICHRNNCYQFFNEEKTWNQSQASCLSQNSSLLKI
 YSKEEQDFLKLKLSYHWMGLVQIPANGSWQWEDGSSLSYNQLTLVEIPKGSCAVYGGSSFKAYTEDCANLN
 TYICMKRAV


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001083322

ORF Size: 657 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001083322.2
RefSeq Size:	3158 bp
RefSeq ORF:	660 bp
Locus ID:	27007
UniProt ID:	O54709
Cytogenetics:	6 63.44 cM
Gene Summary:	Function as an activating and costimulatory receptor involved in immunosurveillance upon binding to various cellular stress-inducible ligands displayed at the surface of autologous tumor cells and virus-infected cells. Provides both stimulatory and costimulatory innate immune responses on activated killer (NK) cells, leading to cytotoxic activity. Acts as a costimulatory receptor for T-cell receptor (TCR) in CD8(+) T-cell-mediated adaptive immune responses by amplifying T-cell activation. Stimulates perforin-mediated elimination of ligand-expressing tumor cells. Signaling involves calcium influx, culminating in the expression of TNF-alpha. Participates in NK cell-mediated bone marrow graft rejection. May play a regulatory role in differentiation and survival of NK cells. Binds to ligands belonging to various subfamilies of MHC class I-related glycoproteins including RAET1A, RAET1B, RAET1C, RAET1D, RAET1E, H60 and MULT1.[UniProtKB/Swiss-Prot Function]