

Product datasheet for **MG222609**

Hap1 (NM_010404) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hap1 (NM_010404) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Hap1
Synonyms:	HAP; HAP-1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG222609 representing NM_010404
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGCGCCCCGAAGAGCAGGTGCAGAGCGGTGCGGGAGACGGGACAGGGTCGGGGACCCAGCAGCAGGCA
 CCCCCACGACCCAGCCTGCAGTTGGTCCCCTCCGGAGCCCTCGGCGGAGCCCAAACCTGCTCCAGCGCA
 GGAACCGGGTCCGGACAAAAATCAGGATCCCAGAACAGACAGGAAGCTTTTGTGCGTCCATGATCATT
 GGTGATTCGGACGCACCATGGACCCGCTACGTATTCCAGGGGCCTTACGGTCCCCGGGCCACTGGCCTGG
 GCACTGGAAGGCGGAGGAATCTGGAAGACACCAGCCGCTACATCGGCCGGAGGCCCGCGTGTCCGG
 CCCTGAGCGTGCGGCTTTATTCGAGAGCTGCAGGAAGCGTTGTGTCTAATCCACCACCCACGAAGAAG
 ATCACCAGATGATGTCAAAGTGTGTTGATTTGCTGGAAGAGAAAGAACGGGACCTGAACACAGCCG
 CCCGATCGCCAGTCCCTGGTAAACAGAACAGTGTCTTGATGGAGGAGAATAATAAGCTGGAACCAT
 GCTGGGCTCAGCCAGGGAGGAGATTTTACATCTCCGGAAGCAGGTGAACCTGCGAGATGACCTTCTTCAG
 CTCTACTCAGACTCTGATGACGATGATGATGAGGAAGACGAGGAAGACGAGGAAGAGGGCGAAGAGGAGG
 AACGAGAAGGACAGAGGGATCAAGACCAGCAGCACGACCACCCCTATGGTGCCCCAAGCCACACCTAA
 GGCTGAGACAGCGCACCGCTGCCACAGCTGGAACCCCTGCAGCAGAAGCTCAGGCTTCTGGAGGAAGAG
 AACGACCACCTGCGAGAGGAGGCTCCACCTTGACAACCTGGAGGACGAAGAGCAGATGCTCATTCTGG
 AATGTGTGGAGCAGTTCTCTGAAGCCAGCCAGCAGATGGCAGAGCTATCGGAAGTGTGTTGAGGCT
 GGAAGGCTATGAGAGGCAGCAGAAAGAGATCACTCAGTGCAGGCCGAGATCACCAGCTACAACAGCGT
 TGTCAGTCTTATGGGGCCAGACAGAGAACTGCAGCAGATGCTGGCCTCAGAGAAGGGGATCCACTCGG
 AGAGCCTGCGAGCTGGCTCTACATGCAGGATTATGGGAGCAGGCTCGTGACCGCCAGGAGGATGGGAA
 GAGTCATCGCCAGCGCTCCTCCATGCCCGCAGGCTCTGTCAACCACTATGGATACAGTGTGCTCTGGAT
 GCACCTTCAAGTTTCCAGAGACTGGCTGAGGAGCTCCGAACATCTCTGAGGAAGTTCATCACTGACC
 CTGCGTATTTATGGAGAGCGTGACACTCACTGCAGGGAGGGCGGAAGAAGGAGCAGAGGGCGATGCC
 ACCCCACCGCTCAAGATCTCAAGCCACCTGAAGATTTTGGGCTCCAGAGGAGTTGGTTCCTGAGGAG
 GAGCTGGGGGCCATCGAAGAGGTGGGACAGCTGAGGATGGGCAGGCAGAGAGAACGAGCAGGCGTCTG
 AGGAGACCGAGGCTGGGAGGAGGTGGAACCGGAGGTGGACGAGACCACAAGGATGAATGTGGTGGTCTC
 CGCCTGAGGCCAGTGGCTGGGCCCTCCACCTGGACATGAAGTATGCTCTCCAGCAACTGTCCAAC
 TGGCAGGACGCCACTCTAAGCGGCAGCAGAAGCAGAAGGTGGTCCCAAAGGTGAGTGTCCCGCAGAG
 GACACCCTCTGCCAGTGGGACAAGCTTTCGATCATCAACCATA

ACGCGTACGCGGCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG222609 representing NM_010404
 Red=Cloning site Green=Tags(s)


MRPKEQVQSGAGDGTGSGDPAAGTPTTQPAVGPAPEPSAEPKPAPAQGTGSGQKSGSRTKTGSFCRSMII
 GDSDAPWTRYVFGPYGPRATGLTGKAEGIWKTPAAYIGRRPGVSGPERAAFIRELQALCPNPPPTKK
 ITEDDVKVMLYLLEEKERDLNTAARIGQSLVKQNSVLMEENNKLETMLGSAREEILHLRKQVNLRDDLLQ
 LYSDSDDDDE
 NDHLREEASHLDNLEDEEQMLILECVEQFSEASQMAELSEVLVLRLEGYERQQKEITQLQAEITKLQQR
 CQSYGAQTEKLQQLASEKGIHSESLRAGSYMQDYGSRRPRDRQEDGKSHRQRSSMPAGSVTHYGYSVPLD
 ALPSFPETLAEELRTSLRKFITDPAYFMERRDTHCREGRKKEQRAMPPPPAQDLKPPEDFEAPEELVPEE
 ELGAIEEVGTAEDGQAEENEQASEETEAWEEVEPEVDETRMNVVVSALASGLGPSHLDKMYVLQQLSN
 WQDAHSKRQKQKQVVPKGECSRRGHPPASGTSFRSSTI

TRTRPLE – GFP Tag – V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_010404

ORF Size: 1794 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010404.3 , NP_034534.1
RefSeq Size:	3774 bp
RefSeq ORF:	1797 bp
Locus ID:	15114
UniProt ID:	O35668
Cytogenetics:	11 D
Gene Summary:	The protein encoded by this gene was first identified as a neuronal protein that binds the HD protein huntingtin. The protein also interacts with kinesin light chain, 14-3-3 proteins, and Abelson helper integration site 1 protein. The protein is involved in intracellular trafficking of vesicles and organelles, and lack of the protein results in neuronal death resembling the hypothalamic degeneration that occurs in Huntington's disease. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2009]