

## Product datasheet for **MG222737**

### **Cadm1 (NM\_001025600) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Cadm1 (NM_001025600) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cadm1
Synonyms:	2900073G06Rik; 3100001I08Rik; AI987920; BI2; Igsf4; Igsf4a; Necl2; RA175; RA175A; RA175B; RA175C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>MG222737 representing NM\_001025600  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGGCGAGTGCTGTGCTGCCGAGCGGATCCAGTGTGCGCGGCAGCGGCTGTGGCGGCGCGCGCGCC  
 CTCCAGGGCTCCGGCTCCGGCTCCTGCTGTTGCTCCTTTTCGGCCGCGCACTGATCCCCACAGGTGATGG  
 ACAGAATCTGTTTACTAAAGACGTGACAGTGATTGAAGGAGAAGTGGAACCATCAGCTGCCAGGTCAAT  
 AAGAGTGACGACTCAGTGATCCAGCTCCTGAACCCCAACAGGCAGACATTTACTTCAGGGACTTCAGGC  
 CTTTGAAGGACAGCAGGTTTCAGCTGCTGAATTTTCTAGCAGTGAACCTCAAAGTGTCACTGACGAATGT  
 CTCAATCTCGGATGAAGGGAGATACTTCTGCCAGCTCTACACGGACCCCCACAGGAGAGTTACACCACC  
 ATCACAGTCTGGTTCTCCACGTAACCTGATGATCGATATCCAGAAAGACACGGCAGTTGAAGGGGAGG  
 AGATTGAAGTCAACTGTACTGCCATGGCCAGCAAGCCAGCGACGACCATCAGGTGGTTCAAAGGGAACAA  
 GGAACCTAAAGGCAAATCAGAGGTGGAGGAGTGGTCCGACATGTACACTGTGACCAGTCAGCTGATGCTG  
 AAGGTGACAAGGAGGACGACGGGTCCCGGTGATCTGCCAGGTGGAGCACCCCTGCGGTCAGTGGAAACC  
 TGCAGACCCAGCGCTATCTAGAAGTGCAGTATAAACCGCAAGTGCATATCCAGATGACTTACCCTCTGCA  
 AGGCCTAACCCGGGAAGGGGATGCATTTGAGTTAACGTGTGAAGCCATCGGGAAGCCCCAGCCTGTGATG  
 GTAACCTGGGTGAGAGTCGATGATGAAATGCCTCAACATGCCGTACTGTCTGGGCCAAACCTGTTTCATCA  
 ATAACCTAAACAAAACAGATAACGGTACTTACCCTGTGAGGCTTCCAACATAGTGGGAAAGGCTCATTC  
 GGACTATATGCTGATGTATACGATTCTCGAGCAGGTGAAGAGGGGACCATTGGGCGAGTGGACCACGCA  
 GTGATTGGTGGCGTCGTAGCCGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGTGGT  
 TTGCCAGACATAAAGGTACATACTTCACTCATGAAGCCAAGGAGCCGATGACGCAGCAGACGCAGACAC  
 AGCTATAATCAATGCAGAAGGAGGACAGAACAACCTCCGAAGAAAAGAAAGAGTACTTCATC

**ACGCGTACGCGGCCGCTCGAG** – GFP Tag – GTTTAA

**Protein Sequence:**

>MG222737 representing NM\_001025600  
 Red=Cloning site Green=Tags(s)

MASAVLPSGSQCAAAAIVAAAAAPPGLRLRLLLLLLAAAALIPTGDGQNLFTKDVTVIEGEVATISCQVN  
 KSDDSVIQLLNPNRQTIYFRDFRPLKDSRFQLLNFSSELKVSLTNVSI SDEGRYFCQLYDPPQESYTT  
 ITVLVPPRNLMIQKDTAVEGEEIEVNCTAMASKPATTIRWFKGNKELKKGKSEVEEWSMYTVTSQLML  
 KVHKEDDGVPVICQVEHPAVTGNLQTRYLEVQYKQVHIQMTYPLQGLTREGDAFELTCEAIGKPQPM  
 VTWVRVDEMPQHAVLSGPNLFINNLNKTDNGTYRCEASNIVGKAHSDYMLVYVDSRAGEEGTIGAVDHA  
 VIGGVVAVVVFAMLCLLIILGRYFARHKGTYF THEAKGADDAADADTAIINAEGGQNNSEEKKEYFI

**TRTRPLE** – GFP Tag – V

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**


**ACCN:** NM\_001025600

**ORF Size:** 1251 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001025600.1](#), [NP\\_001020771.1](#)

**RefSeq Size:** 4270 bp

**RefSeq ORF:** 1254 bp


**Locus ID:** 54725

**UniProt ID:** [Q8R5M8](#)

**Cytogenetics:** 9 A5.3

**Gene Summary:**

Mediates homophilic cell-cell adhesion in a Ca<sup>2+</sup>-independent manner. Also mediates heterophilic cell-cell adhesion with CADM3 and NECTIN3 in a Ca<sup>2+</sup>-independent manner. Acts as a tumor suppressor in non-small-cell lung cancer (NSCLC) cells. Interaction with CRTAM promotes natural killer (NK) cell cytotoxicity and interferon-gamma (IFN-gamma) secretion by CD8<sup>+</sup> cells in vitro as well as NK cell-mediated rejection of tumors expressing CADM3 in vivo. May contribute to the less invasive phenotypes of lepidic growth tumor cells. In mast cells, may mediate attachment to and promote communication with nerves. CADM1, together with MITF, is essential for development and survival of mast cells in vivo. Acts as a synaptic cell adhesion molecule and plays a role in the formation of dendritic spines and in synapse assembly. May be involved in neuronal migration, axon growth, pathfinding, and fasciculation on the axons of differentiating neurons. May play diverse roles in the spermatogenesis including in the adhesion of spermatocytes and spermatids to Sertoli cells and for their normal differentiation into mature spermatozoa.[UniProtKB/Swiss-Prot Function]

**Product images:**

Circular map for MG222737