

## Product datasheet for **MG222955**

### Ebag9 (NM\_019480) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ebag9 (NM_019480) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ebag9
Synonyms:	A1835379; Rcas1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG222955 representing NM_019480 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCATCACTCAGTTTCGGTTATTTAAAGTTTGTACCTGCCTAGCAACAGTGTTCGTTCTAAAGA  
GATTAATATGCAGATCTGGCAGAGGAAGGAAATTAAGTGGAGATCAAATAACTTTGCCACTACAGTTGA  
TTATTCGTCTGCCCTAAGCAGACAGATGTGGAGGAGTGGACTTCCTGGGATGAAGATGCACCCACAAGT  
GTAAGATCGAAGGAGGGAATGGGAATGTGGCTACTCAGCAGAACTCTTTGGAACAACCTGGAACCTGACT  
ACTTCAAGGACATGACACCAACTATACGAAAACTCAGAAAATTGTCATTAAGAAGAGAGAACCATTGAA  
TTTTGGTGTCCCAGATGGTAGCACGGGTTTTCCAGTAGATTGGCAGCTACACAAGACATGCCTTTTATT  
CATCAATCTTCAGAATTAGGTGACTTAGATACATGGCAAGAAAATAGCAATGCGTGGGAAGAGGAAGAAG  
ATGCAGCCTGGCAAGCAGAAGAGGTGCTGAGGCAGCAGAAAATTGCAGACAGAGAGAAGAGAGCAGCAGA  
ACAACAGAGGAAGAAAATGGAGAAGGAAGCGCAGCGCCTGATGAAGAAGGAGCAGAACAAAATCGGTGTG  
AAGCTTTCC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG222955 representing NM\_019480  
 Red=Cloning site Green=Tags(s)

MAITQFRLFKVCTCLATVFSFLKRLICRSGRGRKLSGDQITLPTTVDYSSVPKQTDVEEWTSWDEDAPTS  
 VKIEGGNGNVATQQNSLEQLPDPYFKDMPTPIRKTQKIVIKKREPLNFGVPDGGSTGFSSRLAATQDMPFI  
 HQSSELGDLDTWQENSNAWEEEEAAWQAEVLRQOKIADREKRAAEQQRKKMEKAQRLMKKEQNKIGV  
 KLS


TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_019480

**ORF Size:** 639 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019480.5</a>
<b>RefSeq Size:</b>	1606 bp
<b>RefSeq ORF:</b>	642 bp
<b>Locus ID:</b>	55960
<b>UniProt ID:</b>	<a href="#">Q9D0V7</a>
<b>Cytogenetics:</b>	15 B3.2
<b>Gene Summary:</b>	May participate in suppression of cell proliferation and induces apoptotic cell death through activation of interleukin-1-beta converting enzyme (ICE)-like proteases.[UniProtKB/Swiss-Prot Function]