

Product datasheet for **MG223023**

Braf (NM_139294) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Braf (NM_139294) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Braf
Synonyms:	9930012E13Rik; AA120551; AA387315; AA473386; B-raf; Braf-2; Braf2; C87398; C230098H17
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG223023 representing NM_139294
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGCGCTGAGTGGCGGCGGTGGCAGCAGCAGCGGTGGCGGCGGCGGCGGTGGCGGCGGCGGTGGCG
 GTGGCGACGGCGGCGGCGCGCCGAGCAGGGCCAGGCTCTGTTCAATGGCGACATGGAGCCGGAGGCCGG
 CGCTGGCGCCGCGGCTCTTCGGCTGCGGACCCGGCCATTCTGAAGAGGTATGGAATATCAAGCAAATG
 ATTAAGTTGACACAGGAACATATAGAGGCCCTATTGGACAAATTTGGTGGAGAGCATAACCCACCATCAA
 TATACCTGGAGGCTATGAAGAGTACACCAGCAAGCTAGATGCCCTTCAGCAAAGAGAACAGCAGCTTTT
 GGAATCCCTGGTTTTTCAAACCTCCACAGATGCATCACGGAACAACCCCAAGTACCACAGAAACCTATC
 GTTAGAGTCTTCCTGCCAACAAACAGAGGACAGTGGTACCCGCAAGATGTGGTGTACAGTTTCGAGACA
 GTCTAAAGAAAGCACTGATGATGAGAGGTCTCATCCAGAATGCTGTGCTGTTTACAGAATTCAGGATGG
 AGAGAAGAAACCAATTGGCTGGGACACGGACATTTCTGGCTTACTGGAGAGGAGTTACATGTTGAAGTA
 CTGGAGAATGTCCCACTTACAACACACAACCTTTGTACGAAAACTTTTTTACCTTAGCATTGTTGTGACT
 TTTGCCGAAAGCTGCTTTTCCAGGGTTCCGTTGTCAAACATGTGGTTATAAATTTACCAGCGTTGTAG
 TACAGAGGTTCCACTGATGTGTGTAATATGACCAACTTGATTTGCTGTTTGTCTCCAAGTTCTTTGAG
 CATCACCCAGTACCACAGGAGGAGGCTCCTTCCAGAGACTGCCCTTCCATCTGGATCCTCTCCGCAC
 CCCCTCAGACTCTACTGGGCCCAAATCTCCACAGTCCATCTCCTTCAAAATCCATTCCAATTCACA
 GCCCTTCCGACCAGCAGATGAAGATCATCGCAATCAGTTTGGGCAACGAGACCGGCTCCTCCTCAGCTCCC
 AAAGGGATTTGATTAGAGACCAGGGTTTCGTGGTGATGGAGCCCTTGAACCAACTGATGCGCTGTCT
 TCGGAAATACCAATCCCGGACTCCAGCCCTCCTCCATTCTGTCCCCAGTGAAATAGTGTGTTTTT
 GAGCCTGGCCAGTGTTCAAGGGTCAACCACAGGCTTGTCGCCACCCCGCTGCCTCATTACCTGGCT
 CACTACTAAGCTGAAAGCCTTACAGAAATCTCCAGGTCCTCAGCGGAAAGGAAGTATCTTCTCCTC
 ATCCTCGGAGGACAGAAGTCGGATGAAAACACTTGGTAGAAGAGATTCAAGTGTACTGGGAGATTCT
 GATGGACAGATTACAGTGGGACAGAGAATTGGATCTGGGTCAATTTGAACTGTCTACAAGGGAAAGTGGC
 ATGGTGATGTGGCAGTAAAATGTTGAATGTGACAGCACCCACACCTCAACAGCTACAGGCCTTCAAAAA
 TGAAGTAGGAGTGCTCAGGAAACTCGACATGTGAATATCCTCCTTTTTCATGGGCTATTCTACAAGCCA
 CAACTGGCAATTGTTACACAGTGGTGTGAGGGCTCCAGCTTATATCACCATCTCCACATCATTGAGACCA
 AATTTGAGATGATCAAACCTTATAGATATTGCTCGGCAGACTGCACAGGGCATGGATTACTTACACGCCAA
 GTCAATCATCCACAGAGACCTCAAGAGTAATAATATTTCTTTCATGAAGACCTCAGGTAATAAATAGGT
 GACTTTGGTCTAGCCACAGTGAATCTCGGTGGAGTGGTCCCATCAGTTTGAACAGTTGTCTGGATCTA
 TTTTGTGGATGGCACCAGAAGTAATCAGAATGCAAGATAAAAACCCGTATAGCTTTCAGTACAGCTGTA
 TGCGTTTGGGATTGTTCTGTACGAACCTGATGACCGCCAGCTACCTTATCAAACATCAACAACAGGGAT
 CAGATAATTTTTATGGTGGGACGAGGATACCTATCTCCAGATCTCAGTAAGGTACGGAGTAACTGTCCAA
 AAGCCATGAAGAGATTAATGGCAGAGTGCCTCAAAAAGAAAAGAGACGAGAGACCACTCTTCCCCAAAT
 TCTCGCCTCCATTGAGCTGCTGGCCGCTCATTGCCAAAAATTCACCGCAGTGCAATCAGAACCTTCTCTG
 AATCGGGCTGGTTTCAAACAGAAGATTTTAGTCTGTATGCTTGTGCTTCTCGAAAAACCCCATCCAAG
 CAGGGGATATGGAGAATTTGCAGCCTTCAAG

ACCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG223023 representing NM_139294
 Red=Cloning site Green=Tags(s)


```

MAALSGGGGSSGGGGGGGGGGGGGGGGGGAEQGGQALFNGDMEPEAGAGAAAASSAADPAIPEEVWNIKQM
IKLTQEHEIALLDKFGGEHNPPSIYLEAYEYYSKLDALQQREQLLESVLFQPTDASRNNPKSPQKPI
VRVFLPNKQRTVVPARCGVTVRDSLKALMMRGLIPECCAVYRIQDGEKKPIGWDTDISWLTGEELHVEV
LENVPLTTHNFVRKTFFTLAFCDFCRKLFLQGFRCQTCGYKFHQRCSTEVPLMCVNYDQLDLLFVSKFFE
HHFVPQEEASFETALPSGSSSAPPDSTGPIITSPSPSKSIPQPFPADEDRNQFGQRDRSSSAP
NVHINTIEPVNIDEKFPVELQDQRDLIRDQGFQDGAFLNQLMRCLRKYQSRTPSPLLHVPSEIVDFD
EPGPVFRGTTGLSATPPASLPGSLTNVKALQKSPGQERKSSSSSSSEDRSRMKTLGRRDSSDDWEIP
DGQITVQQRIGSGSFGTYVYKGVKWHGDVAVKMLNVTAPTQQQLQAFKNEVGVLKTRHVNILLFMGYSTKP
QLAIVTQWCEGSSLYHHLHIETKFEMIKLIDIAHQTAQGM DYLHAKSIIHRDLKSNINIFLHEDLTVKIG
DFGLATVKSRWSGSHQFEQLSGSILWMAPEVIRMQDKNPYSFQSDVYAFGIVLYELMTGQLPYSINNRD
QIIFMVGGRGYSPLSKVRSNCPKAMKRLMAECLKKRDERPLFPQILASIELLARSPLKIHRASASEPSL
NRAGFQTEDFSLYACASPKTPIQAGGYGEFAAFK
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI

Plasmid Map:



ACCN: NM_139294

ORF Size: 2412 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in *E. coli* are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139294.5](#), [NP_647455.3](#)

RefSeq Size: 9728 bp

RefSeq ORF: 2415 bp

Locus ID: 109880

UniProt ID: [P28028](#)

Cytogenetics: 6 18.43 cM

Gene Summary: Involved in the transduction of mitogenic signals from the cell membrane to the nucleus. May play a role in the postsynaptic responses of hippocampal neuron (By similarity).
[UniProtKB/Swiss-Prot Function]