

Product datasheet for **MG223140**

Tmed8 (NM_001033475) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmed8 (NM_001033475) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Tmed8
Synonyms: 6430595O10Rik; AI447224; Gm1184; Mem1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG223140 representing NM_001033475
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGGACCGGCAGGCAGCTGAGGGGCCGCCTTCTGGAGCCCCGAGCTCGCCGAGGGTCGGCGGGCG
 GCGTCGGGGACCGACGGGGAGTGGAGGAGAGCCAGGCGGCCCTCAGAGAAGGAAGATCTGAAAGCAC
 CAACGTCAGCTCTCCACTGGCTTCTGCTTCTGATCCAGCAGCCGAGTCTCACCATACAGGCCACAGATG
 GTATCTCCAGCGAGTAAGGATACCAGGAAGATCTACAGAACGTGGCTGGTGCCTCAGAGGGTCAGGCTC
 CGGGGAGCAGGGCGCTGCCAGCAGGCCAGACCCAGGTTCTCAGTGAGATGGCCAAGTACCAAGCTCC
 ACAGCGACCTGAGGACACAGTTATGATCCAGTCTGAGCACACTGGGGCTATAGATGTTCTCCAGCGGAC
 TTGGAATCTGCAGATCTTCTGGGGACCACAGGAAAGTCTCTCCGCTCTGATGGCTCCTCCGTGCTCT
 GGACCTTTGCCAAAGTGAAGGAGTTCAAAGCAAGCTGGGCAAGGAGAAGAACAGCCGTCTGGTAGTGAA
 GCGGGGAGAAGTGGTCACCATCCGGTCCCTACCCACCCGAGGGAAAGCGCGTGTGCTGGGAATTTGCA
 ACTGATGACTATGATATTGGTTTGGAGTTATTTGACTGGACCCAGTAACCAGCACTGATATAACCG
 TGCAGGTCAGTGACTCCAGTGAGGATGAGGAGGAGGAAGAGGATGAAGAGGAGGAGATTGAAGAACCCTG
 CCCGGTGGGAGATGTAGAGAGAGGCTCCAGGAGTCTCTACGGGGCCGCTATGGGAGGTCATGCCTGTG
 TACCGGCGGACAGCCACCGAGATGTGCAAGCTGGCAGCCATGATTACCCTGGGAGGGCATCTACTTGC
 TCAAGTTTGATAACTCTACTCCTTGCTCCGCAACAAGACTCTCTACTTCCATATCTACTATACTAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG223140 representing NM_001033475
 Red=Cloning site Green=Tags(s)

```
MSDRQAAEGPAFWSPAARRGSAGGVGDRRGVEESQAAASEKEDLESTNVSSPLASASDPAAESSPYRPQM
VSPASKDTTDLQNVAGASEGQAPGEQAALPAGQTQVLSSEMAKYQAPQRPEDTVMIQSEHTGAIDLVPAD
LESADLLGDHRKVSPLMAPPVWTFKVKFKSKLKGKNSRLVVKRGEVVTIRVPTHPEGKRVCWEFA
TDDYDIGFGVYFDWTPVTSTDITVQVSDSSEDEEEEEDEEEEEIEEPVPVGDVERGSRSSLRGRYGEVMPV
YRRDShRDVQAGSHDYPGEGIYLLKFDNSYSLLRNKTLVFHIIYYS
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



ACCN: NM_001033475

ORF Size: 978 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033475.3](#), [NP_001028647.1](#)

RefSeq Size: 7179 bp


RefSeq ORF: 981 bp

Locus ID: 382620

UniProt ID: [Q3UHI4](#)

Cytogenetics: 12 D2

Product images:



Circular map for MG223140