

Product datasheet for **MG223185**

Cpeb3 (NM_198300) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cpeb3 (NM_198300) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Cpeb3
Synonyms:	4831444O18Rik; CPE-BP3; mKIAA0940
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG223185 representing NM_198300
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCAGGATGATTACTGATGGACAAAAGCAAAACCCAGCCCCAGTCTCAGCAGCAGCAGCGGCAGCAGC
 AGCAGCAGCAGCAACAGCTCCAGCCCGAGCCCGCGCAGCTGAAGCCCCGTCCACGCCCTCTCCTCAGA
 GATCCCCAAGCCGAAGACAGTAGCGCAGTGCCGGCCCTCAGCCCCGCCTCGGCTCCGCCAGCCCCAAC
 GGCCCCGACAAGATGCAGATGGAGTCGCCGCTCCTGCCGGGCTTGAGTTTCCATCAGCCCCCTCAGCAGC
 CGCCGCCCGCGCAGGAGCCACGGCGCCCGAGCGTCGCTGTGCCGCTCCTTCGGCAGCACCTGGTCCAC
 AGGCACTACAAACGCGGTGGAGGACAGCTTCTCCAGGGGATCACCCAGTCAACGGGACCATGCTCTTC
 CAAAATTTCCCGCACCACTCAACCCGGTCTTCGGAGGCACCTTCTCGCCGAGATCGGCTGGCGCAGA
 CCCAGCACCATCAGCAGCTCCACCGCCCGCGCCAGCCGCCGAGCCCGCGCAGCCCCCGCAGGCGCA
 GCCCTCGCAGCAACGCCGCTCGCCCGCCAGCCCGAGGCGCCCTATGCGCAGAGGAGCGCCGCTGCC
 TACGGCCACAGCCCATCATGACCAGCAAGCCGCTCTATCTCGGCGGTGCGGGCTGCTGCGGCCGCGG
 CCGCTGCCTCTTCGGCCTCGTCCAGCTGGAACACGCACAGAGCGTGAACGCCGCTGGAGTGCTCCGTC
 CAACCCGTGGGGTGGCCTGCAGGCGGGCCGAGACCCTCGCCGCGCGGTGCGCGTGGGCGTGGGTGTAGGT
 GTGGGGGTACCTCCCCGCTTAACCCATCTCGCCGCTCAAAAAGCCCTTCTCCAGCAATGTGATCGCTC
 CACCCAAGTTCCTCGTCAGCCCCGCTCACCTCCAAGTCTGGATGGAGGATAACGCTTTCCGG[AT]C
 CGATAATGGTAACAATCTGTTGCCTTTTCAGGACCGGAGTAGGCCCTATGACACTTTAACTTGACTCC
 TTGGAGAACTCCTTAATGGATATGATAAGGACTGACCATGAGCCTCTGAAAGGTAACACTACCTCCCA
 GTGGCCCAACCAATGAGTTTCGCTGATATAATGTGGAGGAATCATTTTGCAGGACGCATGGGAATAAATTT
 CCACCATCCAGGAACAGATAACATTATGGCACTTAACACCAGGAGCTATGGGCGGAGACGAGGTCCGGTCT
 TCTCTCTTTCCCTTTGAAGACGCTTCTTGATGACAGCCATGGTGATCAGGCCTTATCGTCTGGCTTAA
 GTTCTCTACTCGATGTCAAAATGGGGAACGAGTGGAACGCTACTCTAGAAAGGTGTTTGTGGAGGCT
 TCCTCTGACATCGATGAAGATGAGATCACTGCCAGCTTTCGAGGTTTGGACCCCTTGTTGGTGGACTGG
 CCTCATAAAGCCGAGAGCAAATCGTACTTCCCTCCTAAAGGCTATGCTTTCCTGCTGTTCCAGGAGGAGA
 GCTCAGTGCAAGCTCTGATAGATGCCTGCCTGGAGGAGGATGGGAACTGTACCTGTGTGTGTCAGCCC
 AACCATCAAGGATAAACCAGTGAATCCGACCGTGGAACCTAAGTGACAGTGACTTTGTAATGGATGGC
 TCTCAGCCTTTGGACCCAGAAAACTATCTTTGTTGGGGGAGTCCACGACCCCTTCGAGCTGTTGAAC
 TGGCAATGATAATGGACCGTTTGTACGGTGGTGTGTTGCTATGCTGGCATTGACACAGACCCAGAGCTGAA
 GTACCCGAAAGGTGCTGGGCGTGTGCGTTCTCCAATCAGCAGAGCTACATTGCAGCCATCAGTGCGCGC
 TTTGTGCAACTTCAACACAACGACATTGACAAACGGGTGGAAGTGAAGCCGTACGTGCTGGATGATCAGA
 TGTGTGACGAGTGTGAGGACACGCTGCGGGGGGAAGTTTGCCCATTTCTGTGCCAACGTCACCTG
 TCTGCAGTACTACTGTGAATACTGCTGGCCAGCATCCACTCCCGAGCCGGCCGTGAGTTCCACAAACCG
 CTGGTGAAGGAGGAGGCGACCGCCCCCGGCACGTCCCGTTCCGCTGGAGC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG223185 representing NM_198300
 Red=Cloning site Green=Tags(s)


MQDDLMDKSKTQPQSQQQQRQQQQQQQLQPEPGAAEAPSTPLSSEIPKPEDSSAVPALSPASAPPAPN
 GPDQMOMESPLLPLSFHQPPQPPPPQEPTAPGASLSPSFGSTWSTGTTNAVEDSFFQGITPVNGTMLF
 QNFPHHVNPVFGGTFSPQIGLAQTQHHQPPPPAPQPPQAQPPQAQPSQQRSPASPSQAPYAQRSAAA
 YGHQPIMTSKPSSSSAVAAAAAAAASSSSWNTHQSVNAAWSAPSNPWGGLQAGRDPRRAVGVGVG
 VGVPSPLNPI SPLKKPSSNVIAPPKFPRAAPLTSKSWMEDNAFRXXDNGNLLPFQDRSRPYDTFNLHS
 LENS LMDMIRTDHEPLKKGHYPPSGPPMSFADIMWRNHFAGRMGINFHHPGTDNIMALNTRSYGRRRGRS
 SLFPFEDAFLDDSHGDQALSSGLSSPTRCQNGERVERYSRKVFVGGLPDIDEDEITASFRFGLVVDW
 PHKAESKSYFPPKGYAFLLFQEESSVQALIDACLEEDGKLYLCVSSPTIKDKPVQIRPWNLSDSDFVMDG
 SQPLDPRKTI FVGGVPRPLRAVELAMIMDRL YGGVCYAGIDTDPELKYPKGAGRVAF SNQSYIAAISAR
 FVQLQHNDIDKRVEVKPYVLDDQMCDECQGTTCGGKFAPFFCANVTCLQYYCEYCWASIHSRAGREFHKP
 LVKEGGDRPRHVPFRWS

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2095_g09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:




ACCN: NM_198300

ORF Size: 2151 bp

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_198300.2
RefSeq Size:	5792 bp
RefSeq ORF:	2151 bp
Locus ID:	208922
UniProt ID:	Q7TN99
Cytogenetics:	19 C2

Gene Summary:

Sequence-specific RNA-binding protein which acts as a translational repressor in the basal unstimulated state but, following neuronal stimulation, acts as a translational activator (PubMed:17024188, PubMed:26074072). In contrast to CPEB1, does not bind to the cytoplasmic polyadenylation element (CPE), a uridine-rich sequence element within the mRNA 3' UTR, but binds to a U-rich loop within a stem-loop structure (PubMed:17024188). Required for the consolidation and maintenance of hippocampal-based long term memory (PubMed:26074003). In the basal state, binds to the mRNA 3' UTR of the glutamate receptors GRIA1 and GRIA2 and negatively regulates their translation (PubMed:17024188, PubMed:22153079). Also represses the translation of DLG4, GRIN1 GRIN2A and GRIN2B (PubMed:24155305). When activated, acts as a translational activator of GRIA1 and GRIA2 (PubMed:22153079, PubMed:26074003). In the basal state, suppresses SUMO2 translation but activates it following neuronal stimulation (PubMed:26074071). Binds to the 3' UTR of TRPV1 mRNA and represses TRPV1 translation which is required to maintain normal thermoception (PubMed:26915043). Binds actin mRNA, leading to actin translational repression in the basal state and to translational activation following neuronal stimulation (PubMed:26074072). Negatively regulates target mRNA levels by binding to TOB1 which recruits CNOT7/CAF1 to a ternary complex and this leads to target mRNA deadenylation and decay (By similarity). In addition to its role in translation, binds to and inhibits the transcriptional activation activity of STAT5B without affecting its dimerization or DNA-binding activity. This, in turn, represses transcription of the STAT5B target gene EGFR which has been shown to play a role in enhancing learning and memory performance (By similarity). In contrast to CPEB1, CPEB2 and CPEB4, not required for cell cycle progression (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG223185