

## Product datasheet for **MG223247**

### Ssbp1 (NM\_028358) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Ssbp1 (NM\_028358) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Ssbp1  
**Synonyms:** 2810480P10Rik; G630031O20Rik; mtDBP; MtSSB  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG223247 representing NM\_028358  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGTTTCGAAGACCTGTGTTACAGGATTTTCGTCAGTTTGTAAAGACATGAGTCTGAAGTAGCCAGCAGTT  
 TGGTTCCTGAACGATCTCTGAATCGTGTTCACTTGGACGAGTAGGTCAGGACCCTGTCATGAGACA  
 GGTGGAAGGAAAAACCCAGTCACAATATTTCTCTAGCAACAAATGAGATGTGGCGATCAGGGGATAGT  
 GAAGTATACCAAATGGGTGACGTTAGTCAGAAGACGACGTGGCACAGAATATCAGTGTTCGACCAGGCC  
 TCAGAGATGTGGCATATCAGTATGTGAAAAGGGGCTCGTATATTTGTGGAAGGAAAGTGGACTATGG  
 CGAGTACATGGATAAAAAACAATGTGAGGCGCAAGCAACAACAATCATAGCTGGTAAGAAGCTCGTGGT  
 CACAGTGTCTCTGGTTGTTTCATTGGAAGGTTGGCG

**ACGCGTACGCGGCCGCTCGAG** - GFP Tag - GTTTAA

**Protein Sequence:** >MG223247 representing NM\_028358  
 Red=Cloning site Green=Tags(s)

MFRRPVLQVFRQFVRHESEVASSLVLERSLNRVQLLGRVGDPMRQVEGKNPVTIFSLATNEMWRSRGS  
 EYVQMGDVSQKTTWHRISVFRPGLRDVAYQYVKKGARIFVEGKVDYGEYMDKNNVRRQATTTIAGKLLV  
 HSVSGCSLEGLA

**TRTRPLE** - GFP Tag - V


**Restriction Sites:** Sgfl-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM\_028358

ORF Size: 456 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_028358.2, NP_082634.1</u>
<b>RefSeq Size:</b>	1063 bp
<b>RefSeq ORF:</b>	459 bp
<b>Locus ID:</b>	381760
<b>UniProt ID:</b>	<u>Q9CYR0</u>
<b>Cytogenetics:</b>	6 B1
<b>Gene Summary:</b>	This protein binds preferentially and cooperatively to ss-DNA. Probably involved in mitochondrial DNA replication. Associates with mitochondrial DNA (By similarity). [UniProtKB/Swiss-Prot Function]