

## Product datasheet for **MG223278**

### Ranbp9 (NM\_019930) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ranbp9 (NM_019930) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Ranbp9
Synonyms:	IBAP-1; Ibp1; RanBPM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)



**Protein Sequence:** >MG223278 representing NM\_019930  
 Red=Cloning site Green=Tags(s)

```

MSGQPPPPPPQQQPPPPPPAAGAPGSGVALPAAAAASAGSSPAGSPGGGLGAAAAALLHPPPPPPPPPA
TAAPPPASAAAAPATAPPGLAVGPGPAAGVVPVGLAAGSSAAAPFPHGDSALNEQEKLQRRLKRLYPVAV
DEQETPLPRSWSPKDKFSYIGLSQNNLRVHYKHGKTPKDAASVRATHPIPAACGIYYFEVKIVSKGRDG
YMGIGLSAQGVNMNRLPGWDKHSYGYHGDDGHSFCSSGTGQPYGPTFTTGDVIGCCVNLINNTCFYTKNG
HSLGIAFTDLPPNLYPTVGLQTPGEVVDANFGQHPFVFDIEDYMRERWTKIQAQIDRFPIGDREGEWQTM
IQKMYSSYLVHHGYCATAEAFARSTDQTVLEELASIKNRQRIQKLVLAGRMGEAIEETTQQLYPSLLERNP
NLLFTLKVRFIEMVNGTDSEVRCLGGRSPKSQDSYPVSPRPFSSPSMSPSHGMSIHS LAPKSSTAHFS
GFESCNGVINSKAHQSYCHSKHQLSSLTYPENLSLNVSRSQVNNFTSNDVDMETDHYSNVGETSSNG
FLNGSSKHDHEMEDCDTEMEVDCSQLRRQLCGGSQAAIERMIHFGRELQAMSEQLRRECGKNTANKKMLK
DAFSLLAYSDPWNVSPVGNQLDPIQREPVCALNSAILETHNLPKQPPLALAMQATQCLGLMARSVGVGSC
AFATVEDYLH
  
```

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**ACCN:** NM\_019930

**ORF Size:** 2130 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019930.2](#), [NP\\_064314.2](#)

**RefSeq Size:** 3064 bp

**RefSeq ORF:** 2133 bp


**Locus ID:** 56705

**UniProt ID:** [P69566](#)

**Cytogenetics:** 13 A4

**Gene Summary:** May act as scaffolding protein, and as adapter protein to couple membrane receptors to intracellular signaling pathways. Acts as a mediator of cell spreading and actin cytoskeleton rearrangement. Core component of the CTLH E3 ubiquitin-protein ligase complex that selectively accepts ubiquitin from UBE2H and mediates ubiquitination and subsequent proteasomal degradation of the transcription factor HBP1. May be involved in signaling of ITGB2/LFA-1 and other integrins. Enhances HGF-MET signaling by recruiting Sos and activating the Ras pathway. Enhances dihydrotestosterone-induced transactivation activity of AR, as well as dexamethasone-induced transactivation activity of NR3C1, but not affect estrogen-induced transactivation. Stabilizes TP73 isoform Alpha, probably by inhibiting its ubiquitination, and increases its proapoptotic activity. Inhibits the kinase activity of DYRK1A and DYRK1B. Inhibits FMR1 binding to RNA.[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MG223278