

Product datasheet for **MG223299**

Hoxc5 (NM_175730) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hoxc5 (NM_175730) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Hoxc5
Synonyms: A1181837; Hox-3.4; Hox-6.2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG223299 representing NM_175730
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAGTCCTACGTAGCCAATTCATTCTATAAGCAGAGCCCAATATCCCTGCCTATAACATGCAAACCTT
 GTGGAACTATGGATCGGCCTCAGAGGTGCAGGCATCCAGGTACTGCTACGGCGGATTGGACTTAAGCAT
 CACTTTCCACCCTGCGCCTTCCAACCTCTCCACGGGGTAGACATGGCTGCCAACCCCGGGCTCAC
 CCCGACCGCCCGCTGCAGCGCCGCGCCGCTCCGGGACACGCTCTGGGCAGAGACGAAGCGGCTCCTC
 TGAACCTGGGATGTACAGTCAGAAGGGCGCTCGCCCGGCTGGAGGAGCGAGCTAAGAGCAGTGGGGA
 GATCAAAGAGGAGCAGGCGCAGACAGGGCAGCCTGCCGGACTGAGCCAGCCACCGGCCCGCCACAGATT
 TACCCGTGGATGACCAAACCTGCACATGAGCCACGAGACGGATGGCAAGCGGTCCCGAACAGTTACACGC
 GCTACCAGACCCTGGAACCTCGAGAAAGAATTCATTTAAACCGCTACCTCACTCGCCGAGGCGCATAGA
 GATTGCCAAACCTTGTGTCTCAACGAGAGACAGATCAAGATCTGGTTCCAGAACCAGGATGAAGTGG
 AAGAAAGATTCAAAATGAAGAGCAAAGAAGCTCTT

ACGCGTACGCGGCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG223299 representing NM_175730
 Red=Cloning site Green=Tags(s)

MSSYVANSFYKQSPNIPAYNMQTCGNYGSASEVQASRYCYGGLDLSITFPPPAPSNLSHGVDMAANPRAH
 PDRPACSAAAPGHALGRDEAAPLNPGMYSQKAARPALEERAKSSGEIKEEQAGTGPAGLSQPPAPPQI
 YPWMTKLHMSHETDGKRSRTSYTRYQTLELEKEFHFNRYLTRRRRIEIANNLCLNERQIKIWFQNRMMKW
 KKDSKMKKSKEAL

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_175730

ORF Size: 666 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175730.5](#), [NP_783857.1](#)

RefSeq Size: 2730 bp

RefSeq ORF: 669 bp


Locus ID: 15424

UniProt ID: [P32043](#)

Cytogenetics: 15 58.07 cM

Gene Summary: Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG223299