

Product datasheet for **MG223310**

Vmn1r193 (NM_134225) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Vmn1r193 (NM_134225) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Vmn1r193 |
| Synonyms: | V1ri9 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG223310 representing NM_134225 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCATGTTTACTTCCTTGTGTTTCTGAAAAGGTTTTCTAGTCTGAGCATGAAGATGATTTGGAGTG
ACCTCATCGAGGGGACAATTTTCCTTCACTTGTGGACTTGGAGTTTTAGGAAACAAAATCTTATTTGT
AAGACATTTGTATGCAGTCATCATGGGTCCTGAGAACAAAATTATAGATGTCATTCTCATCCACTTGGCT
TTTGTAACACAATCATTATTTATTGCATAGGGGTCAGAAACATAGCCACAGTTTTTACATCAGAACT
TCCTGGGTGATGTTGGTTGTAACATAATTTATCTAGAAAGGGTTGCCCGTGGCCTCTCCATCTGCAC
CACCTGTCTCCTCAGCGTGGTCCAGGCTGTCCACATCAGTCCCAGGACCACACTGTGGAGAAAGCTCAA
CCACAGACCACATGGCATGTTCTTGTCTTCTCATCCTCTTTGGATATTTAATCCCTGATAAGCTCCA
ACTTGCTTCACTATATCACAGCAGGCAGTAGCATGAACAGGCTGTAGTTGGGATGTTTACTGGGTATTG
CTATATGCTGCCATCCAGGCACACAGTTAAGTGGCTGTTCTTTCTTTCATGGCTCTTCGTGACGTCATT
TTCCAGGGTCTCATGGGCTGCAGCAGTGGGCTATGTCTCTCCGTCTGTACAAACATCATGTCGGTGCC
TCTATCTTACAGCTCAGGGTCTGCAACAATTCAGGCCAGAAATCAGAGCTACACAGAGGGTTCTCAC
TCTCATGACCTGTTTCTTTCTTTTATTTGGCAGATTTTCTCTTATACACAGGTTCCACTGTG
ACACGTGACTCCACAATACTAAATATTAAGCATTTTTAGTACTTAGTTATGCTGGTCTCAGCCCTTTTA
TCCTGATCATCAAGGATATTTGTGTTCTAAGACCTGCTGTGTCC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG223310 representing NM_134225
 Red=Cloning site Green=Tags(s)

MSCLLPVCSEKVFVLVSMKMIWSDLIEGTIFLSLVGLGVLGNKILFVRHLYAVIMGPENKIIDVILIH
 LA FVNTIIICYIGVRNIATVFYIRNFLGDVGCKTIIYLERVARGLSICTTCLLSVVQAVTISPR
 TTLWRKLPQTTHVFLVFLILFWIFNSLISSNLLHYITAGSSMNRSVVMFTGYCYMLPSRHTV
 KWLFLSFMALRDVIFQGLMGCSGSMRLRYKHHVRVLYLYSSGSANNSRPEIRATQVRVLTMT
 CFLFFYLADLIFSLYTGSTVTRDSTILNIKAFVLVSYAGLSPFILIIKDICVPKTCCVP

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_134225

ORF Size: 957 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_134225.1](#), [NP_598986.1](#)


RefSeq Size: 960 bp

RefSeq ORF: 960 bp

Locus ID: 171259

Cytogenetics: 13 A3.1

Product images:



Circular map for MG223310