

Product datasheet for **MG223418**

Pnkd (NM_019999) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pnkd (NM_019999) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pnkd
Synonyms:	2210013N15Rik; 2810403H05Rik; AI854243; Brp17; MNCb-5687; MR-1; Tahccp2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG223418 representing NM_019999 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTTGGCAGGGCTGGCTCGCGCCTTGGTTGTGGGTCTCCGGCTGTTGGCTGCTATTCTTCGCCTTCG
TGCTGCTCCTGAGTCCCCGAGCTGCCAAGAACAGCGGGTTTCCGCGCCTCCTCATGACGCGCAGCCA
GCGGCTGCTTCCGAATCGGGTATAGCCTGTACACCCGAACCTGGCTTGGGTACCTCTTCTACCGACAG
CAGCTTCGGAGGGCTCGGAATCGCTATCCCAAAGGCCACTCCAAAACCTCAGCCCCGCTCTTCAATGGAG
TGAAGGTGCTTCCATCCCTGTCTCTCGACAACACTACAGCTACCTCATCATCGACACCCAGGCTGGCT
GGCAGTGGCTGTGGACCCCTCAGACCCGAGGGCTGTGCAGGCTTCCATTGAAAAGGAAAGAGTAACTTG
GTTGCCATTCTCTGCACCCACAAGCACTGGGACCACAGTGGAGGGAACCGTGACCTCAGCCGGCGGCACA
GGGACTGTGAGTGTATGGGAGCCCTCAGGATGGTATCCCTACCTCACCACCCCTGTGTCATCAAGA
TGTGGTACGCGTGGGACGACTTCAGATTCGAGCCCTAGCCACCCCTGGCCACACTCAAGGCCATCTGGTC
TACCTGCTGGATGGGAGCCCTACAAAGTCTTCTGCTCTTCTCAGGGGACCTGCTCTTCTCTCTG
GCTGTGGACGGACCTTCGAAGGCACAGCGGAGACCATGCTGAGCTCACTGGACACTGTGTTAGACCTGGG
GGATGACACCCCTGCTGTGGCCTGGGCATGAATACGCAGAAGAGAACCTGGGCTTTCAGGGCGTGGTGGAG
CCCGAAGACTTGGCTCGTGAGAGGAAGATGCAGTGGGTACAGAGGCAACGGATGGAACGCAAGAGCAGCAT
GCCCATCCACCCTGGGAGAGGAACGCGCCTACAACCCATTCTGAGAACCATTGCTGGAACACTACAGGA
GGCTCTGGGGCCAGGCCAGGCCACCAGTGATGACGGCTGCTCCCGGGCCAGCTCCTGGAGGAGCTG
CGCCGCTGAAGGACATGCATAAAAGCAAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG223418 representing NM_019999
 Red=Cloning site Green=Tags(s)

MAWQGWLAPWLWVSGCWLLFFAFVLLLSPRSCQEQRGFRGLLMTRSQRLLFRIGYSLYTRTWLGYLFYRQ
 QLRRARNRYPKGHSKTQPRLFNGVKVLPVPLSDNYSYL I IDTQAGLAVAVDPSDPRVQASIEKERVNL
 VAILCTHKHWDHSGGNRDL SRRHRDCRVY GSPQDGI PYL THPLCHQDVVSVGRLQIRALATPGHTQGHLV
 YLLDGEPIYKGPSCLFSGDLLFLSGCGRTFEGTAETMLSSLDTVLDLGD TLLWPGHEYAEENLGFAGVVE
 PENLARERKMQWVQRQRMERKSTCPSTLGEERAYNPFLRTHCLELQEALGPGPGPTSDDGCSRAQLLEEL
 RRLKDMHKSK


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_019999

ORF Size:	1080 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_019999.2 , NP_064383.1
RefSeq Size:	2865 bp
RefSeq ORF:	1083 bp
Locus ID:	56695
UniProt ID:	Q69ZP3
Cytogenetics:	1 C3
Gene Summary:	Probable hydrolase that plays an aggravative role in the development of cardiac hypertrophy via activation of the NF-kappa-B signaling pathway.[UniProtKB/Swiss-Prot Function]