

Product datasheet for **MG223584**

Kcnip2 (NM_145704) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnip2 (NM_145704) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Kcnip2
Synonyms: KChI; KChIP2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG223584 representing NM_145704
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCGGGGCCAAGGCCGAAAGGAGAGTTTGTCCGAATCCCGAGATTTGGACGGCTCCTATGACCAGCTTA
CGGACAGCGTGGAGGATGAGTTTGAAGTATCCACGGTGTGCCACCGGCCTGAGGGTCTGGAACAACCTCA
GGAACAAACCAAGTTCACACGCAGAGAGTTGCAGGTCTGTACAGAGGCTTCAAGAACGAATGTCCAGC
GGAATTGTCAACGAGGAGAACTCAAGCAAATTTATTCTCAGTTCTTTCCCAAGGAGACTCCAGCAACT
ACGCTACTTTTCTTTCAATGCCTTTGACACCAACCATGATGGCTCTGTCAGTTTGGAGACTTTGTGGC
TGGTTTGTCAAGTATTCTTCGGGAACCATAGATGATAGACTGAACTGGGCTTTCAACTTATATGACCTC
AACAAGGATGGCTGTATCACGAAGGAGGAAATGCTCGACATCATGAAGTCCATCTATGACATGATGGGCA
AGTACACCTACCCTGCCCTCCGGGAGGAGGCCCGAGGGAACACGTGGAGAGCTTCTTCCAGAAGATGGA
CAGAAACAAGGACGGCGTGGTGACCATTGAGGAATTCATTGAGTCTTGTCAACAGGACGAGAACATCATG
AGGTCCATGCAACTCTTTGATAATGTCATC

ACGCGTACGCGGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG223584 representing NM_145704
Red=Cloning site Green=Tags(s)

MRGQGRKESLSESRDLGSDYDQLTDSVEDEFELSTVCHRPEGLEQLQEQTKFTRRELQVLYRGFKNECP
 GIIVNEENFKQIYSQFFPQGDSSNYATFLFNAFDTNHDGVSFEDFVAGLSVILRGTIDRLNWFNLYDL
 NKDGCITKEEMLDIMKSIYDMMGKYTPALREEAPREHVESFFQKMDRNDKGVVTIEEFIESCQQDENIM
 RSMQLFDNVI


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_145704

ORF Size: 660 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_145704.2 , NP_663750.1
RefSeq Size:	663 bp
RefSeq ORF:	663 bp
Locus ID:	80906
UniProt ID:	Q9JJ69
Cytogenetics:	19 38.75 cM
Gene Summary:	This gene encodes a member of the voltage-gated potassium channel-interacting protein (KCNIP) family. KCNIP family members are small calcium binding proteins that commonly exhibit unique variation at their N-termini, and which modulate A-type potassium channels. This gene is predominantly expressed in the adult heart, and to a lesser extent in the brain. Disruption of this gene is associated with susceptibility to cardiac arrhythmias and lack of transient outward potassium current in ventricular myocytes, and downregulated expression is associated with cardiac hypertrophy. The encoded protein has also been implicated as a repressor of immune response. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2013]