

Product datasheet for **MG223593**

Asph (NM_001177853) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Asph (NM_001177853) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Asph
Synonyms:	2310005F16Rik; 3110001L23Rik; AI848629; AW261690; AW561901; BAH; C79816; cl-37
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG223593 representing NM_001177853 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTGAAGATAAAGAGGCAAAGCACGGAGGACACAAGAATGGGAGAAGAGGAGGGATTCAGGAGGGT
CCTTTTTACATGGTTCATGGTCATTGCATTGCTCGGCGTCTGGACATCTGTGGCTGTGTTGTTGA
CTTGGTCGATTATGAAGAAGTTCTAGGAAACTAGGAGTCTATGATGCGGATGGCGATGGAGACTTTGAT
GTGGATGATGCCAAAGTTTTATTAGAAGGACCTGGTGGTTAGCCAAGAGGAAAATAAGGCTAAAGGCC
TTAAAGAAAGATCTCCTTCTGAGAGGACATTTCCACCGGAGGAGGAGGCAGAGACTCACGCTGAGCTGGA
GGAGCAGGCCCTGAGGGGGCAGACCTTGAAGCAGATGGGCTGGCAGGAGAACCACAGCCGAGGTTGAG
GACTTCCTCACAGTGACCGACAGTGATGACAGATTTGAGGACCTGGAACCCGGGACAGTTCATGAAGAAA
TCGAGGATACTTACCATGTGGAAGACACAGCATCGCAGAACCATCCAAATGACATGGAAGAGATGACGAA
TGAGCAGGAAAATTCAGATCCCAGTGAAGCAGTGACAGATGCAGGTGTGCTGCTGCCCATGCAGAGGAA
GTAAGACACCAAGACTATGATGAACCAGTATGAACCCTCAGAGCATGAAGGGGTCGGATTTTCAGATA
ACACCATAGATGATTCAGCATAATCTCAGAAGAAATAAATGTCGCCTCTGTTGAAGAACAGCAAGACAC
ACCACCAGATACT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG223593 representing NM_001177853
 Red=Cloning site Green=Tags(s)

MAEDKEAKHGHHKNGRRGGISGGSFFTFWFMVIALLGWVTSVAVVWFDLVDYEEVLGKLGVDADGDGDFD
 VDDAKVLLEGGGLAKRRTKAKGLKERSPSERTFPPEEEAETHAELEEQAPEGADLEADGLAGEPQPEVE
 DFLTVTDSDDRFEDLEPGTVHEEIEDTYHVEDTASQNHNPNDMEEMTNEQENSNDPSEAVTDAGVLLPHAAE
 VRHQDYDEPVYEPSEHEGVAISDNTIDDSIIIEEINVASVEEQDTPPDT

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001177853

ORF Size: 783 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177853.1](#), [NP_001171324.1](#)

RefSeq Size: 2693 bp

RefSeq ORF: 786 bp


Locus ID: 65973

UniProt ID: [Q8BSY0](#)

Cytogenetics: 4 A1

Gene Summary: Isoform 1: specifically hydroxylates an Asp or Asn residue in certain epidermal growth factor-like (EGF) domains of a number of proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG223593