

Product datasheet for **MG223826**

Nxn12 (NM_029173) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Nxn12 (NM_029173) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Nxn12
Synonyms: 4930519N16Rik; Rdcvf2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG223826 representing NM_029173
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGATCGCC

ATGGTGGACGTGCTGGGCGGGCGGCGCCTGGTGACCCGGGAGGGCACGGTGGTGGAGGCCGAGGTGGCGC
TGCAGAACAAAGGTGGTAGCTTTGTACTTTGCGGCGGGCCGGTGCTCGCCAGCCGCGACTTCACGCCGT
GCTCTGCGACTTCTACACGGAGCTGGTGAGCGAGGCGGCGGCCGCTCCCTTCGAGGTGGTTTTCTG
TCGGCAGACGGCAGTGCGGAGGAGATGTTGGACTTCATGCGCGAGCTGCACGGCTCCTGGCTGGCATTGC
CCTTCCACGACCCCTACCGGCATGAAGTGAAGAAGAGGTACGAAATCACCGCCATCCCCAAGCTGGTGGT
CATCAAGCAGAACGGAGCTGTCAACCAACAAAGGGCGGAAGCAGATCCGAGAGCGCGGCTAGCTTGC
TTTCAGAACTGGGTGGAAGCAGCCGATGTTTTCCAAAATTCTCGGGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG223826 representing NM_029173
Red=Cloning site Green=Tags(s)

MVDVLGGRRLVTREGTVVEAEVALQNKVVALYFAAGRCSPSRDFTPLLCDFYTELVSEARRPAPFEVVFV
 SADGSAEEMLDPMRELHGSWLALPFHDPYRHELEKKRYEITAIPKLVVIKQNGAVITNKGKQIRERGLAC
 FQNWVEAADVFQNFSG

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



View online »

Cloning Scheme:



ACCN: NM_029173

ORF Size: 468 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_029173.4](#), [NP_083449.1](#)

RefSeq Size: 1163 bp

RefSeq ORF: 471 bp

Locus ID:	75124
UniProt ID:	Q9D531
Cytogenetics:	13 A5
Gene Summary:	May be involved in the maintenance of both the function and the viability of sensory neurons, including photoreceptors and olfactory neurons. In the retina, isoform 1 may be required for rod function and isoform 2 for cone viability and function.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG223826