

Product datasheet for **MG224123**

Arhgef7 (NM_001113517) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arhgef7 (NM_001113517) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Arhgef7
Synonyms:	betaPix; betaPix-b; betaPix-c; Cool; cool-1; mKIAA0142; p85Cool1; p85SPR; Pak3bp; PIX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG224123 representing NM_001113517
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAATTCGCCGAGCAGACCGTTACGTGGCTCATCACCTGGGCGTGCTGGAGTCGCCAAGAAAACCA
 TCTCGGACCTGAGGTCTTCTTGCAGGCGTCGCTCAAGGACGGGTGGTCTATGCCGGCTGCTGGAGCG
 CCTGCTGCCCGGACCATCGAGAAAGTCTACCCGAGCCGCGGAACGAGAGCGAGTGCCTGAGCAACATC
 CGCGAGTTCCTGCGCGCCTGCGGGCGTCCCTGCGCCTAGAGACTTTTGATGCAAATGATTTGTATCAGG
 GGCAGAACTTTAACAAGTCTCAGCTCCTTGGTGACACTAAAATAAGTGACAGCAGACATTGGACTAGG
 AAGTGATTCTGTGTGTGCTCGGCCTTCGTCCATCGGATAAAGTCTTTTGATTCCCTGGGATCCCAGTCT
 TCACACAGTAGGACTTCAAATTTGCTTCAGAGCCAGTACCGAAGCTTGGACATGACTGATAACCAACA
 GCCAACTGGTAGTACGAGCCAAGTTAACTTCCAGCAGACCAATGAAGATGAACTCTCCTCTCAAAGGG
 TGATGTATCCATGTACACGAGTGGAGGAAGGAGGCTGGTGGGAGGGCACACACAATGGCAGGACCGGC
 TGGTTCGCCAGCAACTAGTTCGAGAGATCAAGCCAAGTGAGAAGCCCGTGTACCCAAATCAGGGACCT
 TGAAGAGCCCTCCCAAAGGGTTCGATACGACTGCCATCAACAAGAGCTATTACAACGTGGTGTACAGAA
 CATCCTGGAAACAGAGCATGAGTATTCCAAGGAGCTGCAGTCTGTGCTGTCCACCTACCTGCGGCCACTG
 CAGACCACTGACAAGTTGAGTTCAGCAAACACTTCATATTTAATGGGAAATCTAGAGGAAATATCTTCTCT
 TCCAGCAAGTGTGTACAGTCCCTAGAAGAATGCACGAAGTCTCCTGAAGCCCAACAGAGAGTTGGTGG
 CTGCTTCTGAGCCTGATGCCGCAGATGAGGACCTGTACCTCGTTACTGTGCCAACCCCATCTGCT
 GTGAGCGTCTCACAGACACAGTGGAGACCTAGGAGAGTTCATGGAAACAAAAGGTGCCAGCAGCCCTG
 GGATCCTGGTGTGACCACCGCCCTGAGCAAGCCCTTCATGCGCCTGGACAAGTACCCCACTGTGATAA
 GGAGCTGGAGAGACACATGGAGGATTATCATCCTGATAGACAAGATATTAGAAGTCTATGACGGCCTTC
 AAAAACCTTTTCAAGTGTCAAGAAGTTCGCAAGAGGAAGGAGCTGGAGCTGCAGATCCTGACGGAGC
 CCATCAGGAGCTGGGAGGGGATGACATAAAGACCCTGGGCAGTGTACATACATGTCCCAAGTCAACCAT
 TCAGTGTGCGGAAGCGAGGAGAAGAATGAGAGATACCTCCTGCTCTTCCCAAACCTTCTGCTCATGTTG
 TCTGCAAGTCCCAGGATGAGTGGTTTCATCTATCAGGGGAAGCTGCCAACACAGGAATGACAATCACAA
 AGCTTGAGGACAGTGAACCATAGGAATGCATTTGAGATATCAGGGAGCATGATCGAGCGGATTCTGGT
 GTCCTGCACCAGCCAGCAGGACTTACACGAGTGGGTGGAACACCTGCAGAAGCAGACGAAGGTACATCT
 GTGAGCAACCCCAACATCAAACCCCACTCGGTGCCATCACACACTTCCCTTCCATCCTCTCACTCCAT
 CCAGCAAACACGCGGACAGCAAGCCCGTGGCACTGACGCCTGCGTACCACACTCCCCACCCCTCTCA
 CCATGGCACCCCAACACACCACATCAGCTGGGGACCCCTGGAGCCTCCGAAGACCCCAAGCCTTGAGC
 CTGAGTTGCCTGCGGCCTGCACCTCCCCTCCGGCCCTCAGCTGCTCTGCTACAAGGAGGATCTCAGTA
 AGAGCCCCAAGACCATGAAAAAGCTGCTGCCAAGCGCAAGCCCGAGCGGAAGCCTTCGGACGAGGAGTT
 CGCTGTGCGCAAGAGCACAGCGCGCTGGAAGAAGACGCTCAGATCCTGAAGGTTATCGAAGCTTATTGC
 ACAAGTGCAAAGACGCGCCAGACCCTGAACTCAACATGGCAAGGCACTGACCTGATGCATAATCACGTCT
 TGGCTGATGACGACCAATCAAGTCTAGACTCCTTGGGTGTCGAGTACGCTTTCTGTTGGAGCCCTC
 AGACCTCTCGGAAGACTCTGAGTATGACAGTATATGGACAGCCCATAGTTACAGAATGGGTTCTGCATCC
 CGTAAGAGCTGCTCATATATCTCTCACCAGAAC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG224123 representing NM_001113517
 Red=Cloning site Green=Tags(s)

```
MNSAEQVTWLITLGVLESPKKTISDPEVFLQASLKDGVVLCRLLELLPGTIEKVYPEPRNESECLSNIR
REFLRACGASLRLETFDANDLYQGQNFNKVLSLVTLNKVTADIGLGSDSVCARPSSSHRIKSFDSLGSQS
SHSRTSKLLQSQYRSLDMTDNTNSQLVVRAKFNFQQTNEDELSFSGGDVIHVTRVEEGGWEGTHNGRTG
WFPSNYVREIKPSEKPVSPKSGTLKSPPKGFDTTAINKSYYNVVLQNILETEHEYSKELQSVLSTYLRPL
QTSKLSANTSYLMGNLEEISSFQQVLVQSLEECTKSPEAQQRVGGCFLSLMPQMRTLYLAYCANHPSA
VSVLTEHSEDLGEFMETKGASSPGILVLTGLSKPFMRDKYPTLLKELEHRHMEDYHPDRQDIQKSMTAF
KNLSAQQCQEVRRKRELELQILTEPIRSWEGDDIKTLGSVTYMSQVTIQCAGSEEKNERYLLLFPNLLLML
SASPRMSGFIYQGKLPPTGMTITKLESENHRNAFEISGSMIERILVSCTSQQDLHEWVEHLQKQTKVTS
VSNPTIKPHSVPSHTLPSHPLTPSSKHADSKPVALTPAYHTLPHPSHHGTPHTTISWGPLEPPKTPKPWS
LSCLRPAAPLRPSAALCYKEDLSKSPKTMKLLPKRKPERKPSDEEFVARKSTAAL EEDAQILKVIEAYC
TSAKTRQTLNSTWQGTDLMHNHVLADDDQSSLDLGRSSLSRLEPSDLSEDSEYDSIWAHSYRMGSAS
RKSCCSYISHQN
```


TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2020_e10.zip

Restriction Sites: Sgfl-Mlul

OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001113517.2
RefSeq Size:	4946 bp
RefSeq ORF:	2349 bp
Locus ID:	54126
UniProt ID:	Q9ES28
Cytogenetics:	8 A1.1
Gene Summary:	Acts as a RAC1 guanine nucleotide exchange factor (GEF) and can induce membrane ruffling. May function as a positive regulator of apoptosis. Functions in cell migration, attachment and cell spreading. Promotes targeting of RAC1 to focal adhesions. Downstream of NMDA receptors and CaMKK-CaMK1 signaling cascade, promotes the formation of spines and synapses in hippocampal neurons (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG224123