

Product datasheet for **MG224273**

Hmga1 (NM_001166537) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Hmga1 (NM_001166537) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Hmga1
Synonyms: AL023995; Hmg; HMG-I(Y; Hmga; Hmga1a; Hmga1b; Hmgj; HMGI(Y; Hmgij; Hmgy
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG224273 representing NM_001166537
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGAGTCGGGCTCAAAGTCCAGCCAGCCCCTGGCCTCCAAGCAGGAAAAGGATGGGACTGAGAAGC
GAGGCCGGGCAGGCCACGCAAGCAGCCTCCGAAAGAGCCAGTGAAGTGCCAACCCGAAGAGACCTCG
GGGCCACCAAAGGGAAGCAAGAATAAGGGCGCCCAAGACCCGAAAGTCACCACAGCTCCAGGGAGG
AAACCAAGGGCAGACCAAGAACTGGAGAAGGAGGAAGAGGGGCATCTCCAGGAGTCTCTGAGG
AGGAGCAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG224273 representing NM_001166537
Red=Cloning site Green=Tags(s)
MSESGSKSSQPLASKQEKDGTEKRGRGRPRKQPPKEPSEVPTPKRPRGRPKGSKNKGAATRKYVTTAPGR
KPRGRPKKLEKEEEEIGISQESSEEEQ

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001166537

ORF Size: 288 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166537.1</u> , <u>NP_001160009.1</u>
RefSeq Size:	1971 bp
RefSeq ORF:	291 bp
Locus ID:	15361
UniProt ID:	<u>P17095</u>
Cytogenetics:	17 14.5 cM
Gene Summary:	This locus encodes a member of the nuclear, non-histone high mobility group protein family. This architectural transcription factor binds to A-T rich DNA sequences and participates in enhanceosome formation, chromatin remodeling and regulation of transcription. This protein functions in many cellular processes, including cell growth and differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Oct 2009]