

Product datasheet for **MG224295**

Hif3a (NM_016868) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hif3a (NM_016868) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Hif3a
Synonyms:	bHLHe17; lpas; MOP7; NEPAS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG224295 representing NM_016868
 Red=Cloning site Blue=ORF Green=Tags(s)

GACGTTGTATACGACTCCTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGGCGCGCC

ATGGACTGGGACCAAGACAGGTGCAACACCCAGCTGCGGAAGGAGAAGTCGCGGGACCGCGCCCGCAGCC
 GGCGCAGCCAGGAGACGGAGGTGCTGTACCAGCTGGCGCACACTTGCCCTTTGCGCGCGGCGTACGCGC
 GCACCTGGACAAGGCCTCCATCATGCGCCTCACAATCAGCTACCTGCGCATGCACCGCCTCTGCGCAGCA
 GGGGAGTGAACAGGTGAAAAAGGGGAGAGCCACTGGACGCCTGCTACCTGAAGGCCCTGGAGGGTT
 TCGTCATGGTACTACCGCCGAGGGAGACATGGCTTACCTGTGCGAAAAATGTCAGCAAGCACCTGGGCCCT
 CAGTCAGCTGGAGCTCATTGGACACAGTATCTTTGATTTTATCCATCCCTGTGACCAAGAGGAACCTCAA
 GACGCCCTGACCCCGAGCCGAACCTGTCAAAGAAGAAGCTGGAAGCCCCAACAGAGCGCCACTTTTCCC
 TGCGAATGAAGAGCACGCTCACCAGCAGAGGGCGCAGCTCAACCTCAAAGCGGCCACCTGGAAGGTGCT
 GCACTGCTCAGGACATATGAGGGCTACAAGCCCCTGCACAGACTTCCCCTGCCGGGAGCCCTCGTCC
 GAGCCTCCCCTGCAATGCCTGGTGTATCTGTGAAGCCATCCCCACCCAGCCAGTCTGGAGCCCCCGC
 TGGGCCGAGGGGCTTTCTCAGTCGCCACAGCCTGGACATGAAGTTCACATACTGCGACGAGAGGATTGC
 AGAAGTTGCTGGCTACAGTCTGATGACCTGATTGGCTGTTCTGCCTATGAATACATCCACGCTTTGGAC
 TCTGATGCGGTACAGCAGGAGCATCCACACTTTGTTGAGCAAGGGCCAGGCAGTAACGGGGCAGTATCGCT
 TCCTGGCCCGACTGGAGGCTATCTGTGGACTCAGACTCAGGCTACAGTGGTGTGAGGGGGCGGGGCC
 CCAGTCGGAAGTATCATCTGCGTCCACTTCTGATCAGCCGTGTAAGAAGACCCGAGTGGTGTGTCT
 CTGGAACAAACGGAGCAACATACTCGCAGACCCCCTCGGCTGAGTGCCTCCTCGCAGAAGGGTATCCCTG
 CCAACAGTGTAGACTCTCCTGCTCCGCGGATCCTGGCCTTCTGCACCCCTCCGGCCCTGAGTGGGCCCTC
 CCTGGCTGCTGACCCTCGCCGTTTCTGTAGTCCAGACCTGCGCCGCTCATGGCACCCATCCTGGATGGA
 CCTCCCCAGCTGCCACGCCAGCACCCACAAGCTACACGAGACCCCAAGTCTCTTCCGGCTGATC
 TCCCAGATAAGTTGGCAGTGGGCTTGGAAAATGCACACAGACTCTCCACTGCCAGAAAAACAAGACCGT
 GGAGACAGATCTAGATATAGCTCAGGACTCTGACACTCTGGACTTGGAGATGCTGGCCCCCTACATCTCC
 ATGGATGATGACTTCCAGCTCACTCCAGTGAAGCAATTGCCAAAGTCCACCGCAGACCTCCAGGGTGG
 CCCGAGGCCCGTGCTCGGAGCTTCCATGGCCTGTGCTCCTATCCCTGAGCCCTCCCTACTGCCCG
 CTGGGGGAGTGATCCACGACTGAACTGTTCCAGTCTTCCAGGGCGATCGCCCCACAGCCTCCCTGATG
 CCTGGAACCTCGGAAGAGGGCCTTGGCCCAGAGCTCAGAGGACAAAGGGTTGGAGCTGCTGGAACCTAAGC
 CTCCCAAGCGGTCCCAAGACTAGAACCTGGAAGCTTCTGCTGCCTCCGCTCAGCCTGAGTTTCTTCT
 GCAAGGTCGACAACCTCTGGAAACCAGCAGGATCCCAGAGCCCCCTCGTGCATTCTCATGAGCCCTTG
 GGCTAGCTCCCTCGCTGCTCTCTCTGCCAGCATGAGGAAACTGTCCAGCCCAGGAACCGCTTCCCGC
 CAGCAGCAGGCTTGGCCAGACCCAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG224295 representing NM_016868
 Red=Cloning site Green=Tags(s)

```
MDWDQDRSNTLNRKEKSRDAARRSRSSQETEVLYQLAHTLPFARGVSAHLDKASIMRLTISYLRMHRLCAA
GEWNQVEKGGPELDACYLKALEGFVMVLTAEGDMAYLSENVSKHLGLSQLELIGHISIFDFIHPCDQEELQ
DALTPRPNLSKKKLEAPTERHFLSRMKSTLTSRGRTLNLKAATWKVLHCSGHMRYKPPAQTSAPGSPRS
EPPLQCLVLICEAIPHPASLEPPLGRGAFLSRHSLDMKFITYCDERIAEVAGYSPDDLIGCSAYEYIHALD
SDAVSRSIHTLLSKGQAVTGQYRFLARTGGYLWTQTQATVVSGGRGPQSESIICVHFLISRVEETGVVLS
LEQTEQHTRRPPRLSASSQKGIPIGNSVDSAPRILAFHPPALSEASLAADPRRFCSPLDRLMAPILDG
PPPAATPSTPQATRRPQSPLPADLPDKLAVGLENAHRLSTAQKNKTIVETDLIAQDSDTLDLEMLAPYIS
MDDDFQLNSSEQLPKVHRRPPRVARRPRARFHLSPPIPEPSLLPRWGS DPRLNCSSPSRGRDPTASLM
PGTRKRALAQSSDKLELLETKPKRSPRLEPGSFLPLSLFLLQGRQLLGNQDPPRAPLVHSHEPL
GLAPSLLSLCQHEETVQPRNRFPPAAGLGQTH
```


TRTRPLE - GFP Tag - V

Restriction Sites: AscI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_016868

ORF Size: 1986 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016868.3](#), [NP_058564.2](#)

RefSeq Size: 6947 bp

RefSeq ORF: 1989 bp

Locus ID: 53417

UniProt ID: [Q0VBL6](#)

Cytogenetics: 7 A2

Gene Summary: Acts as a transcriptional regulator in adaptive response to low oxygen tension. Acts as a regulator of hypoxia-inducible gene expression (PubMed:9840812, PubMed:11734856, PubMed:21546903). Plays a role in the development of the cardiorespiratory system (PubMed:18070924).[UniProtKB/Swiss-Prot Function]