

## Product datasheet for **MG224620**

### Barhl1 (NM\_019446) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Barhl1 (NM\_019446) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Barhl1  
**Synonyms:** Dres115; MBH2  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG224620 representing NM\_019446  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAAGGCTCTAATGGCTTTGGGATTGACTCCATTCTCTCCACCGAGCGGGCAGCCCCGCCCTTCCCA  
 AGGGGGACCCCTTGCTTGGGGACTGCCGTTACCCCTGGAGCTGAGTCCACGCTCAGAGAGCAGTAGCGA  
 CTGCTTTCACCAGCCTCGCCCGAAGAGACTGTCTGGAGACCAGTACCTCGCGGCTGGTGCAGCATCT  
 GGCCAGGTTTGGACTCTCACCTGCAACCGGGCAGCTTTCAGCCCCAGCCAGTCGCGAACTGTCACT  
 CCTCCTTTCTGATCAGGGACATCCTTGCTGACTGCAAACCTCTCGCGCCTGTGCACCCTACTCTAGCAG  
 TGGGCAGCTGCAGCCCTGAGCCTGGGGCCGCCTTGCAGCCAAAGCCGGGGAGGACTTTCGAGACAAG  
 CTGGACAAAAGTGTGAGCAGCGTTCATCAGACTCTGAGTACAAAAGTGAAGGAGGAGGGCAGCCGCGAGA  
 TCTCCAGCTCTCGGGACAGTCCCCCGGTGCGCCTGAAAAAGCCACGCAAAGCACGCCCTTACCGA  
 CCATCAGCTGGCGCAGCTGGAGCGTAGCTTCGAGCGGCAGAAATACCTGAGCGTGAAGACCGCATGGAG  
 CTAGCCGCTCGCTCAACCTCACCGACACGAGGTCAAGACCTGGTACCAGAACCGCAGGACTAAATGGA  
 AGCGACAGACGGCCGTGGGGCTGGAGCTGCTAGCGGAGGCAGGCAATTAATCGGGCTCCAGAGAATGTT  
 CCCGTGCGCTTATTTCTACCCGACAGTCTGGTTTCCAATCTGGACCCCGGTGCCGACTCTATCTGTAC  
 CGCGGACCCAGCGGCCACCTCCTGCCCTCAGAGACCTCTGGTGCCCGCATCTCATCCACGGACTCC  
 AGGGCGCCAGCGAGCCGCCCGCCGCTGCCCGCTGCCCGGCTCCTCCCGCGCCGCGCAGCCCCG  
 G

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

**Protein Sequence:** >MG224620 representing NM\_019446  
 Red=Cloning site Green=Tags(s)

```
MEGSNGFGIDSILSHRAGSPALPKGDPLLGDRCSPLELSRSESSSDCSSPASPGRDCLSTSRPGAAS
GPGLDSHLQPGQLSAPAQSRTVTSFLIRDILADCKPLAACAPYSSSGQPAPEPGGRLAAKAGEDFRDK
LDKSVSSASSDSEYKVKKEEGDREISSSRDSPVRLKKPRKARTAFDHLQAQLERSFERQKYL SVQDRME
LAASLNL TDTQVKTWYQNRRTKWKRQTAVGLELLAEAGNYSALQRMFSPYFYPQSLVSNLDPGAALYLY
RGPSAPPPALQRPLVPRILIHGLQGASEPPPLPPLPGVLPRAAQPR
```

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_019446

**ORF Size:** 981 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_019446.4</a> , <a href="#">NP_062319.1</a>
<b>RefSeq Size:</b>	3674 bp
<b>RefSeq ORF:</b>	984 bp
<b>Locus ID:</b>	54422
<b>UniProt ID:</b>	<a href="#">P63157</a>
<b>Cytogenetics:</b>	2 19.38 cM