

Product datasheet for **MG224791**

H2-K1 (NM_001001892) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	H2-K1 (NM_001001892) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	H2-K1
Synonyms:	H-2K; H-2K(d); H2-K; K-f
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG224791 representing NM_001001892 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGTACCGTGCACGCTGCTCCTGCTGTTGGCGGCCGCCCTGGCTCCGACTCAGACCCGCGGGGCCAC
ACTCGCTGAGGTATTTTCGTACCCGCCGTGTCCCGCCCGCCTCGGGGAGCCCCGGTACATGGAAGTCGG
CTACGTGGACGACACGGAGTTCGTGCGCTTCGACAGCGACCGGAGAATCCGAGATATGAGCCGCGGGCG
CGGTGGATGGAGCAGGAGGGGCCGAGTATTGGGAGCGGGAGACACAGAAAGCCAAGGGCAATGAGCAGA
GTTTCCGAGTGGACCTGAGGACCTGCTCGGCTACTACAACCAGAGCAAGGGCGGCTCTCACAATTCA
GGTGATCTCTGGCTGTGAAGTGGGTCCGACGGGCGACTCCTCCGCGGTACCAGCAGTACGCCTACGAC
GGCTGCGATTACATCGCCCTGAACGAAGACCTGAAAACGTGGACGGCGCGGACATGGCGCGCTGATCA
CCAAACACAAGTGGGAGCAGGCTGGTGAAGCAGAGAGACTCAGGGCTACCTGGAGGGCAGCTGCCTGGA
GTGGCTCCGAGATACCTGAAGAACGGGAACGCGACGCTGCTGCGCACAGATTCCCCAAGGCCCATGTG
ACCCATCACAGCAGACCTGAAGATAAAGTACCCTGAGGTGCTGGGCCCTGGGCTTCTACCTGCTGACA
TACCCCTGACCTGGCAGTTGAATGGGAGGAGCTGATCCAGGACATGGAGCTTGTGGAGACCAGGCTGC
AGGGATGGAACCTTCCAGAAGTGGCATCTGTGTTGCTCCTTGGGAAAGGAGCAGTATTACACATGC
CATGTGTACCATCAGGGCTGCCTGAGCCCTCACCTGAGATGGGAGCCTCCTCCACTGTCTCCA
ACATGGCGACCGTTGCTGTTCTGTTGTCTTGGAGCTGCAATAGTCACTGGAGCTGTGGTGGCTTTTGT
GATGAAGATGAGAAGGAGAAACACAGGTGAAAAAGGAGGGGACTATGCTCTGGCTCCAGGCTCCAGACC
TCTGATCTGTCTCTCCAGATTGAAAAGTATGTTTCATGACCCTCATTCTTAGCG

ACGGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG224791 representing NM_001001892
 Red=Cloning site Green=Tags(s)

MVPCTLLLLLAAALAPTQTRAGPHSLRYFVTAVSRPGLGEPYMEVGYVDDTEFVRFDSDAENPRYEPRARWMEQEGPEYWERETQKAKGNEQSFVRDLRLLGYNQSKGGSHTIQVISGCEVGS DGRLLRQYQYAYD GCDYIALNEDLKTWTAADMAALITKHKWEQAGEAERL RAYLEGT C VEWL RRYL KNGNATLLRTDSPKAHV THHSRPEDKVT LRCWALGFYPADITL TWQLNGEELIQDMELVETRPAGDGT FQK WASVVVPLGKEQYYTC HYYHQGLPEPLTLRWEPPPSTVSNMATVAVLVVLGAAIVTGAVVAFVMKMRRRNTGGKGGDYALAPGSQTSDL SLPDCKVMVHDPHSLA

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2586_e03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001001892

ORF Size: 1107 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001001892.2](#), [NP_001001892.2](#)

RefSeq Size: 1669 bp

RefSeq ORF: 1110 bp


Locus ID: 14972

UniProt ID: [P01901](#)

Cytogenetics: 17 17.98 cM

Gene Summary: Involved in the presentation of foreign antigens to the immune system.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG224791