

Product datasheet for **MG224944**

Prr3 (NM_145487) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Prr3 (NM_145487) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Prr3
Synonyms: 4930540G07Rik; Cat56
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG224944 representing NM_145487
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCGAAACGAAAGAAGCAGGACCAGCCGCCGCTGCCGAGCAGCAACATCTCGCACTGTCTG
 AGCGGGACGAGCCTGGAGATGAGGAGGATGAGCGCCCATGGGACCACCCAGCCTTCTGGGCCCTCCCC
 CATGGCTAACGGGAAGCCTGGTATCCCAAGTCAGCTTCCACAGGGTCTCCAGGATCAAGGGGACGG
 ATGATTCCTCCACTGCTGAGTCTCCCGCTCCTCCCGGGGTAGAGGCTACATTCGGGGAGGCTCGGTC
 CTCGATAGCCCGTATGGTCGTGGTGGTGGGTATCAATGCCGAACCTCCTTTCTGGCCAGGCCA
 TGGGGTCCCTTAGGAAAGCTTTTACAAAGAGGCGAGAAACCCACGAAGGCTCAGAAGCTGGTCTCTT
 GTCAAGAATACCTACCCACCAAGGACAGCCCCAGATGATGGAAGACAAGTCTGACCGCCTGTCTGCC
 GACACTTTTCTAAGAAGGGCCACTGTGCATATGAAGACCATTGTGCCTTCTACCATCCAGGCGTCAACGG
 ACCTCCTCTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG224944 representing NM_145487
 Red=Cloning site Green=Tags(s)

MPKRKKQDQPPPLPQQQHLALSERDEPGDEEDERPMGPPSLLGPPPMANGKPGDPKSAFHRGPPGSRGR
 MIPPLLSLPPPRGRGYIRGGLGPRSSPYGRGWGINAEPPFPGPGHGGPSRESFYKEARNPRRLRSWSL
 VKNTYPPKDSPPQMMEDKSDRPVCRHF SKKGHCYEDHCAFYHPGVNPPPL

TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



Cloning Scheme:



ACCN: NM_145487

ORF Size: 570 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_145487.4](#)

RefSeq Size: 1970 bp


RefSeq ORF: 573 bp

Locus ID: 75210

UniProt ID: [Q811B5](#)

Cytogenetics: 17 B1

Product images:



Circular map for MG224944