

## Product datasheet for **MG225000**

### Cabp7 (NM\_138948) Mouse Tagged ORF Clone

#### Product data:

|                           |                                                                             |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                         |
| Product Name:             | Cabp7 (NM_138948) Mouse Tagged ORF Clone                                    |
| Tag:                      | TurboGFP                                                                    |
| Symbol:                   | Cabp7                                                                       |
| Synonyms:                 | MGC129428; MGC129429                                                        |
| Mammalian Cell Selection: | Neomycin                                                                    |
| Vector:                   | pCMV6-AC-GFP (PS100010)                                                     |
| E. coli Selection:        | Ampicillin (100 ug/mL)                                                      |
| ORF Nucleotide Sequence:  | >MG225000 representing NM_138948<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGCCGTTCCACCCGGTGACGGCGCGTGTGATGTACGGGGCATCTACACCGTGCCCAACCTGCTGTGG  
 AGCAGCGCCCGGTGGATATCCCTGAGGACGAGCTAGAGGAGATCCGAGAGGCCTTCAAGGTATTCGACCG  
 GGATGGCAACGGCTTCATCTCTAAGCAGGAGCTGGGAACAGCCATGCGCTCCCTGGGTTACATGCCCAAC  
 GAGGTGGAGCTGGAGTTATCATCCAGCGGCTAGACATGGATGGTGACGCCAAGTGACTTTGAGGAGT  
 TTGTAACCTCCTGGGACCAAGCTCTCCACCTCGGGATCCCAGAGAAGTTCCATGGCACTGACTTTGA  
 CACGGTCTTCTGGAAGTGTGACATGCAGAAGCTGACAGTGGATGAGCTGAAGCGGCTTCTCTATGACACC  
 TTCTGCGAGCACTTGCCATGAAGGACATAGAGAACATCATCATGACAGAGGAGGAGAGCCACCTGGGCA  
 CTGCTGAGGAGTGTCCCGTGGACGTGGAGACCTGTTCCAACCAGCAGATCCGCCAGACGTGTGTGCGCAA  
 GAGCCTGATCTGCGCCTTCGCCATCGCCTTCATCATCAGTGTCATGCTCATCGCGCCAACCAAGTGTCTG  
 CGCAGCGGCATGAAG

**ACGCGT**ACGCGGCGGCTCGAG - GFP Tag - GTTTAA


[View online »](#)


**Protein Sequence:** >MG225000 representing NM\_138948  
 Red=Cloning site Green=Tags(s)

MPFHPVTAALMYRGIYTPNLLSEQRPVDIPEDELEEIREAFKVFDRDGNFISKQELGTAMRSLGYMPN  
 EVELEVIQRLMDGDGQVDFEEFVTLGPKLSTSGIPEKFHGTDFDTVFWKCDMQKLTVDLKRLLYDT  
 FCEHLSMKDIENIIMTEEEHSLGTAEECPVDVETCSNQIRQTCVRKSLICAFIAFIISVMLIAANQVL  
 RSGMK

TRTRPLE - GFP Tag - V

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**ACCN:** NM\_138948

**ORF Size:** 645 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_138948.4](#)

**RefSeq Size:** 699 bp

**RefSeq ORF:** 648 bp


**Locus ID:** 192650

**UniProt ID:** [Q91ZM8](#)

**Cytogenetics:** 11 A1

**Gene Summary:** Negatively regulates Golgi-to-plasma membrane trafficking by interacting with PI4KB and inhibiting its activity.[UniProtKB/Swiss-Prot Function]

## Product images:



Circular map for MG225000