

Product datasheet for **MG225011**

Zfp784 (NM_001039532) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp784 (NM_001039532) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Zfp784
Synonyms:	C730049F21Rik; Znf784
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG225011 representing NM_001039532 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCCGCCGCGCGCCCGGATCCCCGATTCCGAGCTACCGACCCGGGAGTCGCCATCCCCGGAGCCGC
CGGACCTGGTCCTGGTGCCTGATGGTCGCCCTGTACACCCCCAGGTGACCTCATCGAGATCCAGGTGGT
AAAAGTCACAGATACCACATCGGTCCCTGAGCCCCAGAGCCAGGATCTTTTCATTGTGCCCTGTGTCCT
GCTGCCTTCCGTCTAGTCTCAGAGCTGCTGTTTCATGAGCATGGCCACCTGGCAAGCATGGAGGGGCTTG
GGCAGGATGGGGATCCTAGCCGTTGTCTGTGTGTGGCCACAGCTGTCCAGGTCTGCTAGCCTGCGTGC
CCATTACAGCTTGACACAGGAGAACGGCCTTACCGCTGTTCCCTCTGCCACGGGCTTTTAAGGCCTTG
GCACCTCTGCTCCGTACCAACACCGACATGGGGTAGAGCCCGGCACCTCTGAAAGGCTTCTACCCACAA
CAACAACCGGACAGCCAACTCAAGGGTGGCCCAGGAGAGGTGGAGGTGGTGTGGCTGCAGCAGCTGC
GGGTGCTGTGGTAGGCAAGCCTTTCGCCTGCAGGTTTTGTGCCAAGCCCTCCGACGCTCCTCAGATATG
CGAGACCACGAGCGGGTGCATACAGGGGAGCGGCCCTACCACTGTAGTATCTGTGGCAAAGGCTTACGC
AGTCCTCTGTGCTGAGTGGTCATGCCCGCATCCACTGGCGAGCGCCCTTCCGATGCATGCTCTGTGA
CCGCACTTTCAATAACTCCTCAACTTCCGCAACACAGCGCACCCACTTCCATGGACCAGGGTCTGGA
GTGGGAGAGTCTAGAGGCCAACTGAGATCATCTTCAGTGTCCCAAGATCA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)


Protein Sequence: >MG225011 representing NM_001039532
 Red=Cloning site Green=Tags(s)

MAAAPDPPIPPSSPTRESPSPPEPDLVLPDGRPVTPPGDLIEIQVVKVTDTTSVPEPEPGSFHCALCP
 AAFRLVSELLFHEHGLASMEGLGQDGDPSRCHVCGHSCPGPASLRAHYSLHTGERPYRCSLCPRAFKAL
 APLL RHQHRHGV EPGT SERLLPTTTTGQPN SRVAQERSEVVM AAAAGAVVGKPFACRFCAKPFRRSSDM
 RDHERVHTGERPYHCSICGKGFTQSSVLSGHARIHTGERPFRMCLCDRTFNSSNFRKHQRTHFHGPGSG
 VGESRGQLRSSVSQES

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001039532

ORF Size: 891 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001039532.2](#), [NP_001034621.1](#)

RefSeq Size: 2146 bp

RefSeq ORF: 894 bp


Locus ID: 654801

UniProt ID: [Q8BI69](#)

Cytogenetics: 7 A1

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG225011