

Product datasheet for **MG225339**

Kcnq2 (NM_001006674) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kcnq2 (NM_001006674) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Kcnq2
Synonyms:	HNSPC; KQT2; Nmf134
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG225339 representing NM_001006674
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGTGCAGAAGTCGCGCAACGGTGGCGTGTACCCGGCACCAGCGGGGAAAAGAAGCTCAAGGTGGGCT
 TCGTGGGGCTGGACCCCGCGCGCCGACTCCACACGCGACGGCGGCTACTCATCGCGGCTCCGAGGC
 CCCAAGCGCGGAGCGTTTTGAGCAAGCCGCGACGGCGCGGGAGCCGGGAAGCCCCGAAGCGC
 AACGCCTTCTACCGCAAGCTGCAGAATTTCTCTACAACGTGCTAGAGCGGCCCGCGGCTGGCGTTCA
 TCTACCACGCCTACGTGTTCTTTTAGTCTTCTCTGCTTGTGCTTTCTGTGTTTTCCACCATCAAGGA
 GTACGAGAAGAGCTCTGAGGGGGCCCTACATCTTGGAAATCGTACTATCGTGGTATTCGGTGTGAG
 TACTTTGTGAGGATCTGGGCTGCAGGCTGCTGTTGCCGGTATCGAGGCTGGAGGGCAGGCTCAAGTTG
 CCAGGAAGCCGTCTGTGTGATTGATATCATGGTGTGATTGCCTCCATTGCTGTGCTGGCTGCTGGTTC
 CCAGGGCAATGTCTTTGCCACATCTGCGCTTCGGAGCTTGGGTTCTTGCAAATCTTGGGATGATCCGT
 ATGGACCGGAGGGGTGGCACCTGGAAGCTCTGGGATCGGTAGTCTACGCTCACAGCAAGGAGCTGGTGA
 CTGCCTGGTACATTGGCTTCTCTGCTCATCCTGGCCTCATTTCTGGTGTACTTGGCAGAAAAGGGTGA
 GAATGACCACTTTGACACCTACGCAGATGCACTCTGGTGGGGTCTGATCACCTGACGACCATTTGGCTAC
 GGGGACAAGTACCCTCAGACCTGGAACGGGAGGCTGCTGGCAGCGACCTTACCCTCATTGGTGTCTCGT
 TCTTTGCTCTTCTGCTGGCATTTTGGGATCCGGCTTGGCCCTGAAAGTCCAAGAGCAGCATCGGCAAAA
 ACACCTTGAGAAACGGCGGAACCTGCGGCAGGTCTGATCCAGTCTGCCTGGAGATTCTATGCTACTAAC
 CTCTCACGCCACCGACCTGCACCTCACGTGGCAGTACTACGAGCGGACAGTCACTGTCCCCATGTACAGAC
 TCATCCCACCTCTGAACCAGCTGGAGCTGCTGAGGAATCTCAAGAGCAAATCTGGACTCACCTTCAGGAA
 GGAGCCACAGCCAGAGCCATCACCAAGTCAGAAGGTGAGTTTGAAGATCGTGTCTTCTCCAGCCCCGA
 GGATGGCTGCCAAGGAAAGGGGTCTCCCCAGGCCAGACGGTCCGGCGGTCCCCAGTCCGGATCAGA
 GTCTTGATGACAGCCCGAGCAAGGTGCCAAGAGCTGGAGCTTTGGTGGCCGAGCCGACACGCCAGGC
 TTTCCGCATCAAGGGTGTGATCCCGCAGAATTGAGAAGAAGCAAGCCTCCCTGGGGAGGACATCGTA
 GAGGACAACAAGAGCTGTAAGTGCAGTTTGTGACTGAAGATCTTACCCTGGCCTCAAAGTTAGCATCA
 GAGCTGTGTGTTATGCGGTTCTTGGTATCTAAGCGAAAGTTCAAAGAGAGTCTGCGCCCATATGATGT
 GATGGACGTATCGAACAGTACTCGGCTGGACACTTGGATATGTTGTCCCGCATCAAGAGCTGCAGTCC
 AGGTCCTGCGACTGGAGAGGAGTCTGGCA

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence:

>MG225339 representing NM_001006674
 Red=Cloning site Green=Tags(s)

MVQKSRNGGVYPGTSGEKCLKVGFVGLDPGAPDSTRDGALLIAGSEAPKRGSVLSKPRTGGAGAGKPPKR
 NAFYRKLQNFLYNVLERPRGWAFIYHAYVFLLVFSLVLSVFSSTIKEYESSEGALYILEIVTIVVFGVE
 YFVRIWAAGCCRYRGRWGRGLKFARKPFCVIDIMVLIASIAVLAAGSQGNVFATSALRSLRFLQILRMIR
 MDRRGGTWKLLGSAVYHASKELVTAWYIGFLCLILASFLVYLAEKGENDFDHYADALWWGLITLTTIGY
 GDKYPQTWNGRLLAATFTLIGVSFFALPAGILGSGFALKVQEQHRQKHFEKRRNPAAGLIQSAWRFYATN
 LSRTDLHSTWQYYERTVTVPMYRLIPPLNQLLELLRNLKSKSGLTFRKEPQPEPSPSQKVLKDRVSSPR
 GMAAKGKQSPQAQTVRRSPSADQSLDDSPSKVPKSWSFGDRSRTRQAFRIKGAASRQNSEEASLPGEDIV
 EDNKSCNCFVTEDLTPGLKVSIRAVCMRFLVSKRKFKESLRPYDVMVDVIEQYSAGHLDMLSRIKSLQS
 RSCDWRGVLA

TRTRPLE – GFP Tag – V

Restriction Sites:

Sgfl-MluI

Cloning Scheme:


ACCN: NM_001006674

ORF Size: 1710 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001006674.2](#), [NP_001006675.1](#)

RefSeq Size: 2454 bp

RefSeq ORF: 1713 bp

Locus ID: 16536


UniProt ID: [Q9Z351](#)

Cytogenetics: 2 103.57 cM

Gene Summary:

Associates with KCNQ3 to form a potassium channel with essentially identical properties to the channel underlying the native M-current, a slowly activating and deactivating potassium conductance which plays a critical role in determining the subthreshold electrical excitability of neurons as well as the responsiveness to synaptic inputs. Therefore, it is important in the regulation of neuronal excitability.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG225339