

Product datasheet for **MG225403**

Cldn14 (NM_001165925) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cldn14 (NM_001165925) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Cldn14
Synonyms: AI851731
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG225403 representing NM_001165925
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCCAGCACAGCGGTCCAGCTCCTAGGCTTCCTGCTTAGCTTCCTGGGCATGGTGGGAACGCTCATCA
 CCACCATCTGCCACACTGGCGGAGGACGGCCATGTGGGCACCAACATCCTGACTGCCGTGCTACCT
 GAAGGGACTGTGGATGGAATGTGTGTGGCACAGCACAGGCATCTACCAGTGCAGATCTACCGCTCACTG
 CTGGCGCTGCCCCGGGATCTGCAGGCCGCCGGCGCTCATGGTCATCTCCTGCCTGTTGTCGGGCATGG
 CCTGCGCCTGCGCAGTAGTGGGCATGAAGTGCACACGCTGCGCCAAAGGCACACCCGCCAAGACCACCTT
 TGCAGTGTGGGGGCGCACTCTTCCTGCTGGCCGGCCTGCTGTGCATGGTGGCCGTGTCCTGGACCAGG
 AATGACGTGGTGCAGAATTTTTATAACCCGCTGCTGCCAGTGGCATGAAGTTTGAATCGGCCAGGCC
 TGTACCTGGGCTTCATCTCCTCATCCCTGTCTCTCATCGGGGACCCCTGCTCTGCTTATCCTGCCAGGA
 CGAGGCCCCCTACAGACCCTACCCGCCAGTCCAGGGCCGGAGCTACCACCAGGCTACCGCCCTGCC
 TACCGCCACCAGCGCCTACAAGGACAACCGTGCCCCCTCGGTGACCTCAGCCGCGCACAGTGGGTACA
 GGCTGAATGACTACGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG225403 representing NM_001165925
 Red=Cloning site Green=Tags(s)

MASTAVQLLGFLLSFLGMVGLITTTILPHWRRTAHVGTNILTAVSYLKGLWMECVHSTGIYQCQIYRSL
 LALPRDLQAARALMVISCLLSGMACACAVVGMKCTRCAGKTPAKTTFAVLGGALFLLAGLLCMVAVSWTT
 NDVVQNFYNPLLP SGMKFEIGQALYLGFISSSLSIGGTLCLSCQDEAPYRPPYPPQSRAGATTTATAPA
 YRPPAAAYKDNRAPSVTSAHSGYRLNDYV

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:



ACCN: NM_001165925

ORF Size: 717 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001165925.1 , NP_001159397.1
RefSeq Size:	1283 bp
RefSeq ORF:	720 bp
Locus ID:	56173
UniProt ID:	Q9Z0S3
Cytogenetics:	16 C4
Gene Summary:	This gene encodes a member of the claudin family of tight junction proteins. The encoded protein is an integral membrane protein that may function in maintaining apical membrane polarization in tight junctions located between outer hair cells and supporting cells. Loss of function of this gene is associated with hearing problems. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]