

Product datasheet for **MG225461**

Fbxw7 (NM_001177774) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fbxw7 (NM_001177774) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Fbxw7
Synonyms:	1110001A17Rik; AGO; Cdc4; Fbw7; Fbwd6; Fbx30; Fbxo30; Fbxw6; SEL-10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG225461 representing NM_001177774
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAATCAGGAAGTCTCTCTGTGGCAGCAAAGACGACGAAGTGGAGGCTCTCTGAGAGGGAATGCTT
 CCTCAAGCCAGGTTGATGAGGGACAGATGAATCGCGTGGTTGAGGAGGATCCACAGCAGCAAGCGAGACA
 TCAAGAGGAGGAGCACACTGCGCGGAATGGTGAACCTGTGGGTGCAAACCCTAGGCTGGAGACCAGAAC
 GATACCCAGCAAGGACAAGTGGAGGAAAATAATAACCGCTTTATTTTCAGTAGATGAGGACTCTTCGGGAA
 ATCAGGAAGAGCAAGAGGAAGATGAAGAGCATGCTGGGGAACAGGAGGAGGAAGAGGAGGAAGAGGAAGA
 GGAGGAGGAGATGGACCAGGAGAGTGTGATTTTGTCCGTCTGATGACAGTAGCAGAGAAGATGAACAT
 ACGCACAATAGCAATGTCACAACTGCAGTAGTGTCTCGGACCTGCCCGCTACCCAGCTCTCTCTCCAT
 TCTATACAAAGACAACAAAAATGAAAAGAAAGTTGGACCATGGTCTGAAGTTCGTTCTTTTCTTTGGG
 AAAGAAACCATGCAAAGTCTCAGATTATACCAGTACCCTGGCCTGTACCATGTTCCAGCAACCAACA
 ACTTTTGGGGACCTGAGAGCAGCCAATGGGCAAGGGCAGCAGCGCGGAGGATTACATCTGTCCAACCAC
 CCACAGGCCTTCAAGAGTGGCTGAAAATGTTTCAGAGCTGGAGCGGACCAGAGAAGTTGCTGGCTTTAGA
 TGAGCTCATTGACAGCTGTGAACCAACACAAGTGAAGCATATGATGCAAGTGATAGAGCCCCAGTCCAG
 CGAGACTTCATCTCCTTGCTTCTAAAGAGTTGGCACTCTATGTGCTTTCATTCTGGAACCCAAAGACC
 TGCTGCAAGCGGCTCAGACTTGTGATACTGGAGAATTTTGGCTGAGGATAACCTTCTCTGGAGAGAGAA
 ATGTAAGAAGAGGGGATTGATGAACCGTTGCACATCAAGAGAAGAAAAATAATAAAACAGGTTTCATA
 CACAGCCCATGGAAGAGTGCATATCAGACAGCACAGAATTGATACAACTGGAGACGAGGAGAAGTCA
 AATCTCCTAAGGTGCTGAAAGGGCATGATGACCATGTGATCACATGCCTACAGTTTGTGGCAACCGCAT
 AGTTAGTGGTTCTGATGACAACACTTTAAAAGTTTGGTCAGCGGTACCGGCAAGTGTCTGAGAACGTTA
 GTGGGACATACAGGTGGAGTGTGGTCATCACAGATGAGAGACAATATCATCATCAGTGGATCGACTGACC
 GGACTCTCAAAGTGTGGAATGCTGAACTGGAGAGTGTATACATACTTTATATGGGCACACTTCTACTGT
 ACGGTGTATGCATCTCCATGAAAAAGGGTTGTAAGCGGTTCTCGAGATGCCACTCTCAGGGTTTGGGAT
 ATTGAGACCGCCAGTGTTCACAGTCTGATGGGTACAGTAGCAGCGGTCGCTGCGTTTTCAGTATGATG
 GCAGGAGGTTGTTAGTGGAGCTTATGATTTTATGGTGAAGGTGTGGATCCAGAGACTGAGACCTGTCT
 ACACACGTTACAGGGACACACTAATAGAGTCTATTACAGTTTGTGGCATCCATGTGGTGTGAGTGA
 TCTCTTGATACATCAATCCGAGTCTGGGATGTGGAGACAGGGAATTGTATTACACGCTAACAGGACACC
 AGTCATTAACGAGTGAATGGAAGTCAAAGACAATATTCTGTCTCTGGGAATGCAGATTCTACAGTTAA
 GATCTGGGATATCAAAACAGGACAGTGTTCACAACTTTGCAAGGTCCCAGCAAGCATCAGAGCGTGTG
 ACCTGCTTACAGTTCAACAAGAAGTTCGTAATTACCAGCTCAGACGACGGAACCGTCAAACCTCTGGGACT
 TGAAAACGGGTGAATTTATCCGAAACCTCGTCACATTGGAGAGTGGGGGAGCGGGGAGTTGTGTGGCG
 GATCAGGGCCTCAAACACAAGCTGGTGTGTGCAGTCGGGAGTCGGAATGGAAGTGAAGAAACCAAGCTC
 CTGGTGTGGACTTTGATGTGGACATGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG225461 representing NM_001177774
 Red=Cloning site Green=Tags(s)

```

MNQELLSVGSKRRTGGSLRGNASSQVDEGQMNRVVEEDPQQQARHQEEHTARNGELVGANPRPGDQN
DTQQGQVEENNRFI SVDEDSSGNQEEQEEDDEEHAGEQEEEEEEEEEMDQESDDFDPSDDSSREDEH
THNSNVTNCSSVSDLPAHQ LSSPFYTKTKMKRKL DHGSEVRSFSLGKKPCKVSDYTSTTGLVPCSATPT
TFGDLRAANGQQQRRRITSVQPPTGLQEWLKMFQSWSGPEKLLALDELIDSCEPTQVKHMMQVIEPQFQ
RDFISLLPKELALYVLSFLEPKDLLQAAQTCRYWRILAEDNLLWREKCKEEGIDEPLHIKRRKI IKPGFI
HSPWKSAYIRQHRIDTNWRRGELKSPKVLKGHDDHVITCLQFCGNRIVSGSDDNTLKVWSAVTGKCLRTL
VGHTGGVSSQMRDNIIISGSDRTLKVWNAETGECIHTLYGHTSTVRCMHLHEKRVVSGSRDATLRVWD
IETGQCLHVLMGHVAAVRCVQYDGRRVVSGAYDFMVKVWDPETETCLHTLQGHTNRVYSLQFDGIHVVS
SLDTSIRVWVETGNCIHTLTGHQSLTSGMELKDNILVSGNADSTVKIWDIKTGQCLQTLQGPSKHQSAV
TCLQFNKFNVITSSDDGTVKLWDLKTGEFIRNLVTLESGSGGVVWRIRASNTKLVCAVGSRNGTEETKL
LVLD FVDVDMK
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



ACCN: NM_001177774

ORF Size: 2130 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001177774.1](#), [NP_001171245.1](#)

RefSeq Size: 4092 bp

RefSeq ORF: 2133 bp


Locus ID: 50754

UniProt ID: [Q8VBV4](#)

Cytogenetics: 3 37.7 cM

Gene Summary:

Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins (PubMed:21953459, PubMed:22748924). Recognizes and binds phosphorylated sites/phosphodegrons within target proteins and thereafter bring them to the SCF complex for ubiquitination (PubMed:22748924). Mediates ubiquitination and subsequent degradation of CCNE1 and MYC (PubMed:22748924). Identified substrates include cyclin-E (CCNE1 or CCNE2), DISC1, JUN, MYC, NOTCH1 released notch intracellular domain (NICD), NOTCH2, MCL1 and probably PSEN1 (By similarity). Acts as a negative regulator of JNK signaling by binding to phosphorylated JUN and promoting its ubiquitination and subsequent degradation (By similarity). SCF(FBXW7) complex mediates the ubiquitination and subsequent degradation of NFE2L1 (PubMed:21953459). Involved in bone homeostasis and negative regulation of osteoclast differentiation (PubMed:29149593). Regulates the amplitude of the cyclic expression of hepatic core clock genes and genes involved in lipid and glucose metabolism via ubiquitination and proteasomal degradation of their transcriptional repressor NR1D1; CDK1-dependent phosphorylation of NR1D1 is necessary for SCF(FBXW7)-mediated ubiquitination (PubMed:27238018).[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG225461