

Product datasheet for **MG225469**

Dnajc5b (NM_025489) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Dnajc5b (NM_025489) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Dnajc5b
Synonyms: 1700008A05Rik
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG225469 representing NM_025489
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCATGTAACGCACCCAACCAAGACAGCGGACTTTGTCAACATCAGGAGAAAGTCTATGAAATTC
 TTGGTCTGCATAAGGGAGCATCATGTGAAGAAATTAAGAAAACCTACAGAAAATTAGCCCTGCGACACCA
 TCCAGACAAGAATCCAGATGACCCATCTGCTGCTGAGAAGTTAAAGAAATCAACAACGCCACACAATC
 CTCACCGACACCTCAAAGAGAAACATATGACAAGTATGGATCGCTGGGACTCTACGTGGCTGAACAAT
 TTGGAGATGAAAACGTTAACACATACTTCATGTTGTCAAGCTGGTGGCAAGGACCTGTTTCATCATCAT
 TGGCCTCTTGACCGGCTGCTATTTTTGTTGCTGCCTTTGCTGCTGCAACTGTTGTTGTGGACATTGT
 CGGCCCAAATCATCGACACCAGAAGAGGAGTTCTATGTGTCAACCAGAGGACCTTGAGGAGCAGATCAGAA
 CTGACATGGAGAAAGATATGGACTTTCCAGTTGTTCTCCAGCCTACAAATACAAATGAGAAAACACAGCT
 AATCAGAGAAGGATCTCGAAGTTACTGCACAGACTCT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG225469 representing NM_025489
 Red=Cloning site Green=Tags(s)

MACNAPNQRQRTLSTSGESLYEILGLHKGASCEEIKKTYRKLALRHHPDKNPDDPSAAEKFKENNAHTI
 LTDTSKRNIYDKYGLGLYVAEQFGDENVNTYFMLS~~SS~~WWARTLFIIGLLTGCYFCCLCCCCNCCCGHC
 RPKSSTPEEEFYVSPEDLEEQIRTDMEKMDDFPVVLQPTNTNEKTLIREGSRSYCTDS


TRTRPLE - GFP Tag - V

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_025489

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_025489.2](#)

RefSeq Size: 947 bp

RefSeq ORF: 600 bp

Locus ID: 66326

UniProt ID: [Q9CQ94](#)

Cytogenetics: 3 A2