

Product datasheet for **MG225621**

Sfrp5 (NM_018780) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sfrp5 (NM_018780) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Sfrp5
Synonyms:	A1605071; SARP3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG225621 representing NM_018780 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTGGGTGGCCTGGAGCGCACGGACGGCCGCACTGGCGTTGCTGCTCGGGGCGCTGCATGGGGCGCCAA
CACGCGGCCAGGAGTACGACTACTACGGTTGGCAGGCCGAGCCGCTGCACGGCCGCTCCTACTCCAAGCC
ACCGCAGTGCCTCGACATCCCCGCCGATCTGCCGCTCTGTACACGGTGGGCTACAAGCGCATGCGGCTG
CCCAACCTGCTGGAGCACGAGAGCCTGGCCGAGGTGAAGCAGCAGGCAAGCAGCTGGCTGCCACTGCTGG
CCAAGCGTGGCACTCAGACACCCAGGTCTTCCTCTGCTCGCTCTCGCTCCCGTCTGCCTGGACCGACC
CATCTACCCCTGCCGCTCGCTGTGCAAGCCGTGCGCGCCGGCTGCGCTCCGCTCATGGAGGCCTACGGT
TTCCCTTGGCCCGAGATGCTGCACTGCCACAAGTTCGCCCTGGACAACGACCTCTGCATCGCGGTGCAGT
TCGGGCACCTGCCTGCCACCGCGCTCCAGTGACCAAGATCTGTGCCAGTGTGAGATGGAGCACAGCGC
TGATGGCCTCATGGAACAGATGTGCTCCAGTGACTTTGTGGTCAAGATGCGCATTAAAGGAGATCAAGATA
GACAACGGGGACCGAAAGTTGATTGGAGCCCAGAAGAAGAAGCTGCTCAAGGCAGGCCCTTAAAGC
GCAAGGACACCAAGAAGCTGGTCTGCATATGAAGAACGGGGCAAGCTGCCCTTGCCACAGTTAGACAA
CCTGACGGGCAGCTTCTGGTTCATGGCCGCAAAGTAGAGGGACAGCTGCTGCTCACGGCCGTCTACCGC
TGGGACAAGAAGAATAAGGAGATGAAGTTGCAAGTCAAATTCATGTTCTCCTATCCCTGTTCCCTCTACT
ACCCTTTTTCTATGGGGCAGCTGAACCCAC

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG225621 representing NM_018780
 Red=Cloning site Green=Tags(s)

MWVAWSARTAALALLL GALHGAPTRGQEYDYYGWQAEPLHGRSYSKPPQCLDIPADLPLCHTVGYKRMRL
 PNLL EHESLAEVKQQASSWLPLLAKRCHSDTQVFLCSLFAPVCLDRPIYPCRSLCEAVRAGCAPLMEAYG
 FPWPEMLHCHKFPLDNDLCIAVQFGHLPATAPPVTKICAQCEMEHSADGLMEQMCSDFVVKMRIKEIKI
 DNGDRKLI GAQKKKKLLKAGPLKRKDTKKLVLHMKNGASCPCQLDNL TGSFLVMGRKVEGQLLLTAVYR
 WDKKNKEMKFAVKFMFSYPCSLYYPFFYGAAEPH

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_018780

ORF Size: 942 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018780.3](#), [NP_061250.2](#)

RefSeq Size: 1900 bp

RefSeq ORF: 945 bp

Locus ID: 54612

UniProt ID: [Q9WU66](#)

Cytogenetics: 19 C3

Gene Summary: Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types. SFRP5 may be involved in determining the polarity of photoreceptor, and perhaps, other cells in the retina.[UniProtKB/Swiss-Prot Function]