

Product datasheet for **MG225926**

Rbp4 (NM_011255) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rbp4 (NM_011255) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Rbp4
Synonyms: Rbp-4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG225926 representing NM_011255
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTGGGTGTGGGCGCTCGTGCTGCTGGCGGCTCTGGGAGGCGGCAGCGCCGAGCGCGACTGCAGGG
TGAGCAGCTTCCGAGTCAAGGAGAACTTCGACAAGGCTCGTTTCTCTGGGCTCTGGTATGCCATCGCCAA
AAAGGACCCCGAGGGTCTCTTTTGAAGACAACATCATCGCTGAGTTTTCTGTGGACGAGAAGGGTCAT
ATGAGCGCCACAGCCAAGGACGAGTCCGTCTTCTGAGCAACTGGGAAGTGTGTGCAGACATGGTGGCA
CTTTCACAGACTGAAGATCCTGCCAAGTTCAAGATGAAGTACTGGGGTGTAGCCTCCTTCTCCAGCG
AGGAAACGATGACCACTGGATCATCGACACGGACTACGACACCTTCGCTCTGCAGTACTCCTGCCGCTG
CAGAATCTGGATGGCACCTGTGCAGACAGCTACTCCTTTGTGTTTTCTCGTGACCCCAATGGCCTGAGCC
CAGAGACACGGAGGCTGGTGAAGCAACGGCAGGAGGAGCTGTGCCTAGAGAGGCAGTACAGATGGATCGA
ACACAATGGTTACTGTCAAAGCAGGCCCTCCAGAAACAGTTTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG225926 representing NM_011255
Red=Cloning site Green=Tags(s)

MEWVWALVLLAALGGGSAERDCRVSSFRVKENFDKARFSGLWYAIKKDPEGLFLQDNIIAEFSVDEKGH
MSATAKGRVRLLSNWEVCADMVGTFTDTEPAKFKMKYWGVAFLQRGNDDHWIIDTDYDTFALQYSCRL
QNLDGTCADSYSFVSRDPNGLSPETRRLVQRQEELCLERQYRWIEHNGYCQSRPSRNSL


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_011255

ORF Size: 603 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_011255.2](#)

RefSeq Size: 930 bp

RefSeq ORF: 606 bp

Locus ID: 19662

UniProt ID: [Q00724](#)

Cytogenetics: 19 32.75 cM

Gene Summary: Retinol-binding protein that mediates retinol transport in blood plasma. Delivers retinol from the liver stores to the peripheral tissues. Transfers the bound all-trans retinol to STRA6, that then facilitates retinol transport across the cell membrane.[UniProtKB/Swiss-Prot Function]