

Product datasheet for **MG226194**

Pnpt1 (NM_027869) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pnpt1 (NM_027869) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Pnpt1
Synonyms:	1200003F12Rik; Old35; PNPase; Pnptl1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG226194 representing NM_027869
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGCTGCAGGCTGTGCTGCTTGTGCCCGTGCCTCCGACCGCTGGGCTGCGGCCCCCTCGGCCGGC
 CCGGGCGGAATCGGGCACTCAGCTATTTGCAGATGCGAGCGTTGTGGAGCAGCACCGGATCCAGAGCGGT
 GACTGTAGACCTGGGACACAGAAAATTAGAAAATCCTCTGGGAAACTGGCAAGATTTGCTGATGGCTGT
 GCTGTGATACAGTCAGGCGATACTGCAGTAATGGTCACGGCAGTCAGAAAAAAAAGCTTCGCCTTCCC
 AATTCATGCCGTTGGTGGTGGACTACCGACAGAAGGCTGCTGCAGCAGGCAGAATCCCCACAACTACCT
 TAGAAGGGAGATTGGCTCCTCTGACAGAGAGGTTCTTACAAGTCGAGTAATAGATCGTTCAATTCGCCCT
 CTCTTTCCAGCCGGCTATTTTTATGATACTCAGGTTCTCTGTAATCTGCTAGCAGTAGATGGCATCAATG
 AACCTGATATCCTAGCAGTTAATGGTCTTCTGTAGCCCTCTCCTTATCAGATATCCCTTGAATGGACC
 TGTTGGGGCAGTACGAATAGGAATGATTGATGGAGAATGTGTCGTTAACCCAACAAGGAGAGAAATGTCT
 TCTAGCACTTTAAATTTAGTAGTTGCCGGAGCACCTAAAAGCCAAATTTGTTATGTTGGAAGCCTCTCGAG
 AAAATATTCTACAGCAGGACTTTTGCCATGCTATCAAAGTTGGGGTGAAGTATACACAGCAGATAATTCA
 GGGCATCCAGCAGTTGGTAAAAGAAATCGGTGTTGCCAAGAGGACACCGCAGAAGATATTTACTCCTTCT
 GCAGAGATTGTGAAGTACACGAAGATAATTGCCATGGAGAACTCTATGCGGTTTTTACAGATTATGAAC
 ATGATAAAGTTTTCCAGGGATGAAGCTGTTAACAAGATAAGATTAGATACAGAGGAGCATCAAAGAAAA
 ATTTCCAGAGGTTGACCAATTTGAAATAATAGAATCCTTCAACATTGTTGCAAAGGAGGTTTTCCGAAGT
 ATATTTTTGAATGAATACAAAAGGTGTGATGGAAGAGATCTGACTTCACTTAGGAATATAAGTTGTGAGG
 TCGATATGTTTTAAAACACTTCATGGATCAGCATTATTTCCAGAGAGGACAAACACAGGTACTGTACTGT
 TACATTTGATTCATTAGAATCCAGTATTAAGTCAGATCAAATTATAACAGCTATAAATGGGGTAAAAGAT
 AAAAATTTTCATGCTGCACTATGAGTTCCCTCCTTATGCAACCAATGAAACCGGCAAAAGTTACTGGTGTGA
 ATCGAAGAGAACTTGACATGGTCTCTTGCCGAGAAAGCTTTGTGCTCTGTTATTTCCAAAGATTTTCC
 TTTTACCATAAGAGTTACATCTGAAGTCTCGAATCAAATGGGTCATCTTCTATGGCATCTGCATGTGGT
 GGAAGTTTGGCATTAAATGGATGCAGGGTCCCAATTTCTGCTGTTGCAGGTGTAGCAGTGGGATTGG
 TTACCAAAACCAATCCAGAGAAAGGTGAAATAGAAGACTACCGTTTGCTAACAGATATTCTGGGAATTGA
 AGATTATAATGGTGACATGGATTTCAAATAGCCGGTACAATAAAGGAATAACTGCATTACAGGCTGAT
 ATTAAGTTACCTGGAGTACCAATTAATAATGGAAGCCATCCAACAAGCGTCAGTGGCAAAGAAGG
 AGATACTGCAGATAATGAACAAAACGATTTCAAACCTCGAGCATCAAGAAAAGAAAATGGACCAGTTGT
 AGAAACAGTAAAGGTTCCATTATCAAAACGAGCAAAATTCGTTGGGCCTGGTGGATCACTTAAAAAAA
 CTCCAGGCTGAGACAGGTGTAACAATTAGTCAGGTTGATGAAGAAACCTTCTCCATATTTGCACCAACAC
 CTACTGCAATGCATGAAGCAAGAGATTTTATTACAGAAATTTGCAGAGATGATCAAGAGCAACAATTAGA
 ATTTGGAGCAGTTTATACCGCGACAATAACTGAAATCAGAGACTGGAGTGTGGTAAAACCTGTATCCA
 AACATGACTGCAGTGTGCTTCATAATTCACAACCTTGACCAACGAAAGATTAAACATCCCCTGCCCCTAG
 GACTAGAGGTTGGCCAAGAAATTCAGGTCAAACTTTGGCCGTGATCCAGCTGATGGAAGAATGAGGCT
 TTCTCGTAAAGTACTTCAGTCTCCAGCTACAACCTGCTCTCAAACTCTAAATGATAGAAGCAGCATTGTA
 ATGGGAGAGCCTGTCTCACAGTCATCTAACTCTAACCTT

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG226194 representing NM_027869
 Red=Cloning site Green=Tags(s)


```

MAACRLCCLCPCLRPLGCGPLGRPGRNRALSYLQMRALWSSTGSRVAVTDLGHRKLEISSGKLARFADGC
AVIQSGDTAVMVTAVSKTKASPSQFMPLVVDYRQKAAAAGRIPTNYLRREIGSSDREVLTSRVIDRSIRP
LFPAGYFYDTQVLCNLLAVDGINEPDILAVNGASVALSLSDIPWNGPVGAVRIGMIDGECVVNPTRREMS
SSTLNLVVAGAPKSIIVMLEASAENILQQDFCHAIKVGKVTQQIIQGIQQLVKEIGVAKRTPQKIFTPS
AEIVKYTKIIAMEKLYAVFTDYEHDKVSRDEAVNKIRLDTEEHLKEKFPEVDQFEIIESFNIVAKEVFRS
IILNEYKRCDDGRDLTSLRNISCEVDMFKTLHGSALEFQRGQTQVLCVTFDSLESSIKSDQIITAINGVKD
KNFMLHYEFPPYATNETGKVTGVNRRELGHGALAEKALCPVIPKDFPFTIRVTSEVLESNGSSMASACG
GSLALMDAGVPISSAVAGVAVGLVTKTNPEKGEIEDYRLLTDILGIEDYNGDMDFKIAGTNKGITALQAD
IKLPGVPIKIIIMEAIQQASVAKKEILQIMNKTISKPRASRKENGPPVETVKVPLSKRAKRVGPGGYHLKK
LQAETGVTISQVDEETF SIFAPTPTAMHEARDFITEICRDDQEQQLFEGAVYTATITEIRD TGVMVKLYP
NMTAVLLHNSQLDQRKIKHPTALGLEVGQEIQVKYFGRDPADGRMRLSRKVLQSPATTALKTLNDRSSIV
MGEPVVSQSSNSNP
  
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_027869


ORF Size: 2349 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_027869.2
RefSeq Size:	2704 bp
RefSeq ORF:	2352 bp
Locus ID:	71701
UniProt ID:	Q8K1R3
Cytogenetics:	11 A3.3
Gene Summary:	<p>RNA-binding protein implicated in numerous RNA metabolic processes. Catalyzes the phosphorolysis of single-stranded polyribonucleotides processively in the 3'-to-5' direction. Mitochondrial intermembrane factor with RNA-processing exoribonuclease activity. Component of the mitochondrial degradosome (mtEXO) complex, that degrades 3' overhang double-stranded RNA with a 3'-to-5' directionality in an ATP-dependent manner. Required for correct processing and polyadenylation of mitochondrial mRNAs. Plays a role as a cytoplasmic RNA import factor that mediates the translocation of small RNA components, like the 5S RNA, the RNA subunit of ribonuclease P and the mitochondrial RNA-processing (MRP) RNA, into the mitochondrial matrix. Plays a role in mitochondrial morphogenesis and respiration; regulates the expression of the electron transport chain (ETC) components at the mRNA and protein levels. In the cytoplasm, shows a 3'-to-5' exoribonuclease mediating mRNA degradation activity; degrades c-myc mRNA upon treatment with IFNβ/IFN-beta, resulting in a growth arrest in melanoma cells. Regulates the stability of specific mature miRNAs in melanoma cells; specifically and selectively degrades miR-221, preferentially. Plays also a role in RNA cell surveillance by cleaning up oxidized RNAs. Binds to the RNA subunit of ribonuclease P, MRP RNA and miR-221 microRNA.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MG226194