



**Protein Sequence:** >MG226314 representing NM\_001033405  
 Red=Cloning site Green=Tags(s)

MEPWPLFLLLLLLLLLWLQGCVSGHSNENLYRKVWRREGETLSVQCSYKNRRNLVEAKSWCKVKKKKCDH  
 NFTRSWVRGPGSYSLRDDAKVKVVRITMEALRVQDSGRYWCMRNTAGHFYPLVGFQLEVPALTTERNVPH  
 THLTNTPMDGFVTTGQVHISDPHAPFTSDVTMFTSEVTMFTSGLLTLASGTTTPTPVTGYSFIDTSGTVT  
 EPERNTESQPATLSPSNARSFADPVTSTMSRHQSSSLSTTGCHPLTPNRSQETYIPAMVVVLTFLPA  
 PVVLLVAYGFWKKRHRMGRYNLGSNYAKPWIHLPEGPETPWKPAWSKITQ

TRTRPLE - GFP Tag - V


**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



**Plasmid Map:**



**ACCN:** NM\_001033405

**ORF Size:** 987 bp

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001033405.2</a> , <a href="#">NP_001028577.2</a>
<b>RefSeq Size:</b>	4509 bp
<b>RefSeq ORF:</b>	990 bp
<b>Locus ID:</b>	328833
<b>UniProt ID:</b>	<a href="#">Q2LA85</a>
<b>Cytogenetics:</b>	17 C
<b>Gene Summary:</b>	Cell surface receptor that may play a role in the innate and adaptive immune response. Acts as a counter-receptor for CD276 and interaction with CD276 on T-cells enhances T-cell activation (By similarity).[UniProtKB/Swiss-Prot Function]