

Product datasheet for **MG226317**

Chka (NM_013490) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|-----------------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | Chka (NM_013490) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Chka |
| Synonyms: | Chetk-alpha; Chk; ChoK; CK; CK/EK-alpha; EK; EtnK-alpha |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |
| ORF Nucleotide Sequence: | >MG226317 representing NM_013490 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAACCAAGTTCTGCACCGGGGGCGAGGCCGAGCCGTCGCCTCTCGGGCTGCTGCTGAGCTGCGGTG
GCAACGCTGCCCCGACGCCCGCGTAGGGCAGCAGCGGGACGCCGAGCGGAGCTGGAGTCCAAGCAGCT
TGGTGGCCGAACCAACCTCTCGCGCTGCCGCCGCCGCCGCCGCCCTGCCGCTGCCCCGCCGCCA
TCACCGCCGCTAGCGGACGAACAACCCGAGCCCCGGACCAGGCGCAGGGCCTACCTGTGGTGAAGGAAT
TCCTGCCCGGAGCCTGGCGAGGCCTTCGCGAGGACCAGTCCACATCAGTGTATCAGGGGTGGTCTCAG
TAACATGCTGTTCCAGTCTCCCTGCCAGACTCCATAGCCAGCGTTGGCGATGAGCCTCGGAAAGTGCTC
TTGCGGCTCTATGGGGCAATCTTAAAGATGGGGGCTGAGGCGATGGTTCTGGAGAGTGTTATGTTTGCCA
TTCTTGACAGAGAGTCACTTGGGCCAAAACCTTTTGGCATCTTTCCCAAGGCCGACTGGAGCAGTTTAT
CCCAGCCGGCGATTGGACACTGAAGAGTTACGTTTACCAGATATTTCTGCAGAAATAGCTGAGAAAAATG
GCCACATTTATGGTATGAAAATGCCATTCAATAAGGAACCAAAATGGCTTTTGGAAACATGAAAAAT
ACCTGAATCAAGTACTAAGACTCAAATTCAGCAGGGAGGCCAGAGTTCAGCAACTGCACAAGATCCCTCTC
TTACAACCTGCCTCTTGAGCTCGAGAACCTGAGGTCATTGCTGCAGTATACTAGATCTCCAGTTGTATTT
TGTCATAATGACTGTCAAGAAGTAATATCTTATTGTTGGAAGGCCAGGAGAATCTGAAAGGCCGGAAGC
TGATGCTTATTGACTTTGAGTACAGCAGTTACAATTACAGGGGATTTGACATTGGAATCATTTCTGTGA
ATGGATGTATGATTATACCTATGAGAAGTATCCTTTCTCAGAGCGAACATTCAGAAGTATCCTTCCCGA
AAACAACAGCTCCATTTTATTTCCAGTTACTTGACTACATTCCAAAATGATTTTGAAGCCTCAGTAGTG
AAGAGCAGTTTGTACAAAAGAAGACATGTTGCTTGAAGTCAACAGATTTGCTCTGCCTCTCATTTCTC
CTGGGGACTTTGGTCCATTGTACAGGCCAAGATCTCATCCATTGAATTTGGGTACATGGAATATGCCCAA
GCCAGATTCGAGGCTACTTTGACCAGAAGAGGAAGCTTGGGGTG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online »](#)

Protein Sequence: >MG226317 representing NM_013490
 Red=Cloning site Green=Tags(s)

```


MKTKFCTGGEAEPSPLGLLLSCGGNAAPTPGVGQQRDAAGELESKQLGGRTQPLALPPPPPPPLPLPPP
SPPLADEQPEPRTRRRAYLWCKEFLPGAWRGLREDQFHISVIRGGLSNMLFQCCLPDSIASVGDPRKVL
LRLYGAILKMGAEAMVLESVMFAILAERSLGPFLGIFPQGRLEQFIPSRRLDTEELRLPDISAEIAEKM
ATFHGMKMPFNKEPKWLFGTMEKYLNQVLRKFSREARVQQLHKILSYNLPLELENLRSLLQYTRSPVVF
CHNDQCQEGNILLLEGQENSERRKMLIDFEYSSYNYRGFDIGNHFCEWMDYTYEKYPFFRANIQKYPVSR
KQQLHFISSYLTFQNDFESLSSEEQFATKEDMLLEVNRFALASHFLWGLWSIVQAKISSIEFGYMEYQAQ
ARFEAYFDQKRKLGV
  
```

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM_013490

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ORF Size: | 1305 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_013490.4 , NP_038518.2 |
| RefSeq Size: | 2589 bp |
| RefSeq ORF: | 1308 bp |
| Locus ID: | 12660 |
| UniProt ID: | O54804 |
| Cytogenetics: | 19 3.57 cM |
| Gene Summary: | Has a key role in phospholipid biosynthesis and may contribute to tumor cell growth. Catalyzes the first step in phosphatidylcholine biosynthesis. Contributes to phosphatidylethanolamine biosynthesis. Phosphorylates choline and ethanolamine. Has higher activity with choline (By similarity).[UniProtKB/Swiss-Prot Function] |