

Product datasheet for **MG226419**

Dhx36 (NM_028136) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Dhx36 (NM_028136) Mouse Tagged ORF Clone |
| Tag: | TurboGFP |
| Symbol: | Dhx36 |
| Synonyms: | 2810407E23Rik; AI452301; AU022184; Ddx36; mKIAA1488 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-AC-GFP (PS100010) |
| E. coli Selection: | Ampicillin (100 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>MG226419 representing NM_028136
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAGCTACGACTATCATCAGAGCTGGAGCCGCGATGGGGGCCACGGGGCTCCGGCCAGGGCTCCAGCG
 GCGGCGGCGGAGGAGCCGGGGCTCGGGCGGCGGAGGGGGCCGCGGGCGGGCCGGCATCCCGC
 ACACCTTAAGGGTCGCGAGATCGGCCTGTGGTACGCCAAGAAGCAGACGCAGAAGAACAAGGAGGCTGAG
 AGGCAGGAGAGAGCTGTAGTGCACATGGATGAACGTGCGAGAAGAGCAAATGTGCAGCTGCTGAATTCAG
 TCCAAGCTAAGACTGACAAAGATTCAGAAGCGCAGATATCCTGGTTTGTCTCTGAGGATCATGGGTATGG
 TACTGAAGTTTCTTCAGAGAAAAAATAAACTCAGAGAAGAACTTGACAACCAGGAAAAGAAATGCTA
 AACCAAGAAAAAAGACATTTAGGATCACAGACAAATCATATATTGACCGAGATACTGAGATTTTATTGC
 AAGAAAAAGAGCCAAACCTAAGCTTAGATCAACATTTACTGGAAGATTTACAAGGAAAAAACTGACCC
 TCGATATATAGAGATGCAGCGTTTCAGAAAAAGCTGCCTTCATATGGAATGCAGAAGGAGCTGGTAAAT
 CTAATCAATAACCATCAGGTGACAGTAATAAGTGGTGAAGCTGGTTGTGGCAAAACCACTCAGGTTACGC
 AGTTCATCTTGGATAACTACATCGAAAGAGGAAAAGGGTCTGCCTGCAGAATAGTGTGACTCAGCCAAG
 AAGAATTAGTGCCATTTAGTGCCTGAGAGAGTGGCCACAGAAAGGGCAGAGTCTTGTGGCAATGGTAAT
 AGTACTGGATACCAGATTCGTCTCAAAGTCGGTTGCCAAGGAAACAAGGTTCTATCTTATACTGCACAA
 CAGGAATCATTCTCAGTGGCTCCAGTCCGACTCACGTTTGTCCAGTGTAGTCATATTGTACTTGATGA
 AATTCATGAAAGGAATCTACAGTCGGATGTTTTAATGACTGTTATTAAGATCTTCTCCATTTTCGATCT
 GATCTCAAAGTAATATTGATGAGTGAACCTTGAATGCTGAGAAATTTTCAAGATATCTTTTGAAGATATC
 CAATGATACATATACCCGGGTTACTTTTCCAGTTGTGGAATATCTTTTGAAGATATCATTGAAAAAT
 AAGATATGTTCCAGACCAAAAAAGAACAATAGATCCAGTTCAGAGGGGTTTCATGCAGGGTCATGTAAT
 AGACAAGAAAAAGAAGAAAAAGAGGCCATCTATAAGGAACGCTGGCCAGCGTATATAAGGAACCTGCCGA
 CAAGATACTCTGCAAGTACCGTAGATGTTTTGCAATGATGGATGATGATAAAGTTGATCTGAATTTGAT
 TGCTGCCCTTATTCGATACATTGTTTTGGAAGAAGAGGATGGTGCATATTGGTCTTTCTACCAGGCTGG
 GACAATACAGTACTTTGCATGATCTCTGATGTCACAAGTGTGTTAAATCAGATAAGTTTCTCATT
 TACCTTTACATTCAGTATGCCTACCGTCAACCAGACACAGGATTTAAAAAACTCCTCCCGGTGTTCCG
 GAAAATAGTAATTGCTACCAACATTGCAGAGACTAGCATACCATAGATGATGTGGTTTATGTAATAGAT
 GGAGGAAAAATTAAGAAACACACTTTGACACACAGAACAACATCAGTACCATGTCTGCTGAGTGGGTTA
 GTAAGCTAATGCCAAACAGAGGAAAGGAGAGCTGGAAGAGTTCAGCCTGGTCATTGTTATCATCTGTA
 TAATGGTCTTAGAGCAAGTCTTCTAGATGACTACCAACTACCAGAAATTTGAGAAGTCTTTGGAAGAA
 CTCTGTTTGCAAATAAGATCTTGAGGCTGGGTGGAATTGCTTATTTTCTGAGTAGATTAATGGATCCAC
 CATCAATGAGGCAGTGGTGTCTCCATAAAACACCTAATGGAAGTGAAGTCTTTGGATAAGCAAGAAGA
 ATTAACACCTCTTGGTGTCCATTTAGCCCGACTGCCTGTTGAGCCACATATTGAAAAATGATTCTCTTT
 GGAGCTTTGTTCTGCTGCTTAGATCCAGTTCTCACCATTGCTGCCAGTCTCAGCTTTAAGATCCCTTTG
 TCATTCCATTGGGAAAAGAAAAGATTGCAGATGCAAGAAGAAAAGAACTGGCAAAGGAACTAGAAGTGA
 TCACCTGACAGTTGTAATGCATTTGAGGGCTGGGAAGAAGCCAAACGACGTGGTTTCAGGTATGAAAAA
 GACTACTGCTGGGAATATTTTCTGCTTCAAACACACTACAGATGCTGCATAACATGAAGGGACAGTTTG
 CTGAGCATCTTCTGGAGCTGGATTTGTAAAGCAGTAGAAGTCCCAAAGATCCAAAAGCTAATATAAATTC
 AGATAATGAGAAGATAATTAAGCTGTCTGTGCTGGTTTATATCCCAAAGTTGCTAAAAATCCGACTA
 AATTTGGGTA AAAAGAGGAAAAATGGTAAAAGTTTACACGAAGTCTGACGGCCTGGTTTCTATTATCCTA
 AGTCTGTCAACGTGGAGCAGACAGACTTCCACTACAAGTGGCTTATCTATCACCTGAAGATGAGAACAAG
 CAGTATACTTTGTACGACTGTACAGAAGTGTACCATACTGCCTCCTGTTCTTTGGAGGAGATATTTCC
 ATCCAGAAAGATAAGGATCAGGAAATATTGCTGTAGATGAGTGGATTGTGTTTCAGTCTCCAGAAAGAA
 TTGCCATCTTGTAAAGGACTAAGAAAGAACTGGATAGCCTTCTACAAGAGAAGATTGAAAGCCCTCA
 TCCTGTAGACTGGGACGACACTAAGTCAAGAGACTGTGCGACTGTGCTAGCTATTCTAGACTTGATCAAA
 ACTCAAGAAAAGGCTACTCCAAGAACTTGCCACCACGGTCACAGGACGGATATTATAGC

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG226419 representing NM_028136
 Red=Cloning site Green=Tags(s)


MSYDYHQSWSRDGGPRGSGQGSSGGGGSRGSGGGGGGRGRGRHPAHLKGREIGLWYAKQTQKNKEAE
 RQERAVVHMDERREEQIVQLLNSVQAKTDKDSEAQISWFAPEDHGYGTEVSSEKKINSEKLDNQEKKLL
 NQEKKTRFRITDKSYIDRDEYLLQENEPNLSLDQHLLLEDLQRKKTDPRIEMQRFRRKLLPSYGMQKELVN
 LINNHQVTVISGETGCGKTTQVTQFILDNYIERGKGSACRIVCTQPRRISAI SVAERVATERAESCGNGN
 STGYQIRLQSRLPRKQGSILYCTTGII LQWLQSDSRLSSVSHIVLDEIHERNLQSDVLMTVIKDLLHFRS
 DLKVILMSATLNAEFSEYFGNCPMIHIPGFTFPVVEYLLEDIIEKIRYVPDQKEHRSQFKRGMQGHVN
 RQEKKEKEAIYKERWPAYIKELRTRYASTVDVLQMMDDDKVDLNLIAALIRYIVLEEDGAILVFLPGW
 DNISTLHDLMSQVMFKSDKFLI IPLHSLMPTVNQTQVFKKTPPGVRKIVIAETSIITDDVVYVID
 GGIKIKETHFDTQNNISTMSAEWVSKANAKQRKGRAGRVPQGHYHLYNGLRASLLDDYQLPEILRTPLEE
 LCLQIKILRLGGIAYFLSRLMDPPSNEAVVLSIKHLMEL SALDKQEEL TPLGVHLARLPVEPHIGKMILF
 GALFCCLDPVLTIAASL SFKDPFVIPLGKEKIADARRKELAKETRSDDL TVVNAFEGWEEAKRRGFRIEK
 DYCWEIFLSSNTLQMLHNMKGQFAEHL L GAGFVSSRSPKDPKANINSDNEKI I KAVICAGLYPKVAKIRL
 NLGKKRKMVKVHTKSDGLVSIHPKSVNVEQTDHFYHNL IYHLKMRSSIYL YDCTEVSPYCLLFFGGDIS
 IQDKDKQEI IAVDEWIVFQSPERIAHLVKGLRKELD SLLQEKIESPHVPDWD DTKSRDCAVLSAILDLIK
 TQEKATPRNLPPRSQDGYYS

TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_028136

ORF Size: 3003 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_028136.2](#), [NP_082412.2](#)

RefSeq Size: 4975 bp

RefSeq ORF: 3006 bp

Locus ID: 72162

UniProt ID: [Q8VHK9](#)

Cytogenetics: 3 E1

Gene Summary:

Multifunctional ATP-dependent helicase that unwinds G-quadruplex (G4) structures (PubMed:25611385). Plays a role in many biological processes such as genomic integrity, gene expression regulations and as a sensor to initiate antiviral responses (PubMed:21703541, PubMed:21590736). G4 structures correspond to helical structures containing guanine tetrads (By similarity). Binds with high affinity to and unwinds G4 structures that are formed in nucleic acids (G4-ADN and G4-RNA) (By similarity). Plays a role in genomic integrity (By similarity). Converts the G4-RNA structure present in telomerase RNA template component (TREC) into a double-stranded RNA to promote P1 helix formation that acts as a template boundary ensuring accurate reverse transcription (By similarity). Plays a role in transcriptional regulation. Resolves G4-DNA structures in promoters of genes, such as YY1, KIT/c-kit and ALPL and positively regulates their expression (PubMed:25611385) (By similarity). Plays a role in post-transcriptional regulation (By similarity). Unwinds a G4-RNA structure located in the 3' UTR polyadenylation site of the pre-mRNA TP53 and stimulates TP53 pre-mRNA 3'-end processing in response to ultraviolet (UV)-induced DNA damage (By similarity). Binds to the precursor-microRNA-134 (pre-miR-134) terminal loop and regulates its transport into the synapto-dendritic compartment (By similarity). Involved in the pre-miR-134-dependent inhibition of target gene expression and the control of dendritic spine size (By similarity). Plays a role in the regulation of cytoplasmic mRNA translation and mRNA stability (By similarity). Binds to both G4-RNA structures and alternative non-quadruplex-forming sequence within the 3' UTR of the PITX1 mRNA regulating negatively PITX1 protein expression (By similarity). Binds to both G4-RNA structure in the 5'-UTR and AU-rich elements (AREs) localized in the 3' UTR of NKX2-5 mRNA to either stimulate protein translation or induce mRNA decay in an ELAVL1-dependent manner, respectively (By similarity). Binds also to ARE sequences present in several mRNAs mediating exosome-mediated 3'-5' mRNA degradation (By similarity). Involved in cytoplasmic urokinase-type plasminogen activator (uPA) mRNA decay (By similarity). Component of a multi-helicase-TICAM1 complex that acts as a cytoplasmic sensor of viral double-stranded RNA (dsRNA) and plays a role in the activation of a cascade of antiviral responses including the induction of proinflammatory cytokines via the adapter molecule TICAM1 (PubMed:21703541). Required for the early embryonic development and hematopoiesis (PubMed:22422825). Involved in the regulation of cardioblast differentiation and proliferation during heart development (PubMed:26489465). Involved in spermatogonia differentiation (PubMed:25611385). May play a role in ossification (PubMed:21590736). [UniProtKB/Swiss-Prot Function]