

Product datasheet for **MG226528**

Il1rn (NM_001039701) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Il1rn (NM_001039701) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: Il1rn
Synonyms: F630041P17Rik; IL-1ra
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG226528 representing NM_001039701
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAAATCTGCTGGGGACCTACAGTCACCTAATCTCTCTCTCTCATCCTTCTGTTTCATTAGAGG
 CAGCCTGCCGCCCTTCTGGGAAAAGACCCTGCAAGATGCAAGCCTTCAGAATCTGGGATACTAACCAGAA
 GACCTTTACCTGAGAAACAACCAGCTCATTGCTGGGTACTTACAAGACCAAATATCAAACCTAGAAGAA
 AAGATAGACATGGTGCCTATTGACCTTCATAGTGTGTTCTTGGGCATCCACGGGGCAAGCTGTGCCTGT
 CTTGTGCCAAGTCTGGAGATGATATCAAGCTCCAGCTGGAGGAAGTTAACATCACTGATCTGAGCAAGAA
 CAAAGAAGAAGACAAGCGCTTACCTTCATCCGCTCTGAGAAAGGCCCCACCACCAGCTTTGAGTCAGCT
 GCCTGTCCAGGATGGTTCCTCTGCACAACACTAGAGGCTGACCGTCCCTGTGAGCCTACCAACACACCGG
 AAGAGCCCCTTATAGTCACGAAGTTCTACTTCCAGGAAGACCAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG226528 representing NM_001039701
 Red=Cloning site Green=Tags(s)

MEICWGPYSHLISLLLILLFHSEAACRPSGKRPKMQAFRIWDTNQTFFYLRRNQLIAGYLQGPNIKLEE
 KIDMVPIDLHSVFLGIHGGKLCCLSCAKSGDDIKLQLEEVNITDL SKNKEEDKRFTFIRSEKGPPTTSFESA
 ACPGWFLCTTLEADRPVSLTNTPEEPLIVTKFYFQEDQ

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_001039701

ORF Size: 534 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001039701.3</u> , <u>NP_001034790.1</u>
RefSeq Size:	2480 bp
RefSeq ORF:	537 bp
Locus ID:	16181
UniProt ID:	<u>P25085</u>
Cytogenetics:	2 16.36 cM
Gene Summary:	Inhibits the activity of interleukin-1 by binding to receptor IL1R1 and preventing its association with the coreceptor IL1RAP for signaling. Has no interleukin-1 like activity. [UniProtKB/Swiss-Prot Function]