

## Product datasheet for **MG226812**

### **Rbm3 (NM\_001166411) Mouse Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rbm3 (NM\_001166411) Mouse Tagged ORF Clone  
**Tag:** TurboGFP  
**Symbol:** Rbm3  
**Synonyms:** 2600016C11Rik  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-AC-GFP (PS100010)  
**E. coli Selection:** Ampicillin (100 ug/mL)  
**ORF Nucleotide Sequence:** >MG226812 representing NM\_001166411  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGTCGTCTGAAGAAGGGAAACTCTTCGTAGGAGGGCTCAACTTCAACACCGATGAACAGGCACTTGAAG  
ACCACTTCAGCAGCTTTGGGCCTATCTCTGAGGTGGTTGTTGTCAAGGACCGGGAGACTCAAAGATCCCG  
GGGTTTTGGCTTCATCACCTTCACAAACCCAGAGCATGCCTCAGATGCGATGAGAGCTATGAATGGAGAG  
TCCCTGGATGGGCGCAAATCCGAGTTGATCATGCAGGAAAGTCTGCCAGGGGATCCAGAGGGGGTGCCT  
TTGGTGGCGTGGTTCGAGTTACTCTAGAGGTGGTGGAGACCAGGATATGGAAGTGAAGATATGACAG  
TCGTCTGGAGGATATGGATATGGGTATGGCGGTCTAGAGACTACAGTGGCAGCCAGGGTGGCTATGAC  
CGCTACTCAGGAGAAATTACAGAGACAATTATGACAAC

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

**Protein Sequence:** >MG226812 representing NM\_001166411  
Red=Cloning site Green=Tags(s)

MSSEEGKLFVGGILNFNTDEQALEDHFSSFGPISEVVVKDRETQRSRGFFITFTNPEHASDAMRAMNGE  
SLDGRQIRVDHAGKSARGSRGGAFGGRGRSRSRGGDQGYGSGRYDSRPGGYGYGYSRDRYSGSQGGYD  
RYSGGNYRDNYDN


**TRTRPLE** - GFP Tag - V

**Restriction Sites:** Sgfl-MluI



[View online »](#)

**Cloning Scheme:**

**Plasmid Map:**


ACCN: NM\_001166411

ORF Size: 459 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001166411.2</u>
<b>RefSeq Size:</b>	1267 bp
<b>RefSeq ORF:</b>	462 bp
<b>Locus ID:</b>	19652
<b>UniProt ID:</b>	<u>O89086</u>
<b>Cytogenetics:</b>	X 3.68 cM
<b>Gene Summary:</b>	Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures. Reduces the relative abundance of microRNAs, when overexpressed. Enhances phosphorylation of translation initiation factors and active polysome formation.[UniProtKB/Swiss-Prot Function]