

Product datasheet for **MG226880**

Hrh2 (NM_001010973) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hrh2 (NM_001010973) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Hrh2
Synonyms:	H2R; HH2R
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG226880 representing NM_001010973 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCCCAATGGCACGGTTCATTCTGCTGTTGGACTCTATTGCATTGAAAGTCACCATCAGTGTGG
TCCTCACCACCCTCATCTTCATCACTGTTGCTGGCAATGTGGTCGTCGCTGGCTGTGAGCTTGAATCG
TCGGCTCCGCAGTCTGACCAATTGCTTCATTGTGTCCCTGGCAGCCACTGACCTGCTTCTGGCCTCCTG
GTGATGCCCTTCTGCCATTTACCAGTTGTCTTCAAGTGGAGCTTTGGCCAGGTCTTCTGCAACATCT
ACACCAGCCTGGATGTGATGCTCTGCACAGCTTCCATCCTCAACCTCTTCATGATCAGCTTGACCCTA
CTGCGCTGTACGGACCCACTGCGGTACCCCGTGTGGTACCCCTGTTGCGGTTGCCATCTCTTTGGTC
TTTATTTGGGTCAATTTCCATCACCTCTCCTTCTCTCTATTACCTAGGATGGAACAGCAGAAATGGGA
CCAGAGGGGGCAATGACACCTTCAAGTGCAAAAGTGCAGGTCAACGAGGTATATGGACTGGTGGATGGGAT
GGTCACCTTCTACCTGCCTCTTCTCATCATGTGTGCACCTACTACAGAACTTTAAGATCGCCAGGGAG
CAAGCCAAGAGGATCAACCACATCAGCTCCTGGAAGGCAGCCACCATCAGAGAACACAAAGCCACCGTGA
CCCTGGCTGCCGTCATGGGAGCATTATCGTCTGCTGGTTTCTTACTTCACTGCCTTCGTTTACCGTGG
ACTGAGAGGGGATGACGCTGTGAACGAGGTGGTCGAGGGTATCGTCTATGGTTGGGTTATGCCAATTCA
GCCCTGAACCCCATCCTGTACGCCACCCTCAACAGAGACTTCCGATGGCGTACCAACAGCTCTTCCACT
GCAAGCTTGAAGCCACAACCTCTCACAAAACCTCCCTGAGGTTGAACAACCTCTGCTGTCCAGGAGCCA
AAGCCGAGAAGGTAGATGGCAGGAGGAGAAGCCCTGAAGCTTCAAGTGTGGAGTGGGACAGAACCTCACA
CACCCCGAGGAAAGCCAGTCAGGACACGACTGAGTCACTCAAGTTGTCTCCTGAGCCTGAGCCTCCTGT
CCTTCACTGGAAGCTGGAACTTGATTATCACAGGAGACCCTTCAACCATCACTACATATCAGTGC
T

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA



[View online >](#)

Protein Sequence: >MG226880 representing NM_001010973
Red=Cloning site Green=Tags(s)

MEPNGTVHSCCLDSIALKVTISVVLTTLIFITVAGNVVCLAVSLNRRLRSLTNCFIVSLAATDLLLGLL
 VMPFSAIYQLSFKWSFGQVFCNIYTSLDVMLCTASILNLFMISLDRYCAVDPLRYPVLVTPVRVAISLV
 FIWVISITLSFLSIHLGWSNRNGTRGGNDTFKCKVQVNEVYGLVDGMVTFYLPLLIMCVTYRIFKIARE
 QAKRINHISSWKAATIREHKATVTLAAVMGAFIVCWFPYFTAFVYRGLRGDDAVNEVEGIVLWLVGYANS
 ALNPILYATLNRDFRMAYQQLFHCKLASHNSHKTSLRLNNSLLSRSQSREGRWQEEKPLKLQVWSGTELT
 HPQGSPVRTRLSSHSSCLLSLSSFIWKLGTWIHRRPFQPSLHISA


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001010973

ORF Size:	1191 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001010973.2 , NP_001010973.1
RefSeq Size:	2789 bp
RefSeq ORF:	1194 bp
Locus ID:	15466
UniProt ID:	P97292
Cytogenetics:	13 28.4 cM
Gene Summary:	The H2 subclass of histamine receptors mediates gastric acid secretion. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.[UniProtKB/Swiss-Prot Function]