

Product datasheet for **MG226917**

Foxp1 (NM_053202) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Foxp1 (NM_053202) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Foxp1
Synonyms:	3110052D19Rik; 4932443N09Rik; AI461938; AW494214
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG226917 representing NM_053202
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGATGCAAGAATCTGGGTCTGAGACAAAAAGTAAACGGATCAGCCATCCAGAACGGTCCAGCGGTGGCA
ACCACTTACTAGAGTGCGGGGCACTTCGTGACACTCGGTCCAACGAGAGGCACCAAGCGGTGGACCTGGG
GGCAGCCGACCTTGCCACGTCCAGCAGCAGCAACAGCAGGCCCTGCAGGTGGCAAGGCAGCTCCTCCTT
CAGCAGCAGCAACAGCAGCAGCAGCAACAGCAGCAGCAGCAGCAGCAACAGCAGCAACAGCAGCAGC
AACACAGCAGCAGCAGCAACAGCAGCAGCAGCAGCAAGTTAGTGGATTAAGTCTCCCAAGAGGAA
TGACAAGCAACCAGCTCTTCAGTTCCCGTGTGAGTGGCTATGATGACACCTCAGTTATCACTCCTCAA
CAAATGCAGCAGATCCTCCAGCAGCAGGTGCTGAGCCCCAGCAGCTCCAGGTTCTCTCCAGCAGCAGC
AGGCCCTCATGCTTCAACAGCAGCTTCAAGAATTTTATAAAAAACAACAGGAACAGTTGCAGCTTCACT
TCTCCAACAGCAACATGCTGGAAAACAGCCGAAAGAGCAGCAGGTGGCTACCCAGCAGTTGGCTTCCAG
CAGCAGCTTCTCCAGATGCAGCAGCTGCAGCAGCAGCACCTCTTGTCTCTCCAGCGCCAAGGCCCTCTAA
CAATTCAGCCGGGGCAGCCTGCCCTTCCCTCCAACCCCTCGCTCAAGGCATGATTCCAACGGAAGTGA
GCAGCTCTGAAAAGAAGTGACAAGTGCCACACTGCAGAGGAAACCACAAGCAGCAACCACAGCAGCCTA
GACCTGACCAGCAGATGTGTCTCGTCTCGGCACCTTCCAAGTCTCCCTAATCATGAACCCGCATGCTT
CTACCAATGGACAGCTCTCGGTCCACACTCCCAAAAGGAAAGCTTGTCCACAGGAGCAGCCCCACAG
CCACCCCTCTATGGACATGGCGTATGCAAGTGGCCAGGCTGTGAGGCGGTTTGTGACGACTTCCAGCC
TTTCTAAAACATCTCAACAGTGAAGTGCAGTGGCAGTGAAGAGCAAGAGCGCCTGCAAGCCATGATGACCCACT
GCATGTGAAGTCTACAGAACCCAAAGTGCCTCCAGCCCTGAATCTGGTATCAAGTGTACCCTCTCC
AAGTCTGCCTCAGAGGCTTCTCCACAGAGCTTACCTCATACTCCAACAACCCCCACCGCCCCCTGACTC
CTGTACCCCAAGGCCCTCCGTATCACCACCACAGCATGCACACGGTGGGACCTATCCGAGGCGGTA
CTCAGACAAATACAACGTGCCATTTCTCAGCAGATATTGCGCAGAACCAAGAATTTTATAAGAACGCG
GAAGTTAGACCACATTTACATATGCATCTTAAATCAGGAGGCCATTCTCGAATCTCCAGAAAAGCAGC
TAACACTAAACGAAATCTATAACTGGTTCACACGAATGTTTGCTTACTCCGACGCAATGCAGCCACGTG
GAAGAATGCAGTGCATATAATCTTAGTCTCCACAAGTGTGTTGTGCGAGTAGAGAACGTTAAAGGGGCA
GTATGGACAGTGGATGAAGTAGAGTTCCAAAACGGAGGCCACAAAAGATCAGTGGTAACCTTCCCTTA
TTAAAACATGCAGAGCAGCCACGCCTACTGCACACCTCTCAATGCAGCTTTACAGGCTTCCATGGCTGA
GAATAGTATACCTCTGTACACTACCGCTTCCATGGGAAATCCCACTCTGGGCAGCCTGGCCAGTGCCATC
CGGGAGGAGCTGAACGGGGCCATGGAGCACACCAACAGCAACGAGAGTGACAGCAGTCCAGGCAGATCCC
CTATGCAAGCTGTGCACCCATACACGTCAAAGAAGAACCCTCGACCCCGAGGAAGCTGAAGGCCCTCT
GTCCTTAGTGACAACAGCCAACCAACAGTCCAGATTTTGACCATGCAGAGATTACGAAGACGAACCAAGTA
AATGAGGACATGGAG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG226917 representing NM_053202
 Red=Cloning site Green=Tags(s)

```
MMQESGSETKSNQSAIQNGSSGGNHLLECGALRDTRSNGEAPAVDLGAADLAHVQQQQQQALQVARQLLL
QQQQQQQQQQQQQQQQQQQQQQQQQQQQQQQQVSLKSPKRNDKQPALQVPVSVAMMTPOVITPQ
QMQQILQQVLSPPQLQVLLQQQQALMLQQQLQEFYKKQQEQQLQLLQQQHAGKQPKQQVATQQLAFQ
QQLLQMQLQQQHLLSLQRQGLLTIQPGQPALPLQPLAQGMIPTELQQLWKEVTSAHTAEETSSNHSSL
DLTSTCVSSSAPSKSSLIMNPHASTNGQLSVHTPKRESLSHEEHPHSHPLYGHGVCKWPGCEAVCDDFPA
FLKHLNSEHALDDRSTAQCRVQMVVQLELQLAKDKERLQAMMTHLVKSTEPKAAPQPLNLVSSVTL
KSASEASPQSLPHTPTTAPLTPVTQGPSVITTTSMHTVGPIRRYSKYNVPISSADIAQNQEFYKNA
EVRPPFTYASLIRQAILESPEKQLTLNEIYNWFRMFAYFRNAATWKNVVRHNLHLKCFVRVENVKGA
VWTVDEVEFQRRPQKISGNPSLIKMQSSHAYCTPLNAALQASMAENSIPLYTTASMGNPITGLSASAI
REELNGAMEHTNSNESDSSPGRSPMQAVHP IHVKEEPLDPEEAEGPLSLVTTANHSPDFDHRDYEDEPV
NEDME
```



TRTRPLE - GFP Tag - V


Restriction Sites:

SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shutting:



Plasmid Map:


ACCN: NM_053202

ORF Size: 2115 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_053202.2 , NP_444432.1
RefSeq Size:	2118 bp
RefSeq ORF:	2118 bp
Locus ID:	108655
UniProt ID:	P58462
Cytogenetics:	6 D3
Gene Summary:	<p>Transcriptional repressor. Can act with CTBP1 to synergistically repress transcription but CTPBP1 is not essential (PubMed:11358962, PubMed:14701752). Plays an important role in the specification and differentiation of lung epithelium. Acts cooperatively with FOXP4 to regulate lung secretory epithelial cell fate and regeneration by restricting the goblet cell lineage program; the function may involve regulation of AGR2 (PubMed:11358962, PubMed:22675208). Essential transcriptional regulator of B-cell development (PubMed:16819554). Involved in regulation of cardiac muscle cell proliferation (PubMed:20713518). Involved in the columnar organization of spinal motor neurons. Promotes the formation of the lateral motor neuron column (LMC) and the preganglionic motor column (PGC) and is required for respective appropriate motor axon projections. The segment-appropriate generation of spinal chord motor columns requires cooperation with other Hox proteins (PubMed:18667151, PubMed:18662545). Can regulate PITX3 promoter activity; may promote midbrain identity in embryonic stem cell-derived dopamine neurons by regulating PITX3 (PubMed:20175877). Negatively regulates the differentiation of T follicular helper cells T(FH)s (PubMed:24859450). Involved in maintenance of hair follicle stem cell quiescence; the function probably involves regulation of FGF18 (PubMed:23946441). Represses transcription of various pro-apoptotic genes and cooperates with NF-kappa B-signaling in promoting B-cell expansion by inhibition of caspase-dependent apoptosis. Binds to CSF1R promoter elements and is involved in regulation of monocyte differentiation and macrophage functions; repression of CSF1R in monocytes seems to involve NCOR2 as corepressor. Involved in endothelial cell proliferation, tube formation and migration indicative for a role in angiogenesis; the role in neovascularization seems to implicate suppression of SEMA5B. Can negatively regulate androgen receptor signaling (By similarity).</p> <p>[UniProtKB/Swiss-Prot Function]</p>