

Product datasheet for **MG227017**

Itgb6 (NM_021359) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Itgb6 (NM_021359) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Itgb6
Synonyms:	2210409C20Rik; 4831415H04Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG227017 representing NM_021359
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGGGATTGAGCTGGTCTGCCTGTTTCTGCTACTTCTAGGAAGGAATGATCACGTCCAAGTGCGCTGTG
 CCTGGGGTGGTGCAGAAAGCTGCTCAGACTGCCTGCTCACAGGACCCACTGCGCCTGGTCTCCAGGA
 GAATTTACCCACCTGTCTGGAGCTGGCGAGAGGTGCGACACCCAGCGAATCTTTAGCCAAAGGATGT
 CAACTACCCTTCAATTGAAAACCTGTCTCCCGCATAGAAGTCCTTCAAAATAAGCCTCTCAGCGTGGGCC
 GACAGAAGAACAGTTCTGACATTGTTTCTGATTGCTCCTCAAAGCTTGGTTCTTAAATTGAGACCAGGGCG
 TGAGCAGACTCTGCAAGTGCAGGTCCGCCAACTGAAGATTACCCAGTAGATCTGTATTACCTCATGGAC
 CTCTCCGCTCTATGGATGACGACCTCAACACCATCAAAGAGCTTGGCTCCCGGCTGGCCAAAGAGATGT
 CCAAATAACAGCAACTTTAGACTGGGCTTTGGCTCTTTTGTGAAAAGCCGTTTCTCCTTTTATGAA
 AACACCCAGAGGAAATCACCAACCTTGCAGTAGTATCCCCTATTTCTGCTTACCTACATTTGGATTC
 AAGCACATTTTGGCATTAACTGATGATGCCGAGAGATTCAATGAAATTGTGAGGAAACAGAAAGATTTCTG
 CTAATATTGACACACCTGAAGGTGGATTTGATGCAATCATGCAGGCTGCTGTGTGTAAGGAAAAGATTGG
 CTGGCGCAATGACTCGCTCCACCTCCTGGTTTTTGTGAGTGATGCCGATTCTCATTTTGGAAATGGACAGC
 AAGCTGGCAGGCATTGTCAATCCAATGATGGACTCTGTCACTTGGACCACAGGAATGAATATTCATGT
 CAACTGTCTTGGAAATATCCAATATCGGCCAACTCATTGATAAAGTGGTACAAAACACGTGCTACTGAT
 CTTTGCAGTACCCAAGAACAAGTCCATCTGTATGAGAACTATGCAAACTCATTCTGGAGCAACCGTG
 GGACTGCTTCAGAAGGATTCTGGGAACATCCTTACGTGATCATCTCCGTTATGAAGAATCGCGTCTG
 AGGTGGAAGTGGAAAGTGTAGGGGACACAGAAGGACTCAACCTGTCTTTACAGCTCTCTGTAACAATGG
 TGTCTCTTCCACACCAAAAAGAAATGCTCCACATGAAAGTTGGAGACACAGCATCCTTCAATGTGACT
 GTGAGCGTATCCAAGTGTGAGAAAAGAAGCAGGAACCTCATCATAAAGCCAGTGGGGCTGGGGGACACCC
 TAGAAATACTCGTCAGTGCAGAATGTGACTGCGACTGCCAGAGAGAAATAGAAACTAACAGCTCTAAGTG
 TCACAATGGGAATGGCTCCTTCCAGTGTGGGTGTGTACCTGCAACCCCGGCCACATGGGTCTCACTGC
 GAGTGTGGTGAGGATATGGTGAACAGGATTCCTGCAAGGAGTCCCGGGTACCCTTATGCAGTGGAA
 GGGGTGACTGCTATTGTGGCAGTGCATCTGCCACTTATCGCCGTACGGAAGCATCTACGGACCTTACTG
 CCAGTGTGACAATTTCTCTGTCTGAGACACAAAGGCCTTCTCTGTGGAGATAACGGTGACTGTGACTGT
 GCGGAGTGTGTGCCGGATGGCTGGACAGGTGAATACTGTAAGTGTACAACCAACAGGGACTCGTGTA
 CATCTGAGGATGGAGTGTGTGCAAGTGGAGTGGGGACTGTGTCTGTGGCAAGTGTGTCTGCAGAAACCC
 TGGAGCCTCAGGACCCACCTGTGAACGCTGTCTACCTGTGGCGACCCCTGTAAGTCTAAACGGAGTTGC
 ATCGAGTGTACCTGTCTGCAGATGGCCAGGCCAAGAAGAGTGTGCTGACAAGTGTAAAGCCATTGGTG
 CCACCATCAGTGAAGGAAATTTTCAAAGGATACTTCTGTCTCCTGCTCTCTACAAGGAGAAAATGAATG
 CCTTATTACATTCCTAATAACTACGATAATGAAGGAAAACCATCATTACAACATCAATGAAAAAGAC
 TGCCCCAAACCTCAAACATTCATGATCATGTTGGGGTGTCACTGGCGATCCTGCTCATCGGAGTTG
 TGCTACTGTGATTTGGAAGCTGCTGGTATCATTTTATGACCGGAAGGAAAGTTGCTAAATTTGAAGCAGA
 ACGCTCTAAGGCCAAGTGGCAAACGGGAACCAATCCTCTGTACCGAGGTTCCACCAGCACTTTTAAGAAC
 GTGACCTACAAGCACAGGGAAGCACAAGCAGGCCTTTCCTCAGATGGG

ACGCGTACGCGGCCGCTCGAG – GFP Tag – GTTTAA

Protein Sequence: >MG227017 representing NM_021359
 Red=Cloning site Green=Tags(s)

```


MGIELVCLFLLLLGRNDHVQGGCAWGAESCSDCLLTGPHCAWCSQENFTHLSGAGERCDTPANLLAKGC
QLPFIENPVSRIEVLQNKPLSVGRQKNSSDIVQIAPQSLVLKLRPGREQLQVQVRQTEDYVVDLYLMD
LSASMDDDLNTIKELGSR LAKEMSKLTSNFRLGFGSFVEKPVSPFMKTTPEEITNPCSSIPYFCLPTFGF
KHILPLTDDAERFNEIVRKQKISANIDTPEGGFDAIMQAAVCKEIKGWRNDSLHLLVVFSDADSHFGMDS
KLAGIVIPNDGLCHLDHRNEYSMTVLEYPTIGQLIDKLVQNNVLLIFAVTQEQVHLYENYAKLIPGATV
GLLQKDSGNILQLIISAYEELRSEVELEVLGDTEGLNLSFTALCNNGVLFPHQKCSHMKVGD TASFNVT
VSVSNCEKRSRNLIIKPVGLGDTLEILVSAECD CDCQREIETNSSKCHNGNGSFQCGVCTCNPGHMGPHC
ECGEDMVSTDSCKESPGHPSCSGRGDCYCGQCICHLSPYGSYGPYCQCDNF SCLRHKGLL CGDNGDCDC
GECVCRDGTWTEYCNCTNRDCTSEDGVLCSGRGDCVCGKVCVRNPGASGPTCERCPTCGDPCNSKRSC
IECYLSADGQAQEECADKCKAIGATISEEDFSKDTSVSCSLQGENECLITFLITTDNEGKTIHININEKD
CPKPPNIPMIMLVSLA ILLIGVLLCIWKLLVSFHDRKEVAKFEAERSKAKWQTGNPL YRGSTSTFKN
VTYKHREKHKAGLSSDG
  
```

TRTRPLE - GFP Tag - V


Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:



ACCN: NM_021359

ORF Size: 2361 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_021359.3](#), [NP_067334.1](#)

RefSeq Size: 4788 bp

RefSeq ORF: 2364 bp

Locus ID: 16420

UniProt ID: [Q9Z0T9](#)

Cytogenetics: 2 34.81 cM

Gene Summary:

Integrin alpha-V:beta-6 (ITGAV:ITGB6) is a receptor for fibronectin and cytotactin (By similarity). It recognizes the sequence R-G-D in its ligands (PubMed:10025398). ITGAV:ITGB6 acts as a receptor for fibrillin-1 (FBN1) and mediates R-G-D-dependent cell adhesion to FBN1 (By similarity). Integrin alpha-V:beta-6 (ITGAV:ITGB6) mediates R-G-D-dependent release of transforming growth factor beta-1 (TGF-beta-1) from regulatory Latency-associated peptide (LAP), thereby playing a key role in TGF-beta-1 activation (PubMed:10025398). [UniProtKB/Swiss-Prot Function]