

Product datasheet for **MG227107**

B2m (NM_009735) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: B2m (NM_009735) Mouse Tagged ORF Clone
Tag: TurboGFP
Symbol: B2m
Synonyms: beta2-m; beta2m; Ly-m11
Mammalian Cell Selection: Neomycin
Vector: pCMV6-AC-GFP (PS100010)
E. coli Selection: Ampicillin (100 ug/mL)
ORF Nucleotide Sequence: >MG227107 representing NM_009735
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCTCGCTCGGTGACCCTGGTCTTTCTGGTCTGTCTCACTGACCGCCTGTATGCTATCCAGAAAA
CCCCTCAAATCAAGTATACTCACGCCACCCACCGAGAATGGGAAGCCGAACATACTGAAGTCTACGT
AACACAGTTCCACCCGCCTCACATTGAAATCAAATGCTGAAGAACGGGAAAAAATTCCTAAAGTAGAG
ATGTCAGATATGCCTTCAGCAAGGACTGGTCTTTCTATATCCTGGCTCACACTGAATTCACCCCCACTG
AGACTGATACATACGCCCTGCAGAGTTAAGCATGCCAGTATGGCCGAGCCAAGACCGTCTACTGGGATCG
AGACATG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG227107 representing NM_009735
Red=Cloning site **Green**=Tags(s)
MARSVTLVFLVFLVSLTGLYAIQKTPQIQVYSRHPPEKPNILNCYVTQFHPPHIEIQMLKNGKKIPKVE
MSDMSFSKDSFYILAHTEFTPTETDTYACRVKHASMAEPKTVYWRDM

TRTRPLE - GFP Tag - V

Chromatograms: https://cdn.origene.com/chromatograms/ja2528_e11.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_009735

ORF Size: 357 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009735.3](#), [NP_033865.2](#)

RefSeq Size: 858 bp

RefSeq ORF: 360 bp


Locus ID: 12010

UniProt ID: [P01887](#)

Cytogenetics: 2 60.55 cM

Gene Summary: Component of the class I major histocompatibility complex (MHC). Involved in the presentation of peptide antigens to the immune system.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG227107