

Product datasheet for **MG227244**

Tfap2b (NM_009334) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tfap2b (NM_009334) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Tfap2b
Synonyms:	AI606113; AP-2(beta); E130018K07Rik; Tcfap2b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MG227244 representing NM_009334
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCACTCACCTCCTAGAGACCAGGCTGCGATCATGCTCTGGAACTCGTGGAGAATGTCAAGTACGAAG
 ACATCTATGAGGACCGGCACGATGGCGTCCAAGCCATAGCTCGAGACTCTCTCAACTGGGCTCTGTGTC
 CCAAGGACCCTACTCAAGCGCCCCGCGCTGTCCACACGCCATCATCGACTTCCAGCCACCCTACTTT
 CCGCCCCCTACCAGCCGCTACCCTACCACCAGAGTCAAGACCCGTACTCCACGTCAACGATCCCTACT
 CCCTGAACCTCTGCATCAGCCCCAGCAGACCCGTGGGGCAACGGCAACGCCAAGAAGTGGGCTCAGA
 AGCCGGCTCTCTCTGCCCCAGCCCCGGGCAGCCTTGCCCGAGCTCTCCGGCCTTGATCCCCGAAGGGAC
 TACCCTCTGTCCGCGCCGGACGTGCTGCTGCATTCCGCACATCACGGCTGGACGCCGGCATGGGCG
 ACAGCCTCTCGTTGCACGGCCTTGGGCATCCCGGCATGGAAGACGTGCAGTCAGTTGAAGATGCCAATAA
 CAGCGGCATGAACCTATTGGACCAGTCAGTCATTAATAAAGTTCTGTCCCTCCCAAATCTGTGACTTCT
 CTAATGATGAATAAAGATGGCTTCTTGGGAGGAATGTCAGTCAACACCGCGAGGTATTTTGCTCAGTCC
 CGGGCCGTTTATCTCTACTCAGTTCAACTTCAAAGTACAAAGTCACTGTGGGCGAAGTTCAGAGAAGGCT
 CTCGCCCCCTGAATGCCTCAATGCATCTCTCTGGGCGCGCTCCTCAGAAGAGCCAAATCGAAAAATGGG
 GGGAGATCCTTGAGAGAAAGGCTAGAAAAATCGGTTTGAATTTACCCGCGGGCAGGCGCAAAGCAGCAA
 ATGTCACGTTACTCACCTCACTGGTAGAAGGGGAAGCTGTTCACTTAGCTCGGGATTTCCGGGTATATTTG
 TGAACAGAGTTTCTGCTAAAGCCGTGTCGAGTATTTGAACAGACAGCACACGGACCCAGTACCTG
 CACTCCAGAAAGAATATGCTGCTGGCCACCAAGCAACTTTGTAAGAATTTACGGATCTGCTGGCTCAGG
 ACCGGACACCGATCGGAAACAGCAGGCCAGCCCATCCTGGAGCCGGGCATCCAAAGCTGTCTCACGCA
 CTTCAAGTCTCATCACGCAGGCTTCCGGTCCCGGCCATTTGCGCTGCGCTCACGGCCCTGCAGAACTAT
 CTCACCGAGGCGCTCAAAGGCATGGACAAGATGTTCTTGAACAACACCACTAACAGGCACACGTCTGGG
 AAGGCCAGGTAGTAAAAGTGGCGACAAGGAGGAGAAACACAGGAAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence:

>MG227244 representing NM_009334
 Red=Cloning site Green=Tags(s)

MHSPPRDQAAIMLWKLVENVKYEDIYEDRHDGVPSSHSSRLSQLGSVSVQGPYSSAPPLSHTPSSDFQPPYF
 PPPYQPLPYHQSDPYSHVNDPYSLNPLHQPQQHPWQQRQRQEVGSEAGSLLPQPRAALPQLSGLDPRRD
 YHSVRRPDVLLHSAHHGLDAGMGDSLHLHGLGHPGMEDVQSVEDANNMNLDDQSVIKKVPVPPKSVTS
 LMMNKDGFLLGMSVNTGEVFCVPGRLSLLSSTSKYKVTVEVQRRLLSPPECLNASLLGGVLRRAKSKNG
 GRSLRERLEKIGLNLPAARRKAANVTLLTSLVEGEAVHLARDFGYICETEFPKAVSEYLNQHTDPSDL
 HSRKNMLLATKQLCKEFTDLLAQDRTPIGNSRPSPILEPGIQSCLTHFSLITHGFGAPAICAALALQNY
 LTEALKGMDKMFLLNNTTNRHTSGEGPGSKTGDKKEEHRK

TRTRPLE - GFP Tag - V

Restriction Sites:

SgfI-MluI

Cloning Scheme:


ACCN: NM_009334

ORF Size: 1377 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_009334.3](#), [NP_033360.2](#)

RefSeq Size: 6021 bp

RefSeq ORF: 1380 bp


Locus ID: 21419

UniProt ID: [Q61313](#)

Cytogenetics: 1 A3

Gene Summary:

Sequence-specific DNA-binding protein that interacts with inducible viral and cellular enhancer elements to regulate transcription of selected genes. AP-2 factors bind to the consensus sequence 5'-GCCNNNGGC-3' and activate genes involved in a large spectrum of important biological functions including proper eye, face, body wall, limb and neural tube development. They also suppress a number of genes including MCAM/MUC18, C/EBP alpha and MYC. AP-2-beta appears to be required for normal face and limb development and for proper terminal differentiation and function of renal tubular epithelia.[UniProtKB/Swiss-Prot Function]

Product images:


Circular map for MG227244