

## Product datasheet for **MG227263**

### Sirt1 (NM\_019812) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Sirt1 (NM_019812) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Sirt1
Synonyms:	AA673258; S; Si; Sir2; Sir2a; Sir2alpha; SIR2L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>MG227263 representing NM\_019812  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGGACGAGGTGGCGCTCGCCCTTCAGGCCGCGGCTCCCCTTCGCGGGCGGCCGCATGGAGGCCG  
CGTCGCAGCCGGCGGACGAGCCGCTCCGCAAGAGGCCCGCCGAGACGGGCTGGCTCGGGCGCAGCCC  
GGGCGAGCCGAGCGCAGCAGTGGCGCGCGCGCGGGGTGTGAGGCGGCGAGCGCGGGCCCGGGC  
GCGCTGTGGCGGAGGCGGCGAGGGCGGCGGAGCGGAGCGGGAGGCCCGGCGACGGCCGTGGCCG  
GGGACGGAGACAATGGGTCGGCCTGCGCGGGAGCCGAGGGCGGCTGACGACTTCGACGACGACGAGGG  
CGAGGAGGAGGACGAGGCGGCGGCGCAGCGGCGGCGCAGCGATCGGCTACCGAGACAACCTCCTGTTG  
ACCGATGGACTCCTCACTAATGGCTTTCATTCTGTGAAAGTGATGACGATGACAGAACGTCACACGCCA  
GCTCTAGTGACTGGACTCCGCGGCCGCGGATAGGTCCATATACTTTTGTTCAGCAACATCTCATGTTGG  
CACCGATCCTCGAACAATTCTTAAAGATTTATTACCAGAAACAATTCCTCCACCTGAGCTGGATGATATG  
ACGCTGTGGCAGATTGTTATTAATATCCTTTCAGAACCCAAAGCGGAAAAAGAAAAGATATCAATA  
CAATTGAAGATGCTGTGAAGTACTGCAGGAGTGTAAAAAGATAATAGTTCTGACTGGAGCTGGGGTTTC  
TGTCTCCTGTGGGATTCCGACTTCAGATCAAGAGACGGTATCTATGCTCGCCTTGGCGTGGACTTCCCA  
GACCTCCAGACCCTCAAGCCATGTTTGATATTGAGTATTTAGAAAAGACCCCAAGACCATTCTCAAGT  
TTGCAAAGGAAATATATCCCGGACAGTTCAGCCGCTCTGTGTACAAAATTCATAGCTTTGTGAGATAA  
GGAAGGAAAACACTCTCGAAATTATACTCAAAAATAGATACTTGGAGCAGGTTGCAGGAATCCAAAGG  
ATCCTTCAGTGTATGGTTCCTTTGCAACAGCATCTTGCCTGATTTGTAATAACAAAGTTGATTGTGAAG  
CTGTTCTGGAGACATTTTAAATCAGGTAGTTCCTCGGTGCCCTAGGTGCCAGCTGATGAGCCACTTGC  
CATCATGAAGCCAGAGATTGTCTTCTTTGGTGA AAAACTTACCAGAACAGTTTCATAGAGCCATGAAGTAT  
GACAAAGATGAAGTTGACCTCCTCATTGTTATTGGATCTTCTCTGAAAAGTGAGACCAGTAGCACTAATTC  
CAAGTTCTATACCCCATGAAGTGCCTCAAATATTAATAAATAGGGAACCTTTGCCTCATCTACATTTTGA  
TGTAGAGCTCCTTGGAGACTGCGATGTTATAATTAATGAGTTGTGTATAGGCTAGGTGGTGAATATGCC  
AAACTTTGTTGTAACCCTGTAAGCTTTAGAAAATTAAGAAAACCTCCAGCCCAAAAAGGAATTGG  
TTCATTTATCAGAGTTGCCACCAACACCTTTCATATTTGGAAGACTCAAGTTCACCTGAAAGAAGTGT  
ACCACAAGACTCTTCTGTGATTGCTACACTTGTAGACCAAGCAACAACAACAATGTTAATGATTTAGAA  
GTATCTGAATCAAGTTGTGTGGAAGAAAACCAAGAAGTACAGACTAGTAGGAATGTTGAGAACATTA  
ATGTGGA AAAATCCAGATTTTAAGGCTGTTGGTTCCAGTACTGCAGACAAAATGAAAAGAACTTCAGTTGC  
AGAAAACAGTGAGAAAATGCTGGCCTAATAGACTTGCAAAGGAGCAGATTAGTAAGCGGCTTGAGGGTAAT  
CAATACCTGTTTGTACCACCAATCGTTACATATTCACGGTGCTGAGGTATACTCAGACTCTGAAGATG  
ACGCTTTGTCCTTAGTTCCTGTGGCAGTAACAGTGACAGTGGCAGATGCCAGAGTCCAAGTTTGAAGA  
ACCCTTGAAGATGAAAGTGA AATTGAAGAATTCTACAATGGCTTGAAGATGATACGGAGAGGCCCGAA  
TGTGCTGGAGGATCTGGATTTGGAGCTGATGGAGGGGATCAAGAGGTTGTTAATGAAGCTATAGCTACAA  
GACAGGAATTGACAGATGTA AACTATCCATCAGACAAATCA

**ACGCGT**ACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG227263 representing NM\_019812  
 Red=Cloning site Green=Tags(s)


MADEVALALQAAGSPSAAAAAMEAASQPADEPLRKRPRRDGPGGLGRSPGEPSSAAVAPAAAGCEAASAAAPA  
 ALWREAAGAAAASAEREAPATAVAGDGDNGSGLRREPRAADDFDDDEGEEDEAAAAAIIGYRDNLLL  
 TDGLLTNGFHSCESDDDDRTSHASSDWTTPRPRIGPYTFVQQLMIGTDPRTILKDLLPETIPPELDDM  
 TLWQIVINILSEPPKRRKRKDIINTIEDAVKLLQECKKIIIVLTGAGVSVSCGIPDFRSRDGIYARLAVDFP  
 DLPDPQAMFDIEYFRKDRPFFKFAKEIYPGQFQPSLCHKFIALSDKEGKLLRNYTQNIDTLEQVAGIQR  
 ILQCHGSFATASCLICKYKVDCEAVRGDIFNQVPRCPRCPADEPLAIMKPEIVFFGENLPEQFHRAMKY  
 DKDEVDLLIVIGSSLKVRPVALIPSSIPHEVPQILINREPLPHLHFDVELLGDVCDVIINELCHRLGGEYA  
 KLCCNPVKLSEITEKPPRPQKELVHLSELPTPLHISEDSSSPERTVPQDSSVIATLVDQATNNNVNDLE  
 VSESSCVEEKQEVQTSRNVENINVENPDFKAVGSSTADKNERTSVAETVRKCPWPNRLAKEQISKRLGN  
 QYLFVPPNRYIFHGAEVYSDSEDDVLSSSSCGSNSDSGTCQSPSLEEPLDESEIEEFYNGLEDTERPE  
 CAGGSGFGADGGDQEVVNEAIATRQELTDVNYPSDKS

TRTRPLE - GFP Tag - V

Chromatograms: [https://cdn.origene.com/chromatograms/ja1577\\_h03.zip](https://cdn.origene.com/chromatograms/ja1577_h03.zip)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM\_019812

ORF Size: 2211 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_019812.3](#)

**RefSeq Size:** 3920 bp

**RefSeq ORF:** 2214 bp


**Locus ID:** 93759

**UniProt ID:** [Q923E4](#)

**Cytogenetics:** 10 B4

**Gene Summary:** This gene encodes a member of the sirtuin family of proteins, characterized by their deacetylase activity and proposed role in longevity. The encoded protein regulates gene expression in a wide range of cell and tissue types through its NAD<sup>+</sup>-dependent deacetylation of histones, transcription factors and transcriptional coactivators. Brain-specific overexpression of this gene has been shown to result in increased median lifespan. Viability of homozygous knockout mice for this gene varies with strain background. Homozygous knockout mice of strains that do not exhibit embryonic lethality are sterile and have a reduced lifespan. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]

Product images:



Circular map for MG227263