

Product datasheet for **MG227297**

Elavl4 (NM_010488) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Elavl4 (NM_010488) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Elavl4
Synonyms:	Elav; Hud; PNEM
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)
ORF Nucleotide Sequence:	>MG227297 representing NM_010488 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTGGAAATGGCTTGAAGATGATAATTAGCACCATGGAGCCTCAGGTGTCAAATGGACCGACATCCA
ATACAAGCAATGGACCTCCAGCAACAACAGAACTGTCTTCTCCCATGCAGACAGGCGCCGACAGA
TGACAGCAAAACCAACCTCATCGTCAACTATTTACCCAGAATATGACCCAAGAGGAATTCAGGAGTCTC
TTTGGGAGCATTGGTGAATAGAATCCTGCAAACCTCGTGAGAGACAAAATCACAGGACAGAGTTAGGGT
ATGGATTTGTTAACTATATTGATCCAAAGGATGCAGAGAAAGCCATCAACACTTTAAATGGACTCAGACT
CCAGACAAAACCATAAAGGTGTCTACGCCCGCCAAGCTCGGCCTCGATCAGGGATGCTAACCTGTAT
GTTAGTGGCCTTCCCAAGACCATGACCCAGAAGGAGCTGGAGCAGCTTTTCTCTCAATACGGTCGCATCA
TCACCTCACGCATCCTGGTTGATCAAGTCACAGGGGTCTCTAGAGGGGTGGGATTCATCCGCTTTGATAA
GAGGATTGAGGCAGAAGAAGCCATCAAAGGACTGAATGGCCAGAAGCCCAGCGGTGCTACAGAACCATT
ACTGTGAAGTTTGCCAACAACCCAGCCAGAAGTCCAGCCAGGCCCTGCTCTCCAGCTGTACCAGTCCC
CCAACAGGCGATACCCTGGCCCTCTCCACCACCAGGCTCAAAGATTCAGGCTGGACAATTTGCTTAATAT
GGCCTATGGCGTAAAGAGACTGATGTCTGGACCAGTCCCCCTTCTGCTTGTCCCCCAGGTTCTCCCCA
ATCACCATTGACGGGATGACAAGTCTTGTGGGAATGAACATCCCTGGTCACACAGGACAGGCTGGTGCA
TCTTCGTCTATAACCTGTCCCCTGATTCTGATGAGAGTGTCTCTGGCAGCTCTTTGGCCCTTTGGCGC
AGTGAACAACGTCAAGGTCATCCGTGACTTCAACACCAACAAGTGAAGGATTTCGGCTTCGTCACCATG
ACCAACTACGATGAGGCAGCCATGGCCATCGCCAGCCTCAATGGCTATCGCCTGGGAGACAGAGTGTTC
AAGTTTCTTTAAACCAACAAGCCACAAATCC

ACGGTACGCGGCCGCTCGAG - GFP Tag - GTTAA



[View online »](#)

Protein Sequence: >MG227297 representing NM_010488
Red=Cloning site Green=Tags(s)

MEWNLKMIISTMEPQVSNPSTNTSNGPSSNNRNCPSMQTGAATDDSKTNLIVNYLPQNMTQEEFRSL
 FGSIGEIESCKLVRDKITGQSLGYGFVNYIDPKDAEKAINTLNLRLQTKTIKVSYPSSASIRDANLY
 VSGLPKTMTQKELEQLFSQYGRITTSRILVDQVTGVSIRGVGFIKDFKRIEAEAAIKGLNGQKPSGATEPI
 TVKFANNPSQKSSQALLSQLYQSPNRRYPGPLHHQAQRFRLDNLLNMAYGVKRLMSGPVPSSACPPRFSP
 ITIDGMTSLVGMNIPGHTGTGWCIFVYNLSPDSDSVLWQLFGPFGAVNNVKVIRDFNTNKCKGFGFVTM
 TNYDEAAMAIASLNGYRLGDRVLQVSFKTNKAHKS

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_010488

ORF Size: 1155 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_010488.4](#), [NP_034618.2](#)

RefSeq Size: 4094 bp

RefSeq ORF: 1158 bp


Locus ID: 15572

UniProt ID: [Q61701](#)

Cytogenetics: 4 51.42 cM

Gene Summary: May play a role in neuron-specific RNA processing. Protects CDKN1A mRNA from decay by binding to its 3' UTR. Binds to AU-rich sequences (AREs) of target mRNAs, including VEGF and FOS mRNA (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MG227297