

Product datasheet for **MG227373**

Arnt (NM_009709) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Arnt (NM_009709) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Arnt
Synonyms:	bHLHe2; D3Erttd557e; Drnt; ESTM42; Hif1b; mKIAA4051; W08714
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MG227373 representing NM_009709
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGCGGCGACTACAGCTAACCCAGAAATGACATCAGATGTACCATCGCTGGGTCCACCATTGCTTCTG
GAAACCTTGGACCTGGGATTCAAGGTGGAGGAGCTGTTGTACAGAGGGCTATTAAGCGACCGTCAGGGCT
GGATTTTATGATGAAAGTGAAGTGAACACTAAATTTTTGAGATGCGATGATGACCAGATGTGTAATGAC
AAAGAGCGGTTTCCAGGAAAATCATAGTGAATAGAACGGCGCGACGGAACAAGATGACAGCTTACA
TCACAGAACTGTCAGACATGGTACCTACATGTAGTGCCCTGGCTCGAAAACCAGACAAGCTAACCATCTT
ACGCATGGCCGTTTCTCACATGAAGTCTTGAGGGAACTGGCAACACATCTACTGATGGCTCTACAAG
CCATCTTCTCACTGATCAGGAACTGAAACATTTGATCTTGAGGCAGCAGATGGCTTTCTGTTTATTG
TCTCCTGTGAGACTGGACGGTGGTGTATGTCTCTGACTCAGTACTCCCGTTTTGAACCAGCCACAGTC
TGATGGTTCCGGAGCACACTGTATGATCAGGTGCACCCAGATGATGTGGATAAACTTCGAGAGCAGCTC
TCTACATCAGAAAATGCCCTAACAGGGCGGGTCCCTGGATCTGAAGACTGGAACAGTGAAAAGGAAGGCC
AGCAGTCTTCCATGAGGATGTGCATGGGCTCACGAAGGTCGTTTCACTGCGCCATGAGGTGTGGTACTAG
CTCCGTGGACCCTGTTTCCATGAATAGACTGAGCTTTTTGAGGAACAGATGCAGGAATGGGCTTGGCTCT
GTGAAGGAAGGAGAACCTCACTTTGTGGTAGTCCACTGCACAGGCTACATCAAGGCTGGCCACCAGCAG
GTGTCTCCCTCCCAGATGATGACCCAGAGGCTGGCCAGGGGAGCAAAATTCGCCTAGTGGCCATTGGCAG
GCTGCAGGTAAGTCTCCCACTGTACAGACATGAGTAACATTTGTCAGCCAACAGATTTCACTCC
CGACACAACATTGAAGGGATATTCATTTTGTAGACCATCGTTGTGTGGTACTGTTGGCTACCAGCCAC
AGGAGCTTTAGGGAAGAATATTGTAATTTTGTATCCTCTGAAGACCAACAACCTCTAAGAGACAGCTT
TCAGCAGGTGGTAAAATTAAGGTCAGGTGCTGTCCGTCATGTTCCGATTCCGATCTAAGACCCGAGAA
TGGCTGTGGATGAGAACGAGCTCCTTTACCTTCCAAAACCCTTATTCAGATGAAATTGAGTATATTACT
GCACCAACACCAATGTGAAGAACTTAGCCAGGAACCACGGCCTACACTGTCCAACACCATCCCAAGGTC
ACAGCTAGGTCCGACAGCCAATTTATCCCTAGAGATGGGTACAGGGCAGCTGCCATCCAGGCAGCAGCAG
CAGCAGCACACAGAACTGGATATGGTACCAGGAAGAGATGGGCTGGCCAGCTATAATCATTCCCAGGTTT
CTGTCCAGCCTGTGGCAAGTGCAGGATCAGAACACAGCAAGCCCTTGAGAAGTCAGAAGGTCTCTTTGC
ACAGGACAGAGATCCAAGTTTCCAGAAATCTATCCCAGCATCACTGCAGATCAGAGTAAAGGCATCTCC
TCCAGCACTGTCCCTGCCACCAACAGCTGTTCTCCCAGGCAGCTCATTCCCTCCTAACCCCGGCCGG
CAGAGAATTTAGGAATAGTGGTCTTACCCTCCTGTAACCATTGTCCAGCCATCATCTTCTGCAGGGCA
GATACTGGCCAGATTTACGTCCTCAACCCTGCCAGGGATCAGCGCCGACCTGGACCTCTAGCTCC
CGCCAGGCTTTGCCGCCAGCAGGTGCCACCCAGGCTACAGCCAAGACTCGTTCTTCCCAATTTGGTG
TGAACAACTTTACAGACTTCTTCTCCTTCAAGTGTATGTCTTCCGGGTGCTCCCACTGCCTCATCTGG
TACTGTGCTACCCTGCTCTCCCAACCGTGGCTCCAATTTCTCCTGAGACTGGACAGACCACAGGA
CAGTTCAGGCCCCGACAGCAGAGGGCGTGGGTGTCTGGCCACAGTGGCAGGGCCAGCAGCCCCATCATC
GGTCTAGTTCAGTGAGCAGCATGTTTACAGACACAAGCACAAGCACCTAGCCAGCCTGAGGTCTTTCA
AGAAATGCTGTCCATGCTGGGAGACCAAGCAACACCTACAACAATGAAGAATTTCTGATCTAACTATG
TTTCCCCCTTTTCCGAA

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG227373 representing NM_009709
Red=Cloning site Green=Tags(s)


```

MAATTANPEMTSDVPSLGPITIASGNPQGIQGGGAVVQRAIKRRSGLDFDDEVEVNTKFLRCDDQMCND
KERFARENHSEIERRRRNKMTAYITELSDMVPTCSALARKPKLTLILMAVSHMKSLRGTGNTSTDGSYK
PSFLTQELKHLILEAADGFLFIVSCETGRVVYVSDSVTPVLNQPSQSEWFGSTLYDQVHPDDVDKLRQL
STSENALTRVRLDLKTGTVKKEGQQSSMRMCMGSRSSFICRMRCTSSVDPVSMNRLSFLNRNCRNGLGS
VKEGEPHFVVVHCTGYIKAWPPAGVSLPDDPEAGQGSKFCLVAIGRLQVTSSPNCTDMSNICQPTEFIS
RHNIEGIFTFVDHRCVATVGYQPQELLGKNIVEFCHPEDQQLLRDSFQQVVKLKGQVLSVMFRFRSKTRE
WLWMRTSSFTFQNPYSDEIEYIICTNTNVKNSSQEPRTLNTIPRSQLGPTANLSLEMGTGQLPSRQQQ
QQHTELDMPGRDGLASYNHSQVSVQPVASAGSEHSKPLEKSEGLFAQDRDPRFPEIYPSITADQSKGIS
SSTVPATQQLFSQSSFPNPRPAENFRNSGLTPPVTIVQPSSSAGQILAQISRHSNPAQGSAPTWTSSS
RPGFAAQVPTQATAKRSSQFVNNFQTSSTFSAMSLPGAPTASSGTAAYPALPNRGSNFPETGQTTG
QFQARTAEVGVWPQWQGGQPHRRSSSEQHVQQTQAQAPSQPEVFQEMLSMLGDQSNTYNNEEFPDLTM
FPPFSE
    
```

TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_009709


ORF Size: 2328 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_009709.4, NP_033839.2</u>
RefSeq Size:	4325 bp
RefSeq ORF:	2331 bp
Locus ID:	11863
Cytogenetics:	3 40.74 cM
Gene Summary:	<p>Required for activity of the Ah (dioxin) receptor. This protein is required for the ligand-binding subunit to translocate from the cytosol to the nucleus after ligand binding. The complex then initiates transcription of genes involved in the activation of PAH procarcinogens (By similarity). The heterodimer binds to core DNA sequence 5'-TACGTG-3' within the hypoxia response element (HRE) of target gene promoters and functions as a transcriptional regulator of the adaptive response to hypoxia (PubMed:26245371, PubMed:27782878). The heterodimer ARNT:AHR binds to core DNA sequence 5'-TGCGTG-3' within the dioxin response element (DRE) of target gene promoters and activates their transcription (PubMed:28602820). [UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MG227373