

Product datasheet for **MG227387**

Abcg8 (NM_026180) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Abcg8 (NM_026180) Mouse Tagged ORF Clone
Tag:	TurboGFP
Symbol:	Abcg8
Synonyms:	1300003C16Rik; AI114946; sterolin-2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-AC-GFP (PS100010)
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**ORF Nucleotide
Sequence:**

>MG227387 representing NM_026180
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCTGAGAAAACCAAGAAGAGACCCAGCTGTGGAATGGGACTGTACTTCAGGATGCTTCGGGCCTCC
 AGGACAGCTTGTCTCCTCGAAAAGTGACAACAGTCTGTACTTCACCTACAGTGGTCAGTCCAACACTCT
 GGAGGTACAGATCTCACCTACCAGGTGGACATCGCCTCTCAGGTGCCTTGGTTTGGAGCAGCTGGCTCAG
 TTCAAGATACCCTGGAGGTCTCATAGCAGCAAGACTCCTGTGAGCTGGGCATCCGAAATCTAAGCTTCA
 AAGTGAGGAGTGGACAGATGCTGGCCATCATAGGGAGCTCAGGCTGCGGGAGAGCCTCACTACTCGACGT
 GATCACAGGCAGAGGCCACGGTGGCAAGATGAAATCAGGACAAATTTGGATAAATGGGCAACCCAGTACG
 CCTCAGCTGGTGAGGAAGTGCCTGCGCATGTGCGGCAGCATGACCAACTGCTGCCAACCTGACCGTCA
 GAGAGACCCTGGCTTTCATTGCCAGATGCGCCTGCCAGGACCTTCTCCAGGCCACGCTGACAAACG
 GGTGGAAGACGTAAATCGCCGAGCTGCGGCTGCGGCAGTGCGCCAACACCAGAGTGGGCAACCGTATGTA
 CGTGGGGTGTCCGGGGTGTAGCGCCGACAGTGTGAGCATTGGGGTGCAGCTCCTGTGGAACCCAGGAATCC
 TCATTCTGGATGAACCCACTTCTGGCCTCGACAGCTTACAGCCACAATCTGGTGAACCTTGTCCCG
 CCTGGCCAAGGGCAACAGGCTGGTGTCTATCTCCCTCCACCAGCCTCGCTCTGACATCTTCAGGCTATTT
 GACCTGGTCTTCTGATGACATCTGGCACCCCTATCTACCTGGGGGCGGCGCAGCAAATGGTGCAGTACT
 TCACATCCATTGGCCACCCTTGTCTCGTATAGCAACCTGCGGACTTCTACGTGGACTTGACCAGCAT
 CGACAGACGCAGCAAAGAACGGGAGGTGGCCACCGTGGAGAAGGCACAGTCTCTTGACGCCCTGTTCCCTA
 GAAAAAGTACAAGGCTTTGATGACTTTCTGTGAAAGCTGAGGCAAAGGAACTCAACACAAGCACCCACA
 CAGTCAGCCTGACCCCTCACACAGGACACTGACTGTGGGACTGCTGTTGAGCTGCCGGGATGATAGAGCA
 GTTTTCCACCCTGATCCGTGTCAGATTTCCAATGACTTCCGGGACCTGCCACGCTGCTATTCATGGG
 TCGGAAGCCTGCCTGATGTCCCTCATCATTGGCTTCTTTACTACGGCCATGGGGCAAGCAGCTCTCCT
 TCATGGACACAGCAGCCCTCCTTTCATGATAGGGGCGCTCATTCTTTCAATGTCATCCTGGATGTCGT
 CTCCAAATGTCACCTCGGAGAGGTCAATGCTGTACTATGAGCTGGAAGACGGGCTGTACACTGCTGGTCTC
 TATTTCTTTGCCAAGATCCTAGGAGAATTGCCGGAGCACTGTGCCTACGTATCATCTACCGATGCCCA
 TCTACTGGCTGACAAACCTGCGGCCGTGCTGAGCTCTCCTTCTACACTTCTGCTCGTGTGGTTGGT
 GGTCTTCTGCTGCAGGACCATGGCCCTGGCTGCCTCTGCCATGCTGCCACCTTCCACATGTCCTCCTTC
 TTCTGCAATGCCCTCTACAACCTTCTACCTTACTGCCGGCTTCATGATAAACTTGACAACTGTGGA
 TAGTGCCTGCATGGATCTCCAAGCTGTCGTTCTCCGGTGGTGTCTCCTGGGGCTGATGCAGATTAATT
 TAATGGACACCTTTACACCACACAAATCGGCAACTTACCTTCTCCATCCTCGGAGACACGATGATCAGT
 GCCATGGACCTGAACTCGCATCCACTCTATGCGATCTACCTATTGTCATCGGCATCAGTACGGCTTCC
 TGTTCTGTACTATCTATCCTTGAAGCTCATCAAACAGAAGTCAATCAAGACTGG

ACGCGTACGCGGCCGCTCGAG - GFP Tag - GTTTAA

Protein Sequence: >MG227387 representing NM_026180
 Red=Cloning site Green=Tags(s)


MAEKTKEETQLWNGTVLQDASGLQDSLFSSES DNSLYFTYSGQSNTLEVRDLTYQVDIASQVPWFEQLAQ
 FKIPWRSHSSQDSCELGIRNLSFKVRSQMLAII GSSGCGRASLLDVITGRGHGGMKSGQI WINGQPST
 PQLVRCVAHVRQHDQLLPNLT VRETLAFIAQMRLPRTFSQAQRDKRVEDVIAELRLRQCANTRVGNTYV
 RGVSGGERRRVSIGVQLLWNP GILILDEPTSGLDSFTAHLVTTL SRLAKGNRLVLI SLHQPRSDIFRLF
 DLVLLMTSGTPIYL GAAQOMVQYFTSIGHPCPRYSNPADFYVDLTSIDRRSKEREVATVEKAQSLAALFL
 EKVQGFDDFLWKA EAKELNTSTHTVSLTLTQDTCGTAVELPGMIEQFSTLIRRQISNDFRDLPTLLIHG
 SEACLSLIIIGFLYYGHGAKQLSFMDTAALLFMIGALIPFNVILDV VSKCHSERSMLYEELEDGLYTAGP
 YFFAKILGELPEHCAYV IYAMP IYWL TNLRPVELFLLHFLLVWL VVFCRTMALAASAMLP TFHMSSF
 FCNALYNSFYLTAGFMINL DNLWIVPAWISKLSFLRWCF SGLMQIQFNHLYTTQIGNFTFSILGDTMIS
 AMDLNHPLYAIYLIVIGISYGF LFLYYLSLKL IKQKSIQDW


TRTRPLE - GFP Tag - V

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



Plasmid Map:


ACCN: NM_026180

ORF Size: 2019 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026180.2](#)

RefSeq Size: 3674 bp

RefSeq ORF: 2022 bp

Locus ID: 67470

UniProt ID: [Q9DBM0](#)

Cytogenetics: 17 55.02 cM

Gene Summary: ABCG5 and ABCG8 form an obligate heterodimer that mediates Mg(2+)- and ATP-dependent sterol transport across the cell membrane (PubMed:16352607, PubMed:16867993, PubMed:18402465). Plays an essential role in the selective transport of the dietary cholesterol in and out of the enterocytes and in the selective sterol excretion by the liver into bile (PubMed:12444248, PubMed:14504269, PubMed:14657202, PubMed:25378657). Plays an important role in preventing the accumulation of dietary plant sterols in the body (PubMed:12444248, PubMed:14657202). Required for normal sterol homeostasis (PubMed:12444248, PubMed:14657202). The heterodimer with ABCG5 has ATPase activity (PubMed:16352607, PubMed:16867993).[UniProtKB/Swiss-Prot Function]